

# DENTAL LABORATORY SALE

## FALL 2025

Offers Valid September 28 – December 27, 2025

PROMO CODE **FPD**

### Image APN®

Highly Wear Resistant Resin Teeth



**SAVE OVER \$1,200**  
**& Other Bundles Available!**  
Image APN Bundle Packages  
Pages 54-55

### HENRY SCHEIN®



**SAVE OVER \$350**  
Maxima XL2 Handpiece  
Page 44

### SHINING 3D DENTAL



**EXCLUSIVELY AVAILABLE**  
**by ZAHN DENTAL**  
Aoralscan Lync Shining 3D  
Pages 4-5

### formlabs | dental



**SAVE OVER \$1,000**  
Form 4B Printer Bundles  
Pages 10-11

### myerson®



**NOW AVAILABLE!**  
Trusana1 Denture System  
Page 9

### HENRY SCHEIN®



**SAVE OVER \$75**  
**PLUS**  
**Buy 1, Get 2 Offer**  
Henry Schein Vest  
Page 42

### Zirlux®



**SAVE 15%**  
Zirlux 3D Print Resins  
Page 3

### AMANN GIRRBACH



**Q4 SPECIAL PRICING**  
Amann Girschbach Matron  
Page 24

## WE ARE YOUR TOTAL SOLUTION PROVIDER

[www.zahndental.com](http://www.zahndental.com) | 1.800.496.9500

California customers please use website when ordering for Prop 65 information.

25ZS7855

## INDEX

3D Printing & Accessories.....	9–19
3Shape .....	7
Abrasives .....	59
Abrasive Blasters .....	36
Acrylics & Denture Accessories .....	62
Amann Girrbach.....	24–25
Asiga.....	12–13
AutoFinish.....	57
Coltene .....	47
Curing Units .....	16–17
Dekema .....	29
Dental Creations.....	52
Dental Direkt.....	32
Dentsply Sirona .....	60
Denture Equipment.....	63
Dust Collectors.....	27, 38
Exocad.....	8
Formlabs .....	10–11
Fume Extractors.....	18
Fume Hoods .....	38
Furnaces .....	48
GC America.....	43
General Equipment.....	40
Handpieces.....	48
Henry Schein Brand Products.....	44
Henry Schein Brand Teeth .....	54–55
Henry Schein Financial Services.....	20
Henry Schein Prosthetic Design Services ....	6
Henry Schein ServiceFirst.....	33
Infection Control.....	46
Investments.....	42, 50
Keystone Industries.....	56
Kerr.....	45
Mestra .....	53
Milling Machines.....	22–28
Pac-Dent.....	14
Quatro .....	37
Renfert .....	41, 49, 51
Roland DGShape.....	22–23
Shenpaz .....	30–31
Shining 3D Intra Oral Scanners .....	4–5
Sierra Dental Tools.....	21
Sintering Furnaces.....	29–31
TCS .....	61
Trusana1 .....	9
Vacuum Forming .....	58
Vaniman.....	39
Washing Units.....	19
Zirlux 3D Print Resins .....	3
Zirlux Transitions Zirconia.....	34




**ZAHN**™  
A DIVISION OF  HENRY SCHEIN®

If you think this  
has a lot of  
**dental lab solutions,**  
you should see  
our website:



**zahndental.com**




© 2025 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors.

To guarantee special pricing/offers, where applicable, please include the three-digit promo code when ordering. Henry Schein Dental reserves the right to discontinue this promotion at any time [including, without limitation, if supply or demand for the promotional products is affected due to extraordinary events]. Note: Participating in a promotional discount program (e.g., points, discount redemptions or other special awards) is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the program requirements. For HS promotions: This promotion offers a bundled discount. The cost of any “buy/get” (or “no charge”) items is a discount applied to your other purchases on the same invoice. The reportable value (unit price) of the “no charge” item will be identified on the invoice. To calculate the applicable discount on the other items on the same invoice, divide the total value (unit price(s)) of any no-charge item(s) by the total value (extension price(s)) of all other items purchased on the same invoice, and apply the resulting percentage equally to each item to determine net prices. For manufacturer promotions: The prices stated in the invoice or statement may reflect a discount or be subject to a rebate. You must fully and accurately report this stated discount price, or if applicable, any net pricing, after giving effect to any rebates, to Medicare, Medicaid, Tricare and any other federal or state program upon request by any such program. Accordingly, you should retain these records. It is your responsibility to review any agreements or other documents applicable to these prices to determine if they are subject to a rebate.

**15% Off**  
ALL Zirlux Print Resins\*

# Zirlux® 3D Print Resins



Deliver precision and durability for functional restorations.



Validated on industry leading 3D printers & workflows



Made in the USA

ITEM CODE	DESCRIPTION	SIZE	RETAIL	PROMO PRICE*	SAVINGS
921-0033-FPD	Surgical Guide	1 Kg	\$259.99	\$220.99	\$39.00
921-0034-FPD	Ivory Model	1 Kg	\$79.99	\$67.99	\$12.00
921-0035-FPD	Lt Grey Model	1 Kg	\$79.99	\$67.99	\$12.00
921-0037-FPD	Flex Nightguard & Splint	1 Kg	\$429.99	\$365.99	\$64.50
921-0036-FPD	Rigid Nightguard & Splint	1 Kg	\$429.99	\$365.99	\$64.50
921-0032-FPD	Gingiva Mask	0.5 Kg	\$184.99	\$157.99	\$27.75

\*New customers only. Valid 6/29– 12/27/2025



Scan for list of validated **3D printers & workflows**

Exclusively Available from



# Collaborate Better. Communicate Faster. Digitize Smarter.

## **AoralScan Lync**

Exclusively  
Available From  
Zahn Dental

(312-0268)



## Help Your Clinicians Deliver Fast & Precise Digital Impressions

The Shining 3D AoralScan Lync allows your clinicians to digitize their practice efficiently and affordably, streamlining communication between your lab and their practice.



**Lightest-In-Class  
Ergonomics**



**Fast Scanning with  
Outstanding Accuracy**



**Smarter AI-Driven  
Software**

**Contact your Zahn Representative to learn more and get full workflow support.**



# Shining 3D Intraoral & Facial Scanners

**SHINING 3D**  
DENTAL

Innovative digital solutions for every clinical workflow

## Aoralscan Elite Wireless

### Next-Gen Intraoral Scanner with Photogrammetry (IPG)

Ultra-lightweight (194g) with Wi-Fi 6 connectivity, designed for effortless handling and exceptional accuracy in full-mouth and edentulous cases.

- Two-in-one intraoral + photogrammetry scanning
- AI-powered analysis & oral health reports
- Long battery life with easy charging cradle



## Aoralscan 3 Wireless and Wired

### Intraoral Scanner for Daily Efficiency

Autoclavable, ergonomically designed, and equipped with motion sensing to reduce cross-contamination available as a wireless solution and wired.

- AI scanning with auto data filtering
- Long-lasting battery system (3 included)
- Anti-fog tip heating for uninterrupted scans



## MetiSmile

### Facial Scanner for Digital Smile Design

Captures 3D facial data in seconds and auto-aligns with intraoral scans for precise aesthetic planning.

- Ortho simulation tools
- Mandibular trajectory tracking
- Facial feature measurement & comparison



312-0261-FPD	Aoralscan Elite Wireless
312-0262-FPD	Aoralscan Elite Wireless A Set
312-0263-FPD	Aoralscan Elite Wireless B Set
312-0264-FPD	Aoralscan Elite Wireless C Set
312-0265-FPD	Aoralscan Elite Wireless D Set
147-7884-FPD	Aoralscan Elite A Set
147-7885-FPD	Aoralscan Elite B Set
147-7887-FPD	Aoralscan Elite C Set
147-7888-FPD	Aoralscan Elite D Set
312-0257-FPD	Aoralscan 3 Wireless
312-0258-FPD	Aoralscan 3
312-0259-FPD	Metismile
312-0260-FPD	AutoScan DS EX Pro H

**Don't let a lack of CAD  
Designers hold you back.**

**Try two risk-free designs  
at no charge.**




Our team of digital designers have over 100 years of technical experience and a commitment to consulting and supporting your dental lab through the entire case workflow from case submission to case completion.

**Let us design your next  
Fixed, Removable or Implant case.**

**Fast and Reliable Turnaround**  
Same-Day Designs Available

 **HENRY SCHEIN®**  
PROSTHETIC DESIGN SERVICES | LAB

Explore our full list of design services at  
[www.henryscheindesignservices.com](http://www.henryscheindesignservices.com)





## Novux F8 Scanner

### 3Shape

Engineered for efficient, dual-model scanning workflows requiring minimal scanner interaction, F8 enables the user to do more, in less time, and with fewer steps.

Features and benefits:

- Dual-model scanning
- 4 x 5 MP cameras
- Accuracy (ISO 12836): 4 µm
- Scan speed (full arch): 9 seconds
- Scan speed (full arch impression): 45 seconds
- Texture: color
- Scanning strategy: Die-in-model and all-in-one scanning

#20007940

(143-3935-FPD) .....Ea.

**Call**



## Novux E2 Red Dental System Scan Only (No PC)

### 3Shape

- Cameras: 2–5 MP
- Accuracy: 10 µm
- Scan speed (full arch): 24 seconds
- Scan speed (full-arch impression): 72 seconds
- Texture: B/W
- Scanning strategy: Standard

(142-1606-FPD) .....Ea.

**Call**



## Novux E3 Red Dental System Scan Only (No PC)

### 3Shape

- Cameras: 2–5 MP
- Accuracy: 7 µm
- Scan speed (full arch): 18 seconds
- Scan speed (full-arch impression): 64 seconds
- Texture: Color
- Scanning strategy: Standard

(142-1609-FPD) .....Ea.

**Call**



## Novux E4 Red Dental System Scan Only (No PC)

### 3Shape

- Cameras: 4–5 MP
- Accuracy: 4 µm
- Scan speed (full arch): 9 seconds
- Scan speed (full-arch impression): 45 seconds
- Texture: Color
- Scanning strategy: Die in model

(142-1612-FPD) .....Ea.

**Call**



# SUPERCHARGE YOUR CAD DESIGN SKILLS

Daniel Portal, Dental Technician

**DentalCAD® gives you design power, precision and flexibility like never before.**

Envision taking your CAD skills to superpower levels. Streamline communication across the entire workflow with intuitive tools and new AI services. Be the benchmark and design the future of digital dentistry.

**Imagine the CADabilities!**

**[exocad.com/dentalcad](https://exocad.com/dentalcad)**

Territorial coverage of the registered trademarks is available on our webpage: [exocad.com/imprint](https://exocad.com/imprint)

**exocad**

Explore hands-on learning and master the latest digital dentistry workflow at exocad's North American excellence center.

**Scan the code to see the current course lineup.**





# welcome to the future of dentures

## why print when you can jet!

**trusana<sup>1</sup>**  
DENTURE SYSTEM  
JET PRINTER - LOW VISCOSITY RESINS



Myerson's **breakthrough technology** redefines what's possible for the **3D construction of one print dentures**. The **precise fit and high strength of trusana1** reduces **chair time and recall appointments**.

Your Smile's Future Starts Here



**myerson**  
myersonstooth.com

Every Smile Has a Name – Print  
Your Patient's Name Directly in the  
Denture.



### CHAIRSIDE



*"I've had several long time analog denture patients comment on how comfortable their trusana1 dentures are." - Dr. Daniel Vogl, DDS*

### LAB SIDE



*"It's the best looking resin available with amazing qualities and life like appearance!" - Keith Dixon*

# BLAZING SPEED MEETS UNMATCHED ACCURACY



## FAST

Unprecedented print speeds — print 11 models in nine minutes.

## VERSATILE

Fully validated workflows for 15+ dental applications and complete freedom with Open Material Mode.

## INTUITIVE

An unmatched, seamless workflow with automated, mess-free printing anyone can learn in 15 minutes.

## RELIABLE

99% print success rate measured by an independent global leader in product testing.

## Powered by Formlabs Low-Force Display™

Form 4B and Form 4BL deliver unprecedented speed while maintaining superior dimensional accuracy, surface finish, and reliability across 15+ validated dental applications. With Formlabs' signature ease of use, Formlabs delivers an effortless and mess-free 3D printing workflow for even the busiest dental labs.

*"Consistent print quality with great surface accuracy makes this printer an excellent choice for dentistry. The surface comparison yielded some of the best results I have ever measured on a 3D printer."*

*Stephen Kreimer MDT, Kreimer Dental Lab*

### FORM 4B

Basic Package



Price: \$5,970

Item Number: 670-0077-FPD

### FORM 4BL

Basic Package



Price: \$12,900

Item Number: 670-0065-FPD

## High-Accuracy 3D Printing Materials for Dental Labs and Practices

Formlabs Dental Resins empower dental labs and practices to rapidly manufacture biocompatible surgical guides, splints, fixed patterns and models, clear aligner models, and full dentures.



A DIVISION OF HENRY SCHMIDT

formlabs | dental



# TRANSFORM YOUR DENTAL LAB

## With Blazing Speed and Precision

**Form**

**4B**

Validated for Pac-Dent  
Rodin and Dentsply Sirona  
Lucitone Digital Print



### Form 4B – Blazing Speed Meets Unmatched Accuracy

- Print up to 11 models in 9 minutes with unmatched accuracy
- 95%+ surface accuracy within 50 µm for a perfect fit every time
- Supports 15+ validated dental applications, including surgical guides, dentures, restorations, and orthodontics
- Easy, automated workflow with Build Platform Flex and spill-proof resin dispensing

**Form**

**4BL**



### Form 4BL – Speed That Scales

- Print up to 32 models in just 14 minutes for high-volume dental labs
- Large build volume (13.9 x 7.7 x 13.8 in) for full arches, aligner models, and appliances
- Repeatable accuracy across big batches – scale production without sacrificing quality
- Optimized for large-scale throughput with streamlined operation and maximum uptime

## Why Dental Labs Choose Formlabs

- Simplified post-processing with automated washing and curing solutions
- Wide range of dental resins, including biocompatible, sterilizable options
- Dedicated dental support and dental service plan

# ASIGA


## For high performance aesthetic dentures.



Retail: **\$384.99**

Retail: **\$274.99**

**Call Now for Volume Pricing**



## Max 2

### Asiga

3D printer that features Asiga's pioneering Smart Positioning System (SPS). Allows for precise layer thickness control of every layer to ensure optimal print accuracy and speed. Achieves the highest standard of reliability with the broadest wavelength capability of any imaging technology. Combines the latest generation of high-powered, solid-state, ultraviolet light sources with an integrated radiometer to ensure years of consistent printing performance.

385nm LED/Open Material

(251-0204-FPD) ..... Ea.

**Call**


## Ultra (50)

### Asiga

Asiga's first desktop 4K 3D printer. Offers a very simple user interface, large build sizes, SPS technology for production consistency, an integrated LED power sensor, fast print speeds, open material freedom, and an award-winning, user-centered design. Build plate size: 7½" x 4¼".

(146-2340-FPD) ..... Ea.

**Call**



## PRO 4K80 UV 3D Printer

### Asiga

An open material system that combines the latest 4K imaging technology with Asiga's well-proven Smart Positioning System (SPS) to produce a build volume 3X larger than its desktop 3D printers.

(140-7703-FPD) ..... Ea.

**Call**

Contains: 1 PRO 4K80 UV 3D printer, 1 kg. of DentaMODEL resin, 2 liters of Build Tray, 1 Asiga Flash post-curing device & 1 calibration tool kit (includes software).



**RODIN®**  
3D Resin

A Pac-Dent Brand



Full & Partial  
Denture Bases



All-On-X  
Denture Teeth



Universal Glaze  
Sealant & Repair

Universal Reline for  
all Dentures

## EXPERIENCE THE COMPLETE DENTURE WORKFLOW



### Denture Base 2.0

Denture Base Resin

*Impact-Resistant  
highly aesthetic dentures*

1.0 kg/Bottle ..... \$417.75  
Net Price ..... \$313.13

[146-2358-FPD] Original Pink  
[146-2361-FPD] Classic Pink  
[146-2630-FPD] Light Pink  
[146-2359-FPD] Meharry

### Titan

Ceramic Nanohybrid Resin

*Increased flexibility for long-lasting  
full arch restorations*

1.2 kg/Bottle ..... \$1285.51  
(NO CHARGE 5mL/Bottle, Valued \$99.95)

OM1 [225-0252-FPD]  
OM3 [225-0254-FPD]  
B1 [225-0255-FPD]  
A1 [225-0256-FPD]  
A2 [225-0257-FPD]  
A3 [225-0258-FPD]

### Glaze N<sub>2</sub>-Free

Universal Light-Curing Sealant

*Enables complete polymerization  
without the need of nitrogen curing*

5 mL/Bottle ..... \$99.95  
Net Price ..... \$66.63

15 mL/Bottle ..... \$249.95  
Net Price ..... \$166.66

[225-0312-FPD] Glaze N<sub>2</sub>-Free, 5mL  
[220-0494-FPD] Glaze N<sub>2</sub>-Free, 15mL

### SoftReline

Compatible with 3D, PMMA  
& Acrylic Dentures

*Universal reline with  
effortless application*

300 g/Bottle ..... \$199.95  
Net Price ..... \$133.33

[225-0312-FPD] Soft Reline Kit  
1x 50 mL VPS Reline Cartridge  
1x 10 mL VPS Adhesive Primer  
6x Wide Mouth Mixing Tips  
16x Applicator Swabs  
1x Tapered Cone Finishing Stone  
1x Brown Polishing Wheel

**Buy 1, Get 1  
Soft Reline Kit NO CHARGE**

**Buy 3, Get 1 NO CHARGE**  
(Mix and Match Different Shades)

**Buy 3, Get 1 NO CHARGE**  
(Of the same item)

**Buy 2, Get 1 NO CHARGE**  
(Of the same item)



\* For more information, contact your ZAHN Sales Representative, call (800) 496-9500 or visit [www.zahndental.com](http://www.zahndental.com). Offer valid from 9/28/2025 to 12/27/2025. To qualify and receive your no cost goods, please submit your qualified ZAHN invoice to [rodin@pac-dent.com](mailto:rodin@pac-dent.com). Cannot be combined with any other offer. No cash value. Equal or lesser value item(s) only. Promotions are subject to change or cancellation without notice. Please allow 6-8 weeks for delivery.

670 Endeavor Circle, Brea, CA 92821 · 909.839.0888 · [info@pac-dent.com](mailto:info@pac-dent.com) · [www.pac-dent.com](http://www.pac-dent.com) · [www.rodin-3d.com](http://www.rodin-3d.com)



## Pro 55 S Line 3D Printer

### SprintRay

Desktop 3D printer designed for digital dentistry to improve productivity in any practice. Ultrahigh resolution to meet the requirements of restorative applications.

- 55µm X-Y resolution
- 12cm x 12cm build area

3D

Printer .....Ea.

**Call**

Specify:

Pro55 S (251-0139-FPD) Pro 95 S (251-0138-FPD)

## Midas 3D Printer

### SprintRay

Patented digital press stereolithography (DPS) 3D printer that uses hydrodynamic principles to create an environment where the viscosity of the resin is no longer a barrier to speed or accuracy. This technology allows Midas to print up to 3 definitive crowns in 10 minutes, regardless of filler content. Using resin capsules simply scan, place, press, and remove print.

(146-7515-FPD) .....Ea.

**Call**



## Pro 2 3D Printer


### SprintRay

The latest precision 3D printer from SprintRay built to produce biocompatible dental appliances.

Features innovative Optical Panel Technology and a 385nm light engine, this printer achieves 35µm X-Y resolution for incredible accuracy. Redesigned for superior performance, it supports 15+ high-value workflows, including clear aligners, fixed restorations, night guards, dentures, and many more. It is the most user-friendly SprintRay printer ever, featuring a 12" touchscreen, automatic resin-tank and build-platform detection, and workflow error proofing. Advanced PrintOS software facilitates easy operation, delivering pertinent information to clinician's fingertips.

(146-7321-FPD) .....Ea.

**Call**



## cara® Print Cube

### Kulzer

Seamlessly transition from design to final print with this 3D printer. Its intuitive command center, the cara® Print Cockpit, ensures that each print job is prepared with precision to ensure the reliability of the prints. Features a user-friendly interface and advanced tools, including automatic model repair and support generation. CAM software allows clinician to make smaller changes to the 3D model directly in the CAM step without the need to return to the CAD software. Offers easy print job management of multiple print jobs across multiple printers.

- Build platform: 7.55" x 4.72" x 7.48" (192mm x 120mm x 190mm)
- Light source: 405nm LED

(665-0751-FPD) .....Ea.

**Call**





## Form Cure

### Formlabs

Light-Curing Unit for 3D Printing

Intuitive interface for easy temperature setting and optimal post-cure time.

- Advanced heating: Heat control is the key to precise, successful post-curing, carefully scaling temperatures up to 80°C in the cure chamber
- Balanced light: 13 LEDs use 405 nm of light to trigger the post-curing reaction, working with the heater to post-cure the parts
- Uniform exposure: Rotating turntable, heat flow, and dispersed lights provide uniform exposure during post-curing

(129-3052-FPD) .....Ea. **Call**



## Digital Cure— Small Capacity Unit

### Dentsply Sirona

Curing oven validated for use with all Lucitone Digital Print denture materials.

- Rapid curing of up to 5 arches in 60 minutes
- One-touch presets or custom mode
- Uniquely designed to provide higher temperatures and custom light wavelengths

(565-0112-FPD) .....Ea. **Call**



## Cure 1.5


### Asiga

UV Curing System

A professional UV curing system designed for precise postprocessing of 3D-printed dental and medical components. Key features include:

- Vacuum curing
- Infrared heating
- Inert gas connectivity

(251-0188-FPD) .....Ea. **Call**



## Form Cure L

### Formlabs

All stereolithography (SLA) 3D printed parts require postcuring: it helps them reach their maximum mechanical properties and achieve peak dimensional accuracy. Form Cure L uses two 500-watt heaters and high-powered UV lights alongside a reflective inner cavity and rotating turntable to provide a consistent, high-intensity cure. Whether it is for large prints or high throughput production of smaller parts, this unit can be used with all of Formlabs' SLA 3D printers. 120-Volt Unit

(141-1398-FPD) .....Ea. **Call**



## Digital Cure— Large Capacity Unit

### Dentsply Sirona

Uniquely designed to provide higher temperatures and custom light wavelengths that are required to achieve a validated cure of Lucitone Digital Print materials.

- Cures up to 9 arches in 90 minutes
- Low maintenance, long-life LEDs
- One-touch presets or custom mode
- Simplifies process with less labor

(142-2639-FPD) .....Ea. **Call**



## Otoflash

Optimized for quick curing of medical grade 3D Printed materials. Featuring two bottom mounted flash lamps operating at 10 flashes per second.

(220-0354-FPD) .....Ea. **Call**

## Otoflash-G171 Nitro

Light-curing unit for effective polymerization of light-cured resins under nitrogen gas. Chamber Size: 4.72" x 4.72" x 1.97"

Unit features 2 photo-flash lamps that produce 10 flashes/second with 200 watts of power in a spectral distribution from 280nm to 700nm.

(220-0352-FPD) .....Ea. **Call**





## ProCure 2

### SprintRay

#### Post-Curing System

Unbelievably fast. With its patented Light Motion Drive, this unit delivers outstanding light power and instantaneous heat to parts for same-appointment chairside care. By reducing cure times to 5 minutes or less, it can cut as much as 60 minutes from an in-office 3D printing workflow.

(251-0104-FPD) ..... Ea. **Call**



## cara® Print LEDcure

### Kulzer

Reliable postcuring unit for 3D-printed materials. Features long-life LEDs for both 385nm and 405nm resins, a transparent turntable so there is no need to flip objects, and a fan-assisted heater. Creates a stable thermal environment for better bond formation and excellent mechanical properties in terms of flexural strength, elongation at break, edge strength, and color stability. Also ensures restorations will not come out sticky.

Complete Unit  
(665-0689-FPD) ..... Ea. **Call**



## CURIE Light Cure Box

### Ackuretta

Equipped with a sophisticated, smart UV-light source management system that enables users to do more in less time and space. Maximizes energy efficiency and provides the ideal light strength to prevent any deformation.

- Overall dimensions: 7.9"L x 8.6"W x 6.2"H (20cm x 22cm x 16cm)
- Chamber dimensions: 5.5"L x 5.5"W x 3.1"H (14cm x 14cm x 8cm)
- Weight: 6.6 Lb. (3 Kg.)
- UV wavelength: 365nm, 385nm, 405nm
- Power input: 100–240 V AC, 56–60 Hz

(143-4472-FPD) ..... Ea. **Call**



## CURIE Plus Light Cure Box

### Ackuretta

Can cure up to 8 full arches at once. Uniform light distribution and 365nm, 385nm, and 405nm UV wavelength coverage ensure each print layer is thoroughly cured in a single process, eliminating the need for manual flipping. Compatible with a range of nitrogen sources.

- Overall dimensions: 14.1"L x 10.2"W x 9.2"H (36cm x 26cm x 23.4cm)
- Chamber dimensions: 7.8"L x 5.9"W x 2.4"H (20cm x 15cm x 6cm)
- Weight: 26 Lb. (12.2 Kg.)
- UV wavelength: 365nm, 385nm, 405nm
- LED lights: 65
- Power input: 100–240 V AC, 56–60 Hz

(145-5840-FPD) ..... Ea. **Call**



## Ductless Fume Hood 3D

### Quatro Air Technologies

Wall-mounted, ductless fume hood (DFH) for eliminating 3D-printer and alcohol-cleaning odors. Offers multistage filtration with a HEPA filter.

- Stops odors, harmful particulate matter, and vapors from migrating throughout the workplace
- Collects, neutralizes, and deodorizes local odors from 3D printers and alcohol-based cleaning products
- True deep-bed HEPA filter ensures 99.97% efficiency when capturing airborne particles
- Print below the hood
- Designed for quiet operation
- Comes with a wall-mounting bracket
- Optional: Hand Switch for 6-foot remote Start-Stop (#AA011), DFH Counter Stand (#AA038), 6" diameter Exhaust Kit (#AA088), and Silencer for more noise suppression (#H250)

(270-0665-FPD) ..... Ea. **Call**




## MonoPure 3D Printing and Cleaning Station Kit

### Quatro Air Technologies

Fume extractor kit for desktop 3D printers and alcohol cleaning stations. Equipped with odor and HEPA filters. Includes a MonoPure multistage fume extractor with a remote foot pedal, desktop containment cabinet (30"H x 18"W x 16"D) with hinged Plexiglas® doors, 24" flex-arm with 8" diam. Plexiglas® nozzle and table bracket, and two 8-ft. lengths of 3" diameter hose. Cabinet has a 3" exhaust port that can be installed on either the side or back panel for hose connection to the MonoPure fume extractor.

(270-1121-FPD) ..... Ea. **Call**



## Containment Air Box 3D Printer

### Quatro Air Technologies

New larger containment box for desktop 3D printers: 31"H x 20½"W x 20"D. Powder-coated, heavy-duty steel walls with hinged front doors swing up and away for easy access to the printer. 3" exhaust port can be installed on either the side or back panel for hose connection to Quatro fume extractor.

(270-1123-FPD) ..... Ea. **Call**




## DCS 3D Kit

### Quatro Air Technologies

All-in-one fume extractor and containment cabinet designed for most desktop 3D printers. Has a specially designed odor filter capable of handling harmful fumes associated with 3D printer resins and postproduction alcohol cleaning stations. Tall containment cabinet dimensions: 31½"H x 20½"W x 20"D. Saves heating and air conditioning costs by putting clean air back into the room.

(270-1126-FPD) ..... Ea. **Call**



## MonoPure 3D Printer Kit

### Quatro Air Technologies

Fume extractor kit for desktop 3D printers and alcohol cleaning stations. Equipped with odor and HEPA filters. Includes a MonoPure multistage fume extractor with a remote foot pedal, new larger desktop containment cabinet (31"H x 20½"W x 20"D) with hinged Plexiglas® doors, 24" flex-arm with 8" diam. Plexiglas® nozzle and table bracket, and two 8-ft. lengths of 3" diameter hose. Cabinet has a 3" exhaust port that can be installed on either the side or back panel for hose connection to the MonoPure fume extractor.

(270-1122-FPD) ..... Ea. **Call**



## Form Wash (2nd Generation)

### Formlabs

Automatically cleans printed parts thoroughly and efficiently with IPA or alternative solvents. Allows parts to be kept directly on the build platform as they are transitioned from the printer to the unit, or the parts can be removed and placed into the wash basket. Parts are precisely and thoroughly agitated in solvent with an impeller, getting every nook and cranny perfectly clean, automatically for the appropriate time. When washing completes, unit automatically raises parts out of the IPA for air drying.

(670-0010-FPD) .....Ea. **Call**



## Form Wash L

### Formlabs

All stereolithography (SLA) 3D printed parts require washing to remove residual resin. Form Wash L automates this process, leaving clinician with smooth, clean, high-quality parts every time. Once washing completes, parts are lifted out of the solvent to air dry, avoiding over-soaked, warped prints. Embedded solvent monitor helps maintain peak cleaning performance by notifying clinician when it is time to replace the solvent. Can be used with all Formlabs SLA 3D printers. Wash one Form 3L build platform or two Form 3 build platforms with the included build-platform adapter.

(141-1397-FPD) .....Ea. **Call**



## ProWash S

### SprintRay

The next generation of automatic post-processing equipment for dental 3D printing. It provides a 2-stage wash and thorough dry for printed parts in just 9 minutes.

- Patented mechanical jetting with all-new power-dry capability
- Cloud-connected for automatic job handoff from print to cure
- Dedicated drying hardware and sensors for better quality of life

#SRI-0504012

(145-3466-FPD) .....Ea. **Call**



## cara® Print CleanUp

### Kulzer

Save time and effort with this reliable 3D printer: It efficiently cleans 3D-printed dental models with ease. Designed for versatility, it works seamlessly with multiple cleaning buckets and offers hands-free operation to clean parts directly on the build platform. Audible and light alerts ensure clinician stays in control without constant monitoring. Using isopropyl alcohol (IPA), it delivers thorough, precise cleaning with minimal manual intervention, making it the perfect addition to any workflow.

(665-0752-FPD) .....Ea. **Call**



## CLEANI with Dual Tanks

### Ackuretta

Dual-tank washer for dental 3D printing. Save time and 50% alcohol while removing 90% of resin residue in a gentle, 2-stage process. Clean platforms and prints together, minimizing damage and boosting workflow efficiency.

(146-1076-FPD) .....Ea. **Call**

# 2025 TAX CHANGE GUIDE

## WHAT DENTAL LABS NEED TO KNOW & PLAN FOR NOW

### Section 179 Deduction Calculator



## Turn Your Equipment Purchase Into a Tax Win



### PURCHASE

equipment before December 31st\*



### DEDUCT

Instantly see how much you can deduct – and how much you'll save



### SAVE

on your taxes when equipment is purchased and placed into service by the end of 2025\*

Use this simple tool to estimate your Section 179 tax savings on new equipment purchased and placed into service in 2025

**SEE HOW  
MUCH YOU  
CAN SAVE!**

Equipment Price: 150,000.00 Required field

Estimated Tax Rate/Bracket: .10 Required field

Section 179 Deduction:  
(\$2.5M in 2025) 150,000.00

Tax Savings: 15,000.00

Net Cost of Equipment: 135,000.00

**Amount Paid to the IRS if You  
DO NOT Acquire Equipment!:** 15,000.00

**TALK TO US  
ABOUT BONUS  
DEPRECIATION.**



**EXPLORE EQUIPMENT THAT QUALIFIES FOR SECTION 179**



**DISCLAIMER** Neither Henry Schein, Inc. nor Henry Schein Financial Services provides tax advice. Please consult with a qualified professional tax advisor to determine your eligibility for a Section 179 tax deduction. Must purchase and place into service by December 31, 2025. ©2025 Henry Schein, Inc. Neither Henry Schein, Inc. nor Henry Schein Financial Services provides financial advice. Please consult your financial advisor. Neither Henry Schein, Inc. nor Henry Schein Financial Services is a bank and neither represents itself as such, nor conducts banking activities. Henry Schein Financial Services may receive a marketing fee from the vendor for products/services purchased.





**Sierra Dental Tool**  
POWERED BY ROBBJACK

# TIME TO REPLACE THOSE OLD TOOLS ALREADY?

**Milling and  
Grinding Tools for:**  
**Amann Girrbach**  
**Aidite**  
**DentMill**  
**Imes-Icore**  
**Ivoclar PM7**  
**Sirona**  
**Roland**  
**VHF**  
**Wieland**  
**XTCera**



**SPECIAL OFFER!**

**15% off**

# ALL ROLAND TOOLS!

Mention code **FPD** when ordering!  
Offer expires December 27, 2025

Contact your Henry Schein Sales Representative for additional information




**Sierra Dental Tool**  
POWERED BY ROBBJACK

# Maximize Your In-house Workflow!



**AK-1 Material Adapter**



**DGSHAPE CAM Software**



**DWX-43W  
Wet Milling Solution**

**Get Everything You Need for ONLY \$33,995**

**SAVE \$1,765.40!**

✓ Complete Package ✓ Premium Features Included ✓ Unbeatable Value

**Everything you need is included, so you can start milling  
Titanium Custom Abutments and start growing your business.**

#### **Here's what's included:**

DWX-43W, DGSHAPE CAM Software, (1) Abutment Material Adapter for MEDENTiKA® or DESS®, or GeoMEDI®, (2) ZTR-150D, (2) ZTB-100D, (2) ZTB-50D, 6-Pin Adapter (ZV-42W), & Bonus Tooling: (3) ZGB2-125D, (1) ZGB2-125DS, (1) ZGB2-75D, (3) ZGB2-50D, (1) ZGB2-25D, (3) ZGB-125D, (3) ZGB-50D, and (1) ZGB-25D.

**For more information, contact your Zahn Sales  
representative today!**



Offer valid through 12/27/25. Only available in the USA. Terms and conditions apply. Pricing reflects US MSRP and promo is valid for SKU# DWX-43W-ME-K, DWX-43W-DE-K, and DWX-43W-GE-K: which includes DWX-43W, DGSHAPE CAM Software, (1) Abutment Material Adapter for MEDENTiKA® or DESS®, or GeoMEDI®, (2) ZTR-150D, (2) ZTB-100D, (2) ZTB-50D, 6-Pin Adapter (ZV-42W), & Bonus Tooling: (3) ZGB2-125D, (1) ZGB2-125DS, (1) ZGB2-75D, (3) ZGB2-50D, (1) ZGB2-25D, (3) ZGB-125D, (3) ZGB-50D, and (1) ZGB-25D. Items must be purchased together on one invoice to qualify and cannot be combined with any other offer or prior purchase. Roland DGA reserves the right to modify or discontinue the terms of this promotion at any time, without prior notice.



**DWX-43W****DG Shape/Roland****Wet Milling Solution**


Brings faster milling times for traditional ceramics and custom titanium abutments, while reducing operating costs and increasing ROI.

- Faster milling times
  - No more costly filters
  - Reduced routine maintenance
  - Expanded milling capabilities built-in
  - DGS SHAPE CAM software and Bbnus tooling included
- (595-0607-FPD).....Ea.

**Call****DWX-53D****DG Shape/Roland****Dry Milling Solution**

- Redesigned, more robust spindle
- Faster PMMA and zirconia milling speeds
- Advanced 90°-like milling capabilities
- Built-in, fully integrated webcam allows for real-time monitoring either on-site or remotely


(595-0591-FPD) .....Ea.

**Call****DWX-53DC****DG Shape/Roland****5-Axis Dry Dental Mill with Automatic Disc Changer**

- Redesigned, more robust spindle for 3 times the gripping power with PMMA applications
- Built-in integrated webcam, automatic extractor filter cleaning, and thermohydrometer
- Accelerometer sensors
- Compatible materials: PMMA, zirconia, wax, composite resin, PEEK, gypsum, glass-fiber-reinforced resin, and presintered chromium cobalt

(595-0511-FPD).....Ea.

**Call**



## **ceramill® matron**

Processes Plastic, Metals, Ceramics, Hybrids, Waxes, Composites, Glass Ceramics, and Zirconia

Vmax-Arm (45°) for fastest processing of all dental materials

21.5" curved touch display with guided workflows

## Perfectly paired with



### **zolid bion**

More than just a 3D zirconia

All types of indications  
(single unit to full-arch)  
45-minute\* sintering cycle (crowns)  
1,100 MPa strength throughout



### **zolid fx multilayer**

Beauty inclusive

High esthetic crowns  
up to 3 unit bridge  
2-hour\* sintering cycle  
700 +/- 150 MPa



### **zolid gen-x multilayer**

The esthetic all-rounder

Wide range of indications  
25-minute\* sintering cycle  
for single units  
1000 +/- 150 MPa

With the Ceramill Zolid line, we offer you the right material for every zirconia-based indication. [bit.ly/ZolidZi](http://bit.ly/ZolidZi)



Have Questions? Contact Us.



**Steven Klein**  
**Zahn Product Specialist**  
(800) 496-9500  
[Steven.Klein@henryschein.com](mailto:Steven.Klein@henryschein.com)

\*Can sinter in longer cycles according to indication



### Ceramill Matron Amann Girrbach

Fifth-generation milling machine with superior properties compared to machines of older generations. Provides stronger, faster, and smarter technology than conventional milling machines, regardless of clinician's level of skill. Offers outstanding precision and excellent results at high speed across a wide range of materials. High-performance metal machine with 3mm or 6mm cutter shank, 45° A/B axes, and 5-axis+ milling power. This unit provides maximum operating convenience, maximum versatility, and valuable support. SnapMag magnetic tool tray sold separately.

- Vmax arm for maximum strength and precision
- Retractable performance blank changer
- RFID-supported tool and blank management
- Guided workflows on a 21.50" touch display
- Dimensions: 37.01"W x 33.46"H x 27.36"D

(146-7196-FPD).....Ea.

**Call**



### Ceramill® Motion 3 Amann Girrbach

- Fully AG digital: AG.Live connection, display, and RFID technology
- Lab intelligence to track your business
- Grinding, carving (patented), milling, and detailing capabilities for maximum versatility and the highest quality of restorations
- Self-service operations, automated tools, holder management, and auto-calibration
- Fully automated full-arch restorations with only manual finishing required (C-clamp and strategies)
- 100,000-rpm unit

Dry Package  
(220-0357-FPD).....Ea.

**Call**

Hybrid Package  
(220-0358-FPD).....Ea.

**Call**

Contains: 1 of each: Motion 3 mill, Coolstream & Coolstream cart/trolley.



### Ceramill Motion 3 Air Amann Girrbach

Elevates hybrid milling technology with its integrated air supply. Experience enhanced process reliability, efficiency, and user-friendliness without the need for additional infrastructure.

- Plug and play convenience with integrated compressor
- Interactive 10" display with guided workflows
- Error-prevention with RFID technology
- Time-saving cleaning and calibration
- Versatility through material and indication range

(288-4406-FPD).....Ea.

**Call**



### Ceramill® Matik Amann Girrbach

Fully integrated fabrication unit that offers the laboratory the possibility to design its workflow independently and flexibly. Dramatically reduces the effort and complexity of material and tool management. Features automatic maintenance, so cleaning and changing the operating modes save time and help the technician to focus on the essentials: the dental work itself.  
(288-3759-FPD).....Ea.

**Call**



### XTCERA X-mill 500 Plus

**Shenzhen Xiangton Co., Ltd.**

High-performance milling machine for all-around processing and production of dentures and other materials. Outstanding performance with a highly sealed processing chamber for better dust proofing (accuracy as high as 2 microns). Intelligent monitoring with stronger situational awareness, warning of work status and malfunctions, and convenient remote-control operation. Performs a variety of tasks, including fast cutting of hard materials. Safety ensured with automatic start-stop when opening/closing doors to reduce safety hazards. Outstanding appearance with a simple, sleek design and a glossy, easy-to-clean surface.

(613-0069-FPD) ..... Ea.

**Call**



### XTCERA X-Mill 580 DC

**Shenzhen Xiangton Co., Ltd.**

A 5-axis dry-mill milling machine for zirconia, wax, and PMMA materials. Water cooled, 1.8kW spindle and intelligent automatic disc changing enable continuous (24/7) processing. Large-capacity material storage can hold up to 12 discs with an open-material design. Smart monitoring automatically detects disc lifespan and wear. Remote operation for stronger situational awareness, real-time monitoring of work status, and the ability to control the machine remotely.

(613-0086-FPD) ..... Ea.

**Call**



### X-Mill 300

**XTCERA**

Wet/Dry Milling Solution

- High-accuracy, harmonic A-axis drive of 360° rotation
- 1800-watt Industrial spindle allows fast milling of hard materials like acrylic and aluminum
- Multifunction holder, 6-slot automatic tool changer, and harmonic A-axis drive ensure precision milling

(613-0087-FPD) ..... Ea.

**Call**



## iVAC 2.0 CAD/CAM Dust Collector

### Quatro Air Technologies

The new series of fully automatic, self-cleaning, bagless, dust collectors for bench applications. On-line filter-cleaning function to automatically clean the main filter without operator assistance. HEPA filter provides 99.97% efficiency on submicron particles, including zirconia, alloy, sand, ceramic, acrylic, and die stone. AutoFlow feature allows user to set to the lowest acceptable speed, and the motors automatically speed up to overcome pressure increase. Separate access doors allow user to remove debris without having to access filter section. Silent operation, making it one of the quietest compact units on the market (55 dbA at high speeds). 120 volts. 1-year warranty.

#### iVAC 2.C TWIN CAD/CAM Dust Collector

For heavy-duty milling with any mill. Debris storage of up to 6 gallons. 2 long-lasting INFINITY motors. Dimensions: 26"H x 16"W x 20"L  
Unit with Double INFINITY Motors-1-0052-0001

(131-2003-FPD) ..... Ea.

Contains: 1 of each: unit, auto-start cable & 6-ft. hose.  
(Specify mill model when ordering.)

**Call**

## iVAC 2.C Single Infinity Dust Collector

### Quatro Air Technologies

Fully automatic, self-cleaning, bagless dust collector for CAD/CAM milling. For medium- to heavy-duty milling with any mill. Debris storage of up to 4 gallons. Long-lasting Infinity motor. Dimensions: 25"H x 14"W x 15"L

#### iVAC 2.C Single CAD/CAM Dust Collector

For medium-heavy milling with any mill. Debris storage of up to 4 gallons. Dimensions: 25"H x 14"W x 15"L  
Unit with INFINITY Motor-1-0051-0001

(270-0427-FPD) ..... Ea.

Contains: 1 of each: unit, auto-start cable & 6-ft. hose. (Specify mill model when ordering.)

**Call**

## iVAC 2.C Twin Brushless Dust Collector

### Quatro Air Technologies

CAD/CAM dust collector for the Roland milling machine. Fully automatic online filter-cleaning function: automatically cleans the main filter without operator assistance. HEPA filter: 99.97% efficiency on submicron particles (zirconium dust). Autoflow feature sets unit to lowest acceptable speed; motors automatically speed up to overcome pressure increase.

(270-0560-FPD) ..... Ea.

Contains: 1 of each: unit, auto-start cable for Roland milling machines, 6-ft. hose & remote foot pedal.

**Call**

## iVAC 2.C Twin Brushless Infinity Dust Collector

### Quatro Air Technologies

Fully automatic, self-cleaning, bagless dust collector for CAD/CAM milling. For heavy-duty milling with any mill. Debris storage of up to 6 gallons. Maintenance-free brushless motor. Dimensions: 26"H x 16"W x 20"L

Unit with Brushless Motor-1-0052-4002

(270-0470-FPD) ..... Ea.

Contains: 1 of each: unit, auto-start cable & 6-ft. hose (specify mill model when ordering.)

**Call**

## SILENT compactCAM Extractor

### Renfert USA

Compact, bag-free extraction unit with automatic filter cleaning and TC collector motor. For use with compact, light-workload CAM milling machines. Ultraquiet (56.9 dB) under maximum power. Bidirectional communication with the integrated CAM unit (interface cable sold separately). 3-year warranty included. 120 volts. Volume flow (maximum): 2,500 liters/minute, 1.47 ft.<sup>3</sup>/second, 88 cfm  
Dust drawer capacity: 0.69 gallons  
Filter: Self-cleaning, bagless, 99.9% efficiency  
Noise level (maximum): 56.9 dB(A)  
Motor warranty: 1,000 operating hours  
Motor rating: 800 hours  
Dimensions: 17.3"H x 9.7"W x 19.7"D

(228-4785-FPD) ..... Ea.

Contains: 1-2-meter suction hose & 2 end bushings.

**Call**



## SILENT PowerCAM EC Extractor

### Renfert USA

Extraction unit for CAM systems. Powerful, bag-free extraction unit with automatic filter cleaning and EC (electronically commutated) motor technology for heavily used CAM units.

Advantages:

- Bidirectional communication between extraction unit and CAM unit, thanks to the PLC (programmable logic controller) interface
- Low operating noise of maximum 54.3 dB(A)
- Safe production process through low-maintenance EC motor technology (5000 operating hours)
- High volume flow of approximately 4000 liters per minute
- No filter bag needed

(228-5803-FPD) ..... Ea.

Contains: 1 quick-start guide, 1 quick-reference card, 3 disposal bags (1 piece inserted in the dust drawer) & 1-2-meter suction hose (includes 2 mufflers).

**Call**



## SILENT powerCAM TC Extractor

### Renfert USA

Extraction Unit for CAM Systems

Powerful and inexpensive CAM extraction unit with collector motor, especially suitable for medium to large CAM systems. Consistently full suction power due to automatic filter cleaning. Reliable production process through innovative motor technology (3 times longer service life in comparison to conventional collector motors).

- High suction capacity up to 3,980 lpm
- Volume flow (maximum): 3,980 lpm
- Dust-drawer capacity: 1.85 gallons
- Filter: self-cleaning, bagless, 99.9% efficiency
- Noise level (maximum): 54.3 dB(A)
- Motor warranty: 1,000 operating hours

(288-6237-FPD) ..... Ea.

**Call**





**NO CHARGE Starter Tool Kit PLUS 1 Year Extended Warranty  
VALUED AT \$10,740!**

### E4 Mill

#### VHF, Inc.

4-axis wet milling machine with an optional dry machining set with the following technical data:

- 4 simultaneously working axes
- 3 linear axes (X/Y/Z) with precise ball screw drives
- Rotary axis A with +190° to -10° rotation range
- Synchronous spindle SFE 300E, maximum power output 800 W, continuous power 300 W; rotational speed up to 60,000 rpm
- Automatic tool changer with 6 tool positions and 1 insert for AIRTOOLS
- No compressed air required due to AIRTOOLS and the electrical release mechanism for the collet chuck
- Ethernet port, air extraction interface
- Switching power supply 100–240 V AC
- Footprint: 360mm W x 370mm D x 490mm H
- Weight: 28 kg
- Includes DENTALCAM software with DIRECTMILL function #30645C

(700-6430-FPD) ..... Ea. **Call**



### E5 Mill

#### VHF, Inc.

A dry milling machine for discs that masters even complex indications with the highest precision thanks to its 5 axes. Despite its versatility, the E5 can be operated with ease: it does not require any compressed air and is therefore very cost-effective and particularly sustainable. Equally innovative is its low machine weight. Both aspects reduce the requirements for the installation site, which is particularly helpful in a practice environment. Also features easy serviceability and handling.

- Plug and mill: Unpack, connect, and start milling
- Maximum freedom in terms of indication, scanner, and material
- Airtool eliminates the need for compressed air #30633C

(700-6427-FPD) ..... Ea. **Call**



### R5 Mill

#### VHF, Inc.

A highly automated, 5-axis milling and grinding machine that brings groundbreaking innovations to the lab, along with a 10-position blank changer and the capability for wet and dry processing. On a minimal footprint, it combines maximum precision and stability, and it is capable of handling nearly any material. DirectDisc technology ensures fastest blank loading and clamping without tools. DirectClean technology allows the unit to not only effortlessly switch between wet and dry processing, but also perform intelligent self-drying and self-cleaning programs.

(700-5371-FPD) ..... Ea. **Call**



### S5 Mill

#### VHF, Inc.

An exceptionally versatile and affordable milling and grinding machine for dusk-till-dawn processing. It has 5 simultaneously operating axes and a blank changer for 8 blanks, plus it supports both dry and wet operation. In addition to discs, the machine can process up to 24 blocks or 48 abutment blanks. With the optional wet-grinding module, users can process all common types of glass ceramics and titanium. Three integrated ionizers considerably reduce the cleaning effort in the milling chamber. Optional QuickFrame allows easy loading of blanks without any tools.

(700-5359-FPD) ..... Ea. **Call**



### N4+

#### VHF, Inc.

A wet-process machine for milling and grinding glass ceramic, composite, and zirconia blocks as well as titanium abutments. Combines high precision and fast processing in an extremely compact housing.

- 800-watt spindle capable of speeds of up to 80,000 rpm
- 4-axis milling machine with an integrated camera for easy remote support
- Convenient, removable cooling liquid tank
- Automatic changer for 8 tools
- Webcam for working chamber for remote monitoring and service
- Ethernet interface for stable connection
- Very easy operation via DENTALCAM software with DIRECTMILL Technology (included, without license fees)

(700-6465-FPD) ..... Ea. **Call**



### K5+ Milling Unit

#### VHF, Inc.

Equipped with the innovative DirectDisc technology for securing 98mm discs; separate block and prefab abutment holders are also available. A tool-free blank clamping system enables the machine to be loaded with materials quickly and easily. Unit covers the broadest range of indications, due to the  $\pm 35^\circ$  rotation angle in the fifth axis, and covers up to 40mm blanks.

(700-4288-FPD) ..... Ea. **Call**





# AUSTROMAT™ 674i™ G2

TAKE THE BEST YOU CAN GET.

## SiC TECHNOLOGY



### HOT

SiC technology for **1600 °C MATERIALS**

### CLEAN


SiC technology for **COLOR FIDELITY**

### FAST

SiC technology for **HIGHEST SPEED**

### GREEN

SiC technology for **MAXIMUM EFFICIENCY**



(289-0125-FPD)

**IoT**

### CLEVER

IoT for **EASY PROGRAM STREAMING**

### SMART

IoT for **DRIVING COMPREHENSIVE SERVICE**



Valid from June 29 to December 27, 2025

# Shenpaz

One Step Ahead


## FLEXIBILITY BEYOND LIMITS MAX SPEED, MAX CAPACITY



Sintra Pro Max




P/N: 727-0630 MJR



**Sinter Smarter:** Compatibility, total in-lab production versatility  
**Less Waiting, More Creating:** Efficient milling & sintering in sync

| To place an order, contact your Friendly Henry Schein Rep




## Sintra PRO 90 ZR Furnace

### ShenPaz

- High sintering temperature of 3002°F (1650°C)
- 3 MoSi2 heating elements
- New tray with no beads required
- Maximum capacity of 90 restorations: able to stack 3 trays
- 7" intuitive touch screen
- Unique degassing feature purifies muffle from impurities
- Unlimited programming possibilities
- WiFi connectivity pairing feature to all devices
- 115 volts

Complete Unit

(727-0592-FPD) ..... Ea. **Call**  
Contains: 1 of each: furnace, power cable, sintering table, sintering crucible, sintering crucible top, tweezers, user manual & lift key.



## Sintra PRO 120 ZRF Furnace

### ShenPaz

- 4 MoSi2 heating elements
- Large chamber for FULL large arches
- 7" intuitive touch screen
- Maximum capacity of 120 restorations (3 stackable trays optional)
- Modular sintering table permits adjustment of table height and chamber size
- WiFi connectivity allows pairing with other devices: FULL operational and monitoring control
- Fast sintering within 50 minutes
- Subject to zirconium dioxide material
- 230 volts

Complete Unit

(727-0590-FPD) ..... Ea. **Call**  
Contains: 1 of each: furnace, power cable, firing table, sintering crucible, sintering crucible top, tweezers, USB, user manual & lift key.



## Sintra CS Furnace

### ShenPaz

Fast sintering furnace. Multiprocess: sintering, crystallization, and glazing. Extra-wide table, 70mm capable of sinter 7 individual crowns or a small bridge. Programmable open system can support a variety of existing zirconia materials and can be adapted to support future ones. Rapid cooling in the final stage allows quick handling of the crown by hand. Large touch screen.

230-V, 60-Hz Unit

(727-0603-FPD) ..... Ea. **Call**



## SINTRA PRO Max Sintering Furnace

### ShenPaz

Elevates the dental laboratory's productivity and precision. Offers unparalleled speed, capacity, and compatibility. Achieves rapid results with the ability to sinter up to 10 crowns or small bridges in record time. Handles high-volume workloads effortlessly with the ability to sinter up to 50 crowns and full bridges on two trays. Uses up to 70% less electricity consumption per cycle. Ideal unit for labs looking to maximize output without compromising on quality.

(727-0630-FPD) ..... Ea. **Call**



**High Quality Zirconia, 16mm  
as low as \$79.99**

**Volume Discounts Available  
Call for Details!**

# Discover DD cubeY<sup>®</sup> HL

the Hybrid-layered Zirconia  
for Superior Anterior Esthetics

16 VITA<sup>®</sup> shades  
3 VITA<sup>®</sup> bleach shades  
16mm, H20mm, H25mm size discs

Multi-indicative stability\*  
Incisal area: > 800 MPa  
Body area: 1050 (±150) MPa

\*A dual-gradient zirconia with 3Y in the body for strength and 5Y in the incisal for translucency—designed to mimic the natural balance of durability and esthetics.



**PRODUCT**  
DD cubeY<sup>®</sup> HL



**cube** teeth for tomorrow



# SMARTER EQUIPMENT SUPPORT STARTS HERE

**Your** Equipment. **Your** Support. On **Your** Terms.

## What is ServiceFirst?

ServiceFirst is a monthly support subscription for CAD/CAM equipment purchased through Zahn Dental. With a simple QR scan, your lab can instantly access equipment-specific training, troubleshooting, manuals, chat support, and maintenance tools—anytime, anywhere.

### Subscription

- › For just **\$49.99/month with a 12-month agreement**, get peace of mind and reduce downtime for every CAD/CAM product that Zahn currently supports you with
- › One flat fee. Unlimited devices  
Unlimited access.



### Optimize Uptime, Reduce Downtime

- › Assign tasks by team or equipment to take preventative maintenance steps
- › Track maintenance, open tickets, and prevent issues before they happen



### All Your Support in One Place

- › Centralized access to manuals, training videos, IFUs, and guides
- › Organized by machine, accessible to anyone in the lab



### More Equipment = More Value

- › One subscription covers all CAD/CAM equipment purchased from Zahn
- › Perfect for existing digital workflows or labs adding to their production process



### Scan. Submit. Solve.

- › Every machine gets a custom QR sticker
- › Scan to access a personalized portal for that device



### Always-On Support

- › Skip the phone queue—chat with an agent via chat messaging
- › Upload images, videos, or documents to resolve issues fast



## How It Works

- › Subscribe to ServiceFirst
- › Receive QR stickers for each machine
- › Scan to access your custom support portal
- › Start solving support needs instantly from your smartphone

**SCAN QR CODE  
TO GET STARTED  
TODAY**





# Zirlux® Transitions

## Introducing a Transitional Layered Zirconia that mimics nature and lifelike esthetics without compromising strength

Optimal Translucency Technology  
provides higher chroma in the cervical  
area and decreasing chroma toward  
the incisal area

Perfect balance between  
Strength and Translucency

- 1200 MPa
- 48% Translucency



Available in all 16 classic VITA\*  
shades plus OM2 bleach shade.

- ✓ 98.5 x 14mm
- ✓ 98.5 x 18mm
- ✓ 98.5 x 22mm

One premium zirconia  
for anterior and posterior  
restorations

\*VITA is a registered trademark of Vita Zahnfabrik.

Made in the USA

# MUA

## Why choose DESS® Multi-Units?

A key advantage of the DESS® Multi-Unit system is its standardised prosthetic platform across multiple implant compatibilities. Combined with a consistent treatment protocol, this leads to greater predictability and higher success rates, even in complex cases involving different implant systems.

- Full-arch rehabilitation *with just 4–6 implants*
- Shorter treatment times *with efficient workflows*
- Graftless procedures *in many cases*
- Optimised insertion axis *ensuring a precise prosthesis–implant connection*
- Multiple gingival heights *to adapt seamlessly to different implant depths*

**Buy 2 of the same  
Get 1 of the same at no charge!**



**Ships with Order**

**SAVE OVER 20%**



## Gibraltar® Pure Blast Aluminum Oxide

**HENRY SCHEIN®**

- 50µm Aluminum Oxide
- Machine sized and computer scanned to provide a uniform particle size
- Packaged to prevent moisture absorption
- Tan is designed for non-precious alloys only
- White is for precious alloys, can also be used for non-precious

Tan		
(900-4084-FPD)	25 Lb.	<b>73.99</b>
White		
(900-4085-FPD)	25 Lb.	<b>64.99</b>
Tan		
(900-4086-FPD)	47.5 Lb.	<b>133.99</b>
White		
(900-4087-FPD)	47.5 Lb.	<b>119.99</b>

Buy Basic eco,  
Get 1 - Cobra White 25 Micron 5Kg\*

**SAVE OVER \$275**



## Basic eco Fine Sandblasting Unit

**Renfert USA**

Compact microblaster with 2 tanks. Provides for maximum field of vision, thanks to large viewing panel and LED lighting. Uniform flow of abrasive guarantees a consistent surface conditioning with very low consumption. Allows for quick switching between abrasives without interrupting work. 3-year warranty included. 120 volts.

Basic eco 25 µm/70 µm, 70 µm/250 µm  
(228-2682-FPD) ..... Ea. **899.99**  
\*Ships with order.

Buy Basic Classic,  
Get 1 - Cobra White  
25 Micron 5Kg\*

**SAVE OVER \$325**



## Basic Classic

**Renfert USA**

- 25µm/70µm, 70µm/250µm
- 120 volts
- 3-year warranty

Basic Classic

Fine sandblasting unit with 1-2 tanks. Low-priced entry to professional microblasting technology. Features precise blasting with special mixing chamber technology and tank selector switch for quick switching between abrasives without interrupting work.

(228-2065-FPD) ..... Ea. **1,219.99**  
\*Ships with order.

Buy 2, Rolloblast  
100 micron 5Kg,  
Get 1 of the same\*

**SAVE OVER \$50**



## Rolloblast

**Renfert USA**

Nonabrasive glass beads for smoothing and condensing of surfaces. Guarantees a silky, matte surface with no light reflection. For the finest surface treatment of precious and nonprecious alloys. Does not cause silicosis.

50µm  
(228-9533-FPD) ..... 5 Kg. **45.99**  
100µm  
(228-2461-FPD) ..... 5 Kg. **45.99**

\*Ships with order.

**SAVE \$100**



## Cut & Clean Cabinet XL

**Vaniman Manufacturing**

Dust-cabinet solution for postprocessing CAD/CAM dental mill pucks or other applications where dust needs to be contained and the use of an internal compressed air nozzle is needed. Keeps work area clean and healthy while containing dust created in the workspace. Complete Unit

(787-0460-FPD) ..... Ea. **469.99**



# Healthy Labs Need Clean Air

When you work in a dental lab,  
capturing dust right at the source  
isn't just smart, it's essential.



## THE PROBLEM:



**Health risks**

Trimming, grinding, and polishing are everyday tasks in dental labs and they generate a dangerous mix of airborne contaminants including fine silica and zirconia dust, heavy metal particles, resins, acrylics, and ceramics that pose serious respiratory and long-term health risks to technicians. Without proper capture and containment, these particles can migrate throughout the lab, compromising air quality and worker safety.



**JetStream  
Basic Infinity  
#270-0606-FPD**

- Single-speed
- 6-foot hose

**Price: \$ 659.99**  
**SAVE OVER \$100**

**JetStream  
HealthSmart Infinity  
#270-0624-FPD**

- Variable-speed
- "T" connector
- Two 6-foot hoses

**Price: \$ 749.99**  
**SAVE OVER \$125**

**JetStream  
ASV Infinity  
#270-0607-FPD**

- Intelligent controller displays status and alerts
- ASV cable included
- "Bag-full" indicator
- Variable-speed
- "T" connector
- Two 6-foot hoses

For use with Quatro's  
ASV-Ready Accessories

**Price: \$ 799.99**  
**SAVE OVER \$125**



Clean air is a smart investment  
in health and efficiency.  
We're here to help you find the  
solution that works best for you!

## OUR SOLUTION: SMART SOURCE CAPTURE



The JetStream Series Compact Dust Collectors are purpose-built to tackle these challenges at the source. Connected via hose to specialized accessories, they capture harmful particles at the source, minimizing airborne spread and keeping your workspace cleaner, healthier, and more efficient.

Designed by engineers working closely with dental technicians, these units reflect what matters most in the lab: reliability, performance, and simplicity.

Each unit features a proven multi-stage filtration system, including an antimicrobial bag filter and a certified, medical-grade HEPA filter that is 99.97% efficient at capturing particles 0.3 microns in size.

With a powerful 2,000-hour motor, these units offer long-lasting performance to keep your lab running safely and smoothly, day after day.





Buy SILENT Compact,  
Get 1 Filter Set Replacement\*

**SAVE** **OVER \$375**

## SILENT Compact Single-Workbench Extraction Unit

### Renfert USA

Single-station, extremely quiet extraction unit with a compact design for easy integration into workstation benches. Intuitive operation 1000-hour (guaranteed) motor, 2500 liters/minute with a quick-change motor replacement. Features bagless dust collection with easy calibration to the connected appliance. Mechanical filter cleaning and 99.9% filter efficiency. 120 volts. 3-year warranty included. Noise level (maximum): 56.9 dB. Dimensions: 17.3"H x 9.7"W x 19.7"D. Weight: 30 Lb.

Complete Unit  
(228-4782-FPD) .....Ea. **1,489.99**  
Contains: 1–2-meter suction hose & 2 end bushings.  
\*Ships with order.

**SAVE** **UP TO \$175**



## Vanguard General-Purpose Dust Collectors

### Vaniman Manufacturing

High-speed turbine technology that produces high suction. Unique Accumulator™ Pre-Filter that captures 95% of all dust. Whisper-quiet operation (55 dB) and low-maintenance design. 1-year warranty.

Vanguard 1X  
Single station and single speed.  
#10310–Single Station  
(787-3641-FPD) .....Ea. **679.99**  
Vanguard Gold 2X  
Single or 2 station and multispeed. High suction for heavy-duty use.  
#10320–Single Station  
(787-7281-FPD) .....Ea. **799.99**

**SAVE** **OVER \$250**



## Ductless Fume Hood (DFH)

### Quatro Air Technologies

Ideal selection when counter space is at a premium. This wall-mounted unit is designed to stop odors, harmful particulate matter, and vapors from migrating throughout the workplace.

Acrylic/Monomer Odors Unit–DFH-100-M  
Dimensions: 26"H x 28"W x 23"L  
(270-8400-FPD) .....Ea. **1,159.99**  
Contains: 1 of each: unit & wall-mounting bracket.  
Burnout Furnace Unit–DFH-1000-B  
For conventional burnout applications.  
(270-2933-FPD) .....Ea. **1,159.99**  
Contains: 1 of each: unit, wall-mounting bracket & filters.



**SAVE BIG**

## AIR FILTRATION

SAVE OVER  
**\$175**



**Pure Breeze  
HEPA Fume Hood**  
Fume, smoke,  
and odor  
capture.  
787-8931-FPD  
\$1,761.99

SAVE OVER  
**\$175**



**Pure Breeze  
HEPA Air Purifier**  
Dust, odors,  
& aerosols.  
787-4786-FPD  
\$1,761.99

## REMOVABLES & MODELS

SAVE OVER  
**10%**



**Splash Hood**  
Ideal for lathe  
or grinding.  
787-6045-FPD  
\$144.99

SAVE OVER  
**\$175**

**StoneVac  
Brushless SC**  
Heavy material  
dust collector.  
787-0498-FPD  
\$1,745.99



# VANIMAN

CELEBRATING 40 YEARS

THE **RIGHT** EQUIPMENT  
FOR EVERY AREA OF  
YOUR LAB

**Vanguard  
Platinum 2X**  
Premium  
dental dust  
collector.  
787-0886-FPD  
\$1,199.99



SAVE OVER  
**\$125**

**Downdraft  
Workstation**  
dust capture  
accessory.  
787-0114-FPD  
\$199.99



SAVE OVER  
**10%**

## CROWN & BRIDGE

## MILLING & PROCESSING

**Cut & Clean  
Cabinet**  
Zirconia  
puck post-  
processing.  
787-0446-FPD  
\$356.99



SAVE OVER  
**10%**

SAVE OVER  
**\$200**



**StoneVac  
Mill Sync**  
Dust  
collection  
for dental mills.  
787-0426-FPD  
\$1,939.99

**WE BUILD PRODUCTS THAT SOLVE PROBLEMS.**

Since 1985, we've worked closely with dental labs to build reliable equipment designed to keep work areas clean.

When it comes to dust collection and air filtration, we know our stuff. When you want quality, choose Vaniman.



USA MANUFACTURING



FROM THE WORLD'S  
#1 VANIMAN DEALER



**SAVE** OVER **15%**



### LED Task Lamp

#### Arbe Machinery

An innovative and powerful LED task lamp. Sophisticated LED lighting solution designed to provide premium performance for a number of task lighting applications. Because it provides exceptional lighting quality and ergonomic handling, it is the ideal task lamp for laboratory environments. Offers better directional lighting, produces less glare, and allows for precise illumination with no wasted light. Other features:

- Dimmable
- Fully positionable articulating arm
- 6000–6500°K
- 110–230 V
- USB charging port

(133-7664-FPD) ..... Ea. **199.99**

Buy 1 Labolight DUO,  
Get 1 GRADIA® PLUS Kit  
or OPTIGLAZE SET\*

**SAVE** OVER **\$550**



### Labolight DUO

#### GC America

- Dual-mode - All in one Step light/Cure
- LED light-curing unit for indirect composite resin
- Saves work space by adapting 2 functions into 1 unit

(333-4238-FPD) ..... Ea. **3,899.99**

\*To receive your no additional charge goods, send a copy of your Henry Schein invoice dated 09/28/25 – 12/27/25 to: GC America, 3737 West 127th Street, Alsip, IL 60803 or fax to: 1-800-423-2963, Attn.: Customer Care Group. No additional charge goods must be redeemed within 6 weeks of invoice date. Offers cannot be combined with any other offers or promotions. No substitutions allowed. Offers are valid in the U.S. and Canada only. No additional charge goods must be of equal or lesser value to the lowest selling unit price. Please allow 6–8 weeks for delivery.

**SAVE** OVER **\$925**



### Smartbox X2

#### Amann Girrbach

- Plaster dosing device that eliminates manual weighing and mixing of plaster and water
- Achieve permanent, exact dosing at the push of a button
- Unit is rapid and gram-exact
- Using the new dosage technology, Plaster is sprinkled very finely into the mixing beaker, which noticeably increases the homogeneity of the mixture
- 1-year warranty

Complete Unit

(139-0142-FPD) ..... Ea. **3,119.99**

Contains: 1 of each: Smartbox X2 unit with integrated plaster container, 3-liter water container, power cable, scoop, cover for scale & instruction manual.

**SAVE** OVER **\$150**



### Ultra Waxter 2™

#### Kerr

Electronic waxing unit that features dual spatulas. Comfortable, ergonomic, cushioned, gripped tips and programmable temperatures. Dual spatulas with separate temperature controls. Unit accepts spatulas and tips from the original Ultra Waxter™. Tip holder for convenience and organization. Easy-to-read LED screen. 2-year warranty on unit and 6-month warranty on tips and spatulas. Complete Unit

(236-0163-FPD) ..... Ea. **639.99**

Contains: dual spatula handles & 2 tips (Small PKT #35167 & Large PKT #35168).



100 years making work easy  
3 + 1 year guarantee

# Celebrating 100 years of making work easy

## Faster, Cleaner, Better Model Production Starts Here!

### die:master varnish

Red die varnish  
Item No. 228-2407-FPD

Blue die varnish  
Item No. 228-2406-FPD

Silver die varnish  
Item No. 228-2409-FPD

Gold die varnish  
Item No. 228-2408-FPD

~~\$21.00~~

**\$19.49**

SAVE OVER 15%



Isofix 2000 (1 liter, 2pk)

Item No. 228-6000-FPD

Plaster separating agent  
ideal for working models

~~\$106.00~~

**\$98.79**

SAVE OVER 15%

### Prepping

### Pinning



Bi-V-Pins

1,000 pk with Plastic Sleeve  
Item No. 228-6159-FPD

~~\$145.00~~

**\$135.99**

SAVE OVER 15%

1,000 pk with Metal Sleeve  
Item No. 109-3077-FPD

~~\$196.00~~

**\$182.99**

SAVE OVER \$30

### Mixing



Mixing Bowls & Paddles

200ml bowl with paddle  
Item No. 228-3558-FPD

~~\$255.06~~

**\$216.99**

500ml bowl with paddle  
Item No. 228-9363-FPD

~~\$265.61~~

**\$225.99**

700ml bowl with paddle  
Item No. 228-8112-FPD

~~\$271.13~~

**\$229.99**

SAVE OVER \$35



GO 2011 speed

Item No. 228-5559-FPD

Cleaning solvent for mixing bowls,  
impression trays, castings and more!

~~\$62.00~~

**\$57.79**

SAVE OVER 15%

### Vibrating



Vibrating Ball

Item No. 228-0452-FPD

~~\$136.86~~

**\$127.99**

SAVE OVER 15%

### Trimming



Renfert Trimmer Discs

80 Grit, 5 Pack  
Klettfix: 228-6617-FPD

~~\$92.83~~

**\$86.79**

SAVE OVER 15%

Partially Diamond Coated  
Marathon Disc: 228-2209-FPD

~~\$879.97~~

**\$819.99**

SAVE OVER \$150

Fully Diamond Coated  
Infinity Disc: 228-2157-FPD

~~\$1,722.73~~

**\$1,589.99**

SAVE OVER \$300

Offers valid 9/28/2025 - 12/27/2025.



Buy 1 Henry Schein Vest Powder,  
Get 2- Liquid Only 1 liter\*

**SAVE OVER \$75**



## Henry Schein Vest

**HENRY SCHEIN®**

Investment for all alloys and pressable ceramic;  
for use in speed and overnight technique.

- High-speed bench set in 15 minutes
- 5–7 minutes normal working time
- Casting within 1 hour
- Single liquid for all applications
- Easy divesting

Powder Only, 100 Gm. Pouches ..... 100 per Box **169.99**

Liquid Only  
(900-2957-FPD) ..... 1 Liter **29.99**

\*Ships with order.

**SAVE OVER 10%**



## Bellovest® SH

**BEGO USA**

- Crown-and-bridge investment material that offers versatility and flexibility.
- Mixing liquid: BegoSol® HE
- Processing time at 68°F (20°C): approximately 4½–5 minutes
- Shelf life in unopened bag: 24 months
- Beginning of solidification (Vicat time): approximately 10 minutes
- Compressive strength after 2 hours: 4.2–5.1 MPa
- Linear thermal expansion: 0.85%
- Flowability: approximately 140mm–145mm

Bellovest® SH Powder, 100 Gm. Pouches ..... 100 per Pkg. **96.99**

BegoSol® HE Liquid  
(412-2242-FPD) ..... 1 Liter **30.99**

(412-0151-FPD) ..... 5 Liters **118.99**

Buy 3,  
Get 1 any type of equal or less value\*

**SAVE OVER \$250**



## GC FujiVest® II

**GC America**

Rapid-fire, carbon-free, phosphate-bonded investment that produces ultrasmooth, precise castings quickly and easily. It is easy to divest and resistant to fracture, making it suitable for all crown-and-bridge dental alloys. Once it has set (20 minutes after investing), it can be placed directly into a preheated furnace at final temperature. There is no need for progressive heat-up: the investment can be fired immediately after bench setup time.

90 Gm.  
(333-6790-FPD) ..... 67 per Pkg. **128.99**

Contains: 67–90 Gm. powder pouches.

2 Kg.  
(333-3455-FPD) ..... 4 per Pkg. **169.99**

GC FujiVest II Liquid  
(333-6417-FPD) ..... 900 mL **39.99**

\*To receive your no additional charge goods, send a copy of your Henry Schein invoice dated 09/28/25 – 12/27/25 to: GC America, 3737 West 127th Street, Alsip, IL 60803 or fax to: 1-800-423-2963, Attn.: Customer Care Group. No additional charge goods must be redeemed within 6 weeks of invoice date. Offers cannot be combined with any other offers or promotions. No substitutions allowed. Offers are valid in the US and Canada only. No additional charge goods must be of equal or lesser value to the lowest selling unit price. Please allow 6–8 weeks for delivery.

**SAVE OVER 10%**



## Formula 1

**Ransom & Randolph**

Universal Investment for All Alloys and  
Pressable Ceramics

This high-performance, universal investment works equally well for all alloys and pressable ceramics. Can be used in rapid or standard burnout due to its wide expansion rate. Robust bimodal nature strengthens molds and promotes easier divesting when working on more fragile materials such as all-ceramic restorations. Using only a 15-minute bench set and 30-minute burnout, Formula 1 reduces processing time and increases lab productivity.

(898-0768-FPD) ..... 100 x 100 Gm. **187.99**

(898-0773-FPD) ..... 1 Liter **39.99**



# GC LiSi PressVest

Phosphate Bonded Investment for Press Ceramic Techniques

## Offering you an easier, faster and safer process

- Easy removable of the reaction layer - reduction of operating time, reproducible fit due to no use of hydrofluoric acid
- Easy divesting - reduction of operating time
- Super high fluidity and long working time - accurate and detailed reproduction
- Precise expansion control - optimal fit
- Optimized for quick heating technique
- Adapted to multiple press ceramic technique

### GC LiSi PressVest

[3334326-FPD] Powder Refill x60 (100 g): \$199.99

[3334328-FPD] Liquid (900 mL): \$34.99



Buy 3 GC LiSi PressVest Powder Packages, **Get 1!\***



800.323.7063 • [www.gc.dental/america](http://www.gc.dental/america) • [www.gc.dental/america/training](http://www.gc.dental/america/training) • © 2025 GC America Inc.

\*Promo code 3067. To receive your no charge goods, send a copy of your Zahn Dental invoice dated 9/28/25 – 12/27/25 to GC America, 3737 West 127th Street, Alsip, IL 60803 or fax to 800-423-2963 Attn: Customer Care Group. Terms and Conditions of ALL Promotions: **Limit three redemptions per customer without prior written authorization. All offers expire 12/27/2025.** Promotions cannot be combined with any other promotional offers. No substitutions allowed. Offers are valid in the US only. No charge goods must be redeemed within six weeks of Zahn Dental invoice date. No charge goods will be of equal or lesser value to the lowest selling unit price. Please allow 6-8 weeks for No charge goods delivery. All no charge goods will be fulfilled by GC America Inc. except where noted. Please indicate requested no charge goods on Zahn Dental invoice. Only no charge good requests that come directly from the Dental Office/Dentist and/or Staff will be honored. Qualifying purchases must appear on a single Zahn Dental invoice.

California customers please use website when ordering for Prop 65 information.

**SAVE OVER 20%**



## Dispos-A-Trap

**HENRY SCHEIN®**

- Reliable, sanitary alternative to metal plaster traps
- Made from high-density polyethylene (HDPE) plastic
- Disposable translucent container allows for easy monitoring of waste level
- Easy installation under most sinks
- Flexible PVC connection hoses simplify installation and require only 9" between top of container and bottom of sink pipe

Complete Unit—3.5 Gallons

(101-9233-FPD) ..... Ea. **49.99**

Contains: 1 plaster trap & 2 connection hoses & cuffs.

Complete Unit—5 Gallons

(101-0868-FPD) ..... Ea. **51.99**

Contains: 1 plaster trap & 2 connection hoses & cuffs.

Refill—3.5 Gallons

(101-3463-FPD) ..... Ea. **25.99**

Contains: 1 plaster trap & 2 end caps.

Refill—5 Gallons

(101-2412-FPD) ..... Ea. **29.49**

Contains: 1 plaster trap & 2 end caps.

**SAVE OVER \$350**



## Maxima® XL2 Micromotor

**HENRY SCHEIN®**

Brushless micromotor that features a streamlined, modern design with a small footprint to fit in most work places. Quieter and less vibration than previous model. Rotation speed can be effortlessly controlled between 1000 and 50,000 rpm. Smooth, stable, and powerful torque is delivered across the entire speed range, which is displayed digitally. An easy-to-hold handpiece design ensures precision. Brushless motor does not require replacing the carbon brushes. Built-in torque meter provides real-time feedback and error-code display for overload protection. Easy-to-read digital display. Dust-proof mechanism. 1-year warranty.

Complete System

(570-0624-FPD) ..... Ea. **959.99**

Contains: 1 control unit, 1 handpiece, 1 foot control, 1 handpiece stand & power cables.

**SAVE OVER \$125**



## Ultrasonic Cleaner with Touch Pad

**HENRY SCHEIN®**

One-touch, press-set digital pad timer in beige, vinyl-clad housing. Solid-state circuitry. Stainless steel drain with reversible right or left hookup. 43 kHz of cleaning power. Includes an FCC-approved EMI filter and a 60-minute timer with infinity hold. 110-volt unit. 2½-year warranty.

Tank dimensions: 9⅞" x 5⅞" x 4"

Tank capacity: 3⅔ Qt.

Complete Unit

(104-4057-FPD) ..... Ea. **459.99**

Contains: cover & drain hose assembly. Basket sold separately.

**SAVE OVER 20%**



## The Vibrator

**HENRY SCHEIN®**

- Improved 3-speed switch and power circuitry: provides more power and reliability
- Durable 4"-diameter rubber platform: provides more work area and quieter operation
- Soft, quiet suction-cupped feet: quiet and stable operation
- Small size: convenient

120 Volt

(900-4053-FPD) ..... Ea. **115.99**

## NTI® Diamond Discs

NTI diamond discs are precision instruments which are manufactured by applying the UniMatrix diamond coating technology.

The discs are firmly riveted to the mandrel. The various thicknesses of the disc blanks which range from 0.08 mm to 0.3 mm together with the grit sizes of the diamond grains allow for a variety of applications.

Available in individual packs.



FOR A LIMITED TIME  
**BUY 4, GET 1  
 AT NO COST!\***  
 (MIX & MATCH)

Valid: Sept 28 - Dec 27, 2025 | Kerr Code: ZNHSLB925

## NTI® Universal Cutters

NTI® Universal Cutters incorporate sharp, precision ground tungsten carbide blades with stainless steel shanks to ensure concentricity.

NTI® Universal Cutters are color coded to indicate blade configurations, and will cut smoothly on stone, acrylic, metals and porcelain without clogging or chattering.

Available in individual packs.



FOR A LIMITED TIME  
**BUY 4, GET 1  
 AT NO COST!\***  
 (MIX & MATCH)

Valid: Sept 28 - Dec 27, 2025 | Kerr Code: ZNHSLB925

\* Promotion valid for standard accounts only. Promotional goods must be of equal or lesser value, cannot be combined with other offers. Purchase must be made on one invoice between 9/28/25 and 12/27/25, unless otherwise noted. Redeem by 1/27/26. To receive your promotional goods, email or fax a copy of your authorized Henry Schein invoice noting promo code **ZNHSLB925** and product(s) desired to: [kerrpromo@kerrdental.com](mailto:kerrpromo@kerrdental.com) or FAX: 888.727.2614. Limit to 3 redemptions. Incomplete submissions will not be processed. Allow 8 weeks for delivery. Offer valid in the 50 United States and the District of Columbia only. Promotions are subject to change or cancellation without notice. Offer void if purchased product is returned. Note that you may have an obligation under federal, state or local law to reflect discounts on product given pursuant to this promotion on any cost report forms submitted to a federal or state government or private payer who provides reimbursement for that product. ©2025 Kerr Corporation. All trademarks are property of Kerr Corporation. All Rights Reserved. MKT-25-1307 Rev 0

**SAVE** OVER **15%**



## CRITERION™ Latex Polymer Coated

**CRITERION**

Latex Powder-Free Exam, Non-Sterile

Finished with a smooth interior polymer coating for fast and easy donning. Comfortable fit and feel for a superior range of motion. Textured finish for an enhanced grip. Beaded cuff. 3.15-mil thickness at the palm. White.

Gloves .....100/Box **8.49**  
10 Boxes (Mix & Match)@ **7.99**

Specify:

XS (570-3412-FPD) L (570-3415-FPD)  
S (570-3413-FPD) XL\* (570-3416-FPD)  
M (570-3414-FPD)

\*90/Box

**SAVE** OVER **15%**



## CRITERION™ Nitrile N100

**CRITERION**

Formerly known as Criterion Nitrile. Same Great Product, New Criterion™ Packaging Nitrile Powder-Free Exam

- Formulated with 100% synthetic nitrile, which is odor-free and not made with natural rubber latex, eliminating potential allergic reactions to natural latex protein
- Superior puncture-, tear-, and chemical resistance, with excellent tactile sensitivity and dexterity
- Textured surface provides a strong grip for both wet and dry applications
- Eco-friendly box design

Gloves .....100 per Box **7.49**  
10 Boxes (Mix & Match)@ **6.99**

Specify:

XS (111-8533-FPD) L (111-8537-FPD)  
S (111-8535-FPD) XL\*\* (111-8539-FPD)  
M (111-8536-FPD)

\*\*90/Box

**SAVE** OVER **15%**



## CRITERION™ Earloop Masks Level 3

**CRITERION**

- Light and comfortable for extended wear
- Filter media is surrounded by soft, nonirritating inner and outer facing
- Full-width flexible nosepiece provides fit and security
- Not made with natural rubber latex

Masks.....50 per Box **12.49**  
6 Boxes (Case)@ **11.99**

Specify:

Blue (570-2044-FPD) White (570-2043-FPD)  
Teal (570-2042-FPD)

**SAVE** OVER **15%**



## CRITERION™ Earloop Masks, Level 1

**CRITERION**

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

Masks.....50 per Box **9.49**  
6 Boxes (Case)@ **8.99**

Specify:

Blue (104-3809-FPD) White (104-6611-FPD)

**SAVE** OVER **15%**



## Ultrasonic Cleaning Solution

**HENRY SCHEIN®**

General Purpose Non-ammoniated Cleaner

(100-3836-FPD) .....1 Gal. **24.99**  
4 x 1 Gal. **23.99**

Plaster and Stone Remover ♦  
(101-3277-FPD) .....1 Gal. **24.99**  
4 x 1 Gal. **23.99**

♦ Hazardous material



# Lab Essentials for Less

Stock up on the supplies you love from a name you can trust.



## GI-MASK Automix New Formula

GI-MASK Automix New Formula provides accurate reproduction of soft tissue structure. The A-Silicone based gingival mask provides all information for implant supra-structures and frameworks.

**BUY 2** GI-MASK Automix New Formula refill 2 x 50 ml

**GET 1** of the same!  
promo code: CS2521

GI Mask Automix Refill C8065 Ea  
Item # 888-9675-FPD  
Price per unit after promotion: \$ 121.99



## Lab Putty

A silicone putty especially designed for the use in the dental laboratory. Lab Putty impresses with its outstanding properties and various final hardnesses.

**BUY 2** Standard Lab Putty

**GET 1** of the same!  
promo code: CS2523

Lab Putty Kit 900mL/Bt  
Item # 888-8192-FPD  
Price per unit after promotion: \$ 97.99


Lab Putty Base Only 7500mL  
Item # 888-8165-FPD  
Price per unit after promotion: \$ 319.99



 **COLTENE**

**Terms & Conditions:** Promotions valid September 28, 2025 through December 27, 2025. To qualify for **no charge** goods on individual promotions, full quantity must be purchased at one time. Eligibility is restricted to the United States and Canada (excludes US territories). May not be combined with other Coltene promotions unless otherwise noted. **No charge** goods may be limited to most popular shades or configurations within a product family. **Redemption Information:** Upload invoice and complete **no charge** good request form at [rebates.coltene.com](https://rebates.coltene.com) by January 16, 2026 to qualify. Incomplete submissions will not be processed. Invoices with backordered product will not be honored. Orders must be shipped in full. Promotions can be fulfilled with invoice at the following website: [rebates.coltene.com](https://rebates.coltene.com). Or, via email at [free.goods@coltene.com](mailto:free.goods@coltene.com)

**SAVE** **OVER \$1,250**




### FOCUS 4010 Porcelain Furnace with Pump ShenPaz

- Precise vacuum calibration: Vacuum precision ensures translucency and vitality of the restoration; eliminates overbrightness or opaqueness
- Preprogrammed tempering process adjusts the right coefficient for zirconia, metal, and ceramic materials
- 200 freely programmable memory locations individually preprogrammed on request
- Intuitive full-color digital display
- 3 preprogrammable cooling steps
- 4 fully programmable drying steps
- 200 programs
- Display setup features choice of 8 languages
- USB connection
- 115 volts
- White

(141-2155-FPD) ..... Ea. **4,749.99**

**SAVE** **OVER \$325**



### Vulcan® Three-Stage Burnout Furnace Ransom & Randolph, LLC

Each 3-stage digital programmable model has 2 standard muffle sizes, a 3-series controller, 9 multistage programs, a hold-temperature program, and a delay-start option, which allows a program's start time to be delayed for up to 7 days. Unit is constructed with a heavy-gauge steel cabinet and is finished with a durable, baked-enamel coating. Vertical lift door directs heat away from the operator and is standard for safety and space conservation.

- High-performance hybrid muffle
- Longer life and more durable than fiber-constructed furnaces
- Faster heating and cooling than full-firebrick furnaces
- Integrated door safety switch
- Durable, stainless-steel front panel
- Operating temperature range: 77°F–2012°F (25°C–1100°C)
- Programmable hold time: 0:00–99:59 hours
- 1-year limited warranty
- Made in the USA

Model #3-130

120-V Unit

(140-7694-FPD) ..... Ea. **2,949.99**

Model #3-550

120-V Unit

(140-7696-FPD) ..... Ea. **3,259.99**

**SAVE** **OVER \$100**



### Presto II Air Turbine Handpiece System NSK

- Lubrication-free
- 320,000 min-1
- Water coolant spray
- Designed specially for precise cutting of ceramic
- Decreases dispersal of debris and chips
- Twist-free swivel
- Unique dust-proof mechanism
- NSK LED is equivalent to daylight

Complete Unit

(272-3197-FPD) ..... Ea. **799.99**

Contains: 1 of each: unit, handpiece, handpiece stand, primary air hose, secondary air hose & foot control.

**SAVE** **OVER \$150**



### Presto Aqua II Air Turbine Handpiece System NSK

- Lubrication-free air-turbine handpiece with water spray
- Designed for precise cutting of ceramic
- Coolant water spray enables shorter operation time and longer bur life
- Two water spray options (droplets or fine-mist spray)
- Unique dust-proof mechanism to protect internal components from debris. 320,000 rpm

Complete Unit

(272-3102-FPD) ..... Ea. **1,159.99**

Contains: 1 of each: unit, handpiece, primary air hose, secondary air hose, coupler & foot control.

**Renfert**

100 years making work easy  
3 + 1 year guarantee

# Celebrating 100 years of making work easy



## LIGHT 1

Experience light in its ideal form

SAVE  
OVER  
\$150

Table mount version  
Item No.: 228-6627-FPD

**\$1,439.99**

SAVE  
OVER  
\$140

Wall mount version  
Item No.: 228-6631-FPD

**\$1,258.99**



3 + 1 year guarantee


Register now:  
[renfert.com/three-plus-one](https://renfert.com/three-plus-one)

### Mounting Accessories

1. Base stand, LIGHT 1  
Item No.: 228-6635-FPD

**\$219.99**

SAVE OVER \$30



2. Table clamp, LIGHT 1  
Item No.: 228-6633-FPD

**\$105.99**

SAVE OVER 10%



## SILENT XS

The ultimate mobile dust collector for All-on-X procedures.

Item No.: 228-6638-FPD

Designed specifically for the demands of chairside adjustments during All-on-X restorations, this compact, cordless, and powerful unit captures airborne dust at the source, ensuring a cleaner, healthier environment for both patient and practitioner.

**\$499.99**

SAVE OVER \$50



Offers valid 9/28/2025 - 12/27/2025.



Buy 3,  
Get 1 any type of equal or less value\*

**SAVE OVER \$175**



## GC Fujivest® Premium

### GC America

Suitable for all types of dental alloys. Accommodates ringless techniques in both quick-heating and slow-heating procedures. Optimal expansion control by liquid dilution ensures a highly accurate fit. Carbon-free phosphate-bonded investment. 90 Gm.

(333-2130-FPD) .....67 per Pkg. **126.99**

Liquid Only

(333-2134-FPD) .....900 mL **39.99**

\*To receive your no additional charge goods, send a copy of your Henry Schein invoice dated 09/28/25 – 12/27/25 to: GC America, 3737 West 127th Street, Alsip, IL 60803 or fax to: 1-800-423-2963, Attn.: Customer Care Group. No additional charge goods must be redeemed within 6 weeks of invoice date. Offers cannot be combined with any other offers or promotions. No substitutions allowed. Offers are valid in the US and Canada only. No additional charge goods must be of equal or lesser value to the lowest selling unit price. Please allow 6–8 weeks for delivery.

**SAVE OVER 15%**



## Vitex™ Investment

### Garreco

Ethyl-silicate, partial-denture framework investment material. Developed especially for producing chrome cobalt castings. Offers precise details due to an ultrafine particle size and unique 2-part blend of liquids. May be used for all high-heat techniques. Formulated to mix and pour easily, it maintains excellent integrity when used in conjunction with Vitex Model Hardener.

- Precise detail
  - Easy mix/pour
- Powder & Liquid Kit, Small

(575-0332-FPD) ..... Ea. **109.99**

Contains: 1–50 Lb. carton of powder investment, 4–28 oz. bottles of “A” liquid & 4–3.7 oz. bottles of “B” liquid.

“B” Liquid

(575-0154-FPD) .....4 oz. **13.99**

(575-0465-FPD) .....14.8 oz. **46.99**

**SAVE OVER 10%**



## Wirovest®

### BEGO USA

Phosphate-Bonded Cobalt-Chromium Refractory Investment  
Part of the Bego line of accurate refractory investment for cobalt-chrome alloys. Produces consistent, ultrasmooth model surfaces and smoother castings. Also yields predictable and controllable expansions when used with Begosol® mixing liquid.

Powder, 400 Gm. Pouches  
(412-1159-FPD) .....45 per Pkg. **129.99**

Begosol® Liquid, Clear

(412-8935-FPD) .....1 Liter **29.99**

**SAVE OVER 10%**



## WiroFine

### BEGO USA

Universal investment material for all indications in the partial denture and combination technique, for gel or silicone duplication. Can be heated rapidly or conventionally to 1,050°F with ideal expansion values. Offers the level of flexibility essential for the modern dental laboratory. Ideal flow properties make for reliable, fatigue-free working; even the finest areas are precisely reproduced. Excellent deflasking properties due to the minimal reaction between the investment material and the alloy.  
Powder Envelopes  
400 Gm.

(623-7215-FPD) .....50 per Pkg. **188.99**

Begosol® K Mixing Liquid

(623-3020-FPD) .....1 Liter **29.99**

(412-0111-FPD) .....5 Liters **115.99**

Begosol® Liquid, Clear

(412-8935-FPD) .....1 Liter **29.99**

**SAVE 10%**



## Multi-Vest®

### Ransom & Randolph

- A phosphate-bonded dental investment for casting cobalt-chromium and other high fusing precious and non-precious alloys
- Phosphate-bonded dental investment designed for casting cobalt-chromium and other high fusing precious and non-precious alloys
- Provides a smooth casting finish with a short setting time, reducing the risk of flashing
- Ideal for casting crowns, bridges, and partial denture frames
- Color: Pink
- Distilled water is required for optimal results. Universal Expansion Liquid can be used
- to adjust expansion as necessary, but is not required for use with for Multi-Vest investment.

Powder

(140-8529-FPD) .....33 Lb. **82.99**



**Renfert**

100 years making work easy  
3+1 year guarantee

# Celebrating 100 years of making work easy



**POWER steamer 1**  
(Manual filling)

Item No.: 228-6145-FPD

~~\$2,475.00~~

**\$2,089.99**

SAVE OVER \$300



**POWER steamer 2**  
(Automatic filling)

Item No.: 228-6144-FPD

~~\$3,410.00~~

**\$2,879.99**

SAVE OVER \$400



\*Stand sold separately



**Twister**  
Vacuum mixing unit

Item No.: 228-5266-FPD

~~\$1,978.00~~

**\$1,659.99**

SAVE OVER \$300

*Renfert*



**Top Spin**

Pin drilling unit

Item No.: 228-2627-FPD

~~\$1,705.00~~

**\$1,399.99**

SAVE OVER \$250



**MT3 Klettfix**

(80 grit, Corundum disc)

Item No.: 228-2064-FPD

~~\$1,375.00~~

**\$1,149.99**

SAVE OVER \$175



**MT3 Marathon**

(Partially diamond-coated disc)

Item No.: 228-2071-FPD

~~\$2,145.00~~

**\$1,844.99**

SAVE OVER \$300

Offers valid 9/28/2025 - 12/27/2025.

# Master **EVERY** Detail **No Bubbles! No Voids!**



SAVE OVER 20%

**\$22.49**

(510-0095-FPD)

- ⌚ Saves time
- ⌚ Noncontaminating
- ⌚ Reduces surface tension
- ⌚ Produces bubble-free castings
- ⌚ Produces smooth and accurate models

**dang good  
debubbler**

SAVE UP TO 20 MINUTES!

**wonderfill®**  
THE EZ IN & OUT TONGUE  
& VOID FILLER!

- ⌚ Fills in bubbles, air pockets, missing teeth and undercuts
- ⌚ Enhances the fit of appliances and bleaching trays
- ⌚ Excellent in creating a temporary dam when making partial molds



SAVE  
OVER 20%

**\$47.99**


Net Wt. 1.13 kg  
(510-0210-FPD)



from the makers of Wonderfill®

**dentalcreations** Ltd.

Effective dates: 9/28/25 – 12/27/25



Cool Down with  
**0% Interest**  
All Month Long

Take advantage of refreshing savings on everyday  
lab supplies, CAD/CAM equipment, and more.

**\$0 Down | 0% Interest | Up To 12 Months Financing**

*Limited time offer—don't miss out!*



www.zahndental.com | (800) 496-9500

California customers please use website when ordering for Prop 65 information.

Some exclusions may apply. Contact your Zahn Sales Consultant for more details. 25ZS7855








### Workbox Orion

907-6560-FPD

- High luminosity LED lamp.
- Absence of shadows inside.
- Ergonomic work position.
- Optimal absorption of powder by means of external suction (not included).
- Easy cleaning.
- Includes a wedge in the interior with support for micromotor.



#### Specifications

Height 11.8 in • Length 14.17 in  
Width 14.2 in • Weight 5.7 lb • Interior light

**\$ 479.99**  
Save over **\$ 135.00**

### Sandblaster Eclipse II

685-0258-FPD

Cabin entirely made of ABS of reduced dimensions maintaining the amplitude in the work area and with a substantial improvement in lighting.

#### Specifications (685-0257-FPD)

Height 19.29 in • Width 14.17 in  
Depth 17.71 in • Weight 17lb  
Maximum air consumption 26.4 gal/min.  
Operating pressure 2 a 5 atm.  
Interior light

### Sandblaster Eclipse II-T

685-0257-FPD

They have two mobile micro-blasting nozzles and, in the hopper version, a large sand tank and a 0.11 in static gun.

#### Specifications (685-0257-FPD)

Height 19.29 in • Width 14.17 in  
Depth 17.71 in • Weight 17lb  
Maximum air consumption 26.4 gal/min.  
Operating pressure 2 a 5 atm.  
Interior light



Eclipse II

**\$ 1,029.99**  
Save over **\$ 289.00**



Eclipse II-T

**\$ 1,129.99**  
Save over **\$ 339.00**

### Polymerizer for flasks M9

685-0204-FPD

#### Stainless steel

Rust-free frame. This machine can hold up to 9 flasks for polymerizing. Synoptic and intuitive programming control, with delayed start time in hours and minutes and with warning buzzer at programs end.

High heating power (1750 W). Easy draining through a lateral tap. Power supply: 115 V.



**\$ 2,549.99**  
Save over **\$ 287.00**

### Washing and polymerizing machine E2

685-0201-FPD

#### Stainless steel

This machine can hold up to 6 flasks for dewaxing and up to 12 flasks for polymerizing. Once the polymerization is over, just lay the flasks on the upper tray, and warm water shower will dewax them. Different programs are available at your choice. Power supply: 220 V.



**\$ 5,249.99**  
Save over **\$ 577.00**

### Hydraulic Press

685-0197-FPD

Made of high resistant aluminum. Robust design (up to 15 ton force). It admits up to three flasks with clamps. Astonishing smooth-spinning spindle. Very well balanced workbench.



**\$ 1,029.99**  
Save over **\$ 84.00**

### Gel mixer Aneto

685-0269-FPD

#### Stainless steel

Attractive and reliable gel mixer featuring state-of-the-art technology: tank heating by lateral and bottom resistors, which prevents gel burning from high temperatures; oxide-free casing, made of stainless steel and polyurethane; high capacity

(Approx. 1.85 Gal tank); state-of-the-art electronic controls allowing multiple adjustments... A complete machine designed to work without problems for many years.  
Power supply: 115 V.



**\$ 4,019.99**  
Save over **\$ 453.00**

### Polishing Box with suction

685-0218-FPD

#### Stainless steel

Rust-free frame, made of stainless steel. Supplied with powerful Lathe 160 W. Two speed and with a great flow suction motor. Minimum noise. Simple access to the suction filters system, for an easy maintaining and cleaning. The unit is equipped with LED lighting and with two removable rubber hoods. With manual safety reset switch to prevent unexpected start-ups. Power supply: 115 V.



**\$ 4,939.99**  
Save over **\$ 368.00**

### Box with suction for polishing

685-0205-FPD

Robust and rust-free frame. It features and exhauster and filter unit, lighting system, polycarbonate protective screens and a Schuko-type power connection with manual safety reset to prevent unexpected start-ups. Supplied with two rubber screens. Power supply: 115 V.



**\$ 2,709.99**  
Save over **\$ 307.00**

### Automatic polymerizing pot Geisser

685-0253-FPD

Polymerizer pot with high capacity tank (up to 1,58 gallon). Electronic controls to regulate temperature and time. Useful for the vast majority of materials (polymerizable by pressure and temperature) present in the market. The quality of the result is simply perfect. Power supply 115V.



#### Specifications

Dimension 13.77 x 15.75 x 14.17 in  
Weight 11lb • Max. pressure 2 bar  
Max. Temperature 212°F  
Volume tank 1.6 gal.  
Heating power 600 W

**\$ 1,649.99**  
Save over **\$ 200.00**

### Alginate Mixer

685-0282-FPD

Useful tool when making alginate paste in the dental clinic. The user can control visually the process at every stage, and this guarantees that the paste made is in optimal conditions for application. In order to make cleaning easier the bowl is detachable. It is factory supplied with a bowl and a special plastic spatula.



#### Specifications

Dimensions 9.06 x 7 x 9.85 in  
Weight 6.6 lb  
Power 60 W

Bowl 750 cc

Diameter 5.12 in  
Height 4.33 in  
Weight 4.4 lb

**\$ 529.99**  
Save over **\$ 97.00**

### Steam cleaning box Altamira

685-0198-FPD

The indispensable complement for any steam cleaning device. Made of polypropylene, very strong and rust-free material. Wide chamber, LED lighting, powerful ventilation system, integrated water tank at the bottom.

With an easy and practical tank-draining system.  
Power supply: 115 V.



**\$ 1,239.99**  
Save over **\$ 230.01**

# Image APN®

Highly Wear Resistant Resin Teeth

## The Smarter Choice in Carded Teeth – Image APN®

Henry Schein® Brand **EXCLUSIVE**



**Image APN® Tooth Stock Bundle: \$1,699.99**

**SAVE OVER \$1,200**

**Bundle Includes;**

- 216 most popular moulds in shades A1, A2, A3, B1 & C1
  - > Upper Anterior: LD21, LE12, LE22, LH11, LX21, LG42, LG12, LF42, LH45, LD42
  - > Lower Posterior: 29M, 31M, 32M, 33M, 10/32, 10/30, 10/34, 0/30, 0/32, 10/36
  - > 72 Upper & Lower Posteriors: 29M, 31M, 32M, 33M, 10/32, 10/30, 10/34, 0/30, 0/32, 10/36
- 10-drawer cabinet with inserts
- 5 shade guides
- USB key chain scanner for easy restocking
- Barcode labels for inventory management



> 3-layer posterior and anteriors.

> Extended mould and shade range in the most popular U.S. moulds and shades.

> APN® (Advanced Polymer Network) formulation backed by lifetime warranty.

Promotions cannot be combined with other offers or discounts. Valid until 12/27/2025. All transactions are subject to the satisfaction of underwriting guidelines, credit approval, and documentation requirements, and not all applicants will qualify. Certain other restrictions and additional terms and conditions apply.

**HENRY SCHEIN®**  
BRAND PRODUCTS

**Claim Your Bundle Today!**  
Limited Time Offer





# Pick Your Mould. We'll Handle the Rest.

## Dual Form V Tooth Stock Bundle: \$999.99

**Save over \$500**

### Bundle Includes;

- 360 sets of Dual Form V moulds
  - > Choice of any Dual Form V mould and VITA shade
- 10-drawer cabinet with inserts
  - > Holds 44 anteriors or 28 posteriors
- 5 shade guides
- USB key chain scanner for easy restocking
- Barcode labels for inventory management

**DUAL FORM V®**



**ZDENT™**  
DENTURE TEETH

## Z-Dent Tooth Stock Bundle: \$729.99

**Save over \$300**

### Bundle Includes;

- 385 sets of Z-Dent moulds
  - > Choice of any Z-Dent mould in VITA and Classic Shades
- 10-drawer cabinet with inserts
  - > Holds 44 anteriors or 33 posteriors
- 5 VITA or Classic shade guides
- USB key chain scanner for easy restocking
- Barcode labels for inventory management



**HENRY SCHEIN®**  
BRAND PRODUCTS

**Claim Your Bundle Today!**  
Limited Time Offer

# DIAMOND D®

ULTRA IMPACT DENTAL ACRYLIC

## Leading the Way in Trusted Denture Solutions



Advanced natural-looking dentures require artistry as well as technical expertise. Diamond D® was specially formulated for exceptional working and handling characteristics for anatomic sculpting.



### BUY POWDER, GET LIQUID!

#### Save Over \$75!

Buy **5 lbs** of Diamond D® heat cure denture acrylic, get 1 qt. of liquid at no charge!\*

#### Save Over \$275!

Buy **25 lbs** of Diamond D® heat cure denture acrylic, get 4 qts. of liquid at no charge!\*

\*No charge liquid ships with order



Original

25 lb.....242-4291-FPD.....**\$774.99**  
5 lb.....242-5462-FPD.....**\$233.99**



Light

25 lb.....242-4112-FPD.....**\$774.99**  
5 lb.....242-9520-FPD.....**\$233.99**



Light  
Reddish Pink

25 lb.....242-0574-FPD.....**\$774.99**  
5 lb.....242-5109-FPD.....**\$233.99**



Chroma  
Essence

25 lb.....242-5288-FPD.....**\$774.99**  
5 lb.....242-0424-FPD.....**\$233.99**



Dark  
Veined

25 lb.....242-0115-FPD.....**\$774.99**  
5 lb.....242-0035-FPD.....**\$233.99**

keystone®  
industries

To order, visit [www.zahndental.com](http://www.zahndental.com) or call 800.496.9500  
California customers please use website when ordering for Prop 65 information.

This offer cannot be combined with any other promotion or discount. No limit on purchases. Other restrictions may apply. Customer will be invoiced for purchased goods from Zahn Dental and no charge liquid will ship with the same order. Promo is valid from 09/28/2025-12/27/2025 and in the US only. Not responsible for typographical errors.

# SAVE TIME. REDUCE LABOR.

## *Increase Consistency*

The AutoFinish™ Finishing and Polishing Machines streamline denture production — whether you're working with digital or traditional workflows.

### WHY LABS CHOOSE AUTOFINISH™

- **Finish dentures faster**—automated processing cuts manual labor dramatically
- **Supports both digital and traditional dentures** with precision and repeatability
- **Free up technician time** for higher-value tasks—less grinding, more productivity
- **Consistent, high-quality results** that reduce remakes and improve turnaround times
- **Pays for itself in under 12 months** with daily use on dentures and other restorations

#### 3D PRINTED DIGITAL DENTURE

Before

After



#### ANALOG DENTURE

Before

After



VALIDATED FOR CROWN & BRIDGE, REMOVABLES, IMPLANT WORKFLOWS



We Are Your Total **Solutions** Provider  
**zahndental.com** | 1.800.496.9500

California customers please use website when ordering for Prop 65 information.

© 2025 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. 25ZS2655



**SAVE** **OVER 15%**



## Vacuum-Forming Materials

**HENRY SCHEIN®**

5" x 5" Thermoforming Sheets

- Used with vacuum-forming machines for various dental appliances
- Consistent thickness: uniform appliances

Base Plate Material—5" x 5"

Pink—Hard—0.060"

(100-4623-FPD) ..... 25 per Box **16.99**

Pink—Hard—0.060"

(100-2306-FPD) ..... 625 per Box **262.99**

Coping Material—5" x 5"

White—Soft—0.020"

(100-3513-FPD) ..... 50 per Box **13.99**

0.030"

(101-5448-FPD) ..... 50 per Box **20.79**

Mouth Guard Material—5" x 5"

Clear—0.040"

(101-8910-FPD) ..... 12 per Box **13.49**

Clear—0.060"

(101-7962-FPD) ..... 12 per Box **14.49**

Clear—0.080"		
(100-7718-FPD) .....	12 per Box	<b>14.49</b>
Clear—0.120"		
(101-2123-FPD) .....	12 per Box	<b>14.99</b>
Clear—0.150"		
(100-1199-FPD) .....	12 per Box	<b>17.99</b>
Surgical Tray Material—5" x 5"		
Clear—0.060"		
(100-1330-FPD) .....	25 per Box	<b>22.79</b>
Clear—0.080"		
(100-4239-FPD) .....	25 per Box	<b>25.49</b>
Clear—0.100"		
(101-4413-FPD) .....	25 per Box	<b>31.79</b>
Temporary Splint Material—5" x 5"		
Clear 0.020"		
(100-3822-FPD) .....	50 per Box	<b>24.49</b>
Tray Material—5" x 5"		
Pink—0.080"		
(100-2184-FPD) .....	25 per Box	<b>20.79</b>
Pink—0.080"		
(100-9127-FPD) .....	625 per Box	<b>330.99</b>
Pink—0.100"		
(100-6962-FPD) .....	25 per Box	<b>19.99</b>
Blue—0.125"		
(100-2195-FPD) .....	25 per Box	<b>28.79</b>
Blue—0.150"		
(100-6997-FPD) .....	25 per Box	<b>29.49</b>
Blue—0.150"		
(100-3064-FPD) .....	312 per Box	<b>333.99</b>

**SAVE** **OVER 15%**



## Block Out Compound

**HENRY SCHEIN®**

Blocks out remaining teeth and undercuts on models when making appliances. Heat resistant. Won't run or stick.

- Heat resistant: will not run or stick
- Reusable: economical

(100-3302-FPD) ..... 1 Lb. **16.99**

**SAVE** **OVER 15%**



## Insta Cool™

**Keystone Industries**

Chilling spray is ideal for instantly cooling all types of plastics and materials processed in a vacuum or pressure-forming machine, ensuring maximum thermoforming adaptation.

- Easy to use: simply point can and spray
- Eliminates costly downtime while waiting for materials to cool, thus increasing productivity and profits
- Environmentally friendly, CFC-free

(201-0168-FPD) ..... 10 oz. Can **24.99**

**SAVE** **OVER 15%**



## Vacuum-Forming Machine

**HENRY SCHEIN®**

- Top designed for better heat dispersion
- Slotted rear post for perfect position over model
- Precision vacuum adapter
- Replacement gasket fits underneath lower plate gasket

Complete Unit—110 Volt  
(100-6234-FPD) ..... Ea. **389.99**

**SAVE** **OVER 15%**



## Wheel Saw

**Keystone Industries**

Cutting wheel used for hard sheet resin and dual laminate trimming.

(255-5974-FPD) ..... Ea. **19.99**



**SAVE UP TO 20%**



## Laboratory Pumice

**HENRY SCHEIN®**

- For smoother polishing without scratching
- Consistently fast, smooth-cutting, long-lasting abrasive
- Suitable for polishing all full or partial dentures

Jar .....5 Lb. **15.99**

Specify:

Flour (100-4143-FPD) Medium (100-8129-FPD)

Fine (100-9881-FPD) Coarse (100-1454-FPD)

Medium-Fine (101-7313-FPD)

Carton .....25 Lb. **33.99**


Specify:

Flour (100-3821-FPD) Medium (100-2455-FPD)

Fine (100-1916-FPD) Coarse (100-2796-FPD)

Medium-Fine (101-9319-FPD)

**SAVE OVER \$100**



## Diamond Model Trimmer Wheels

**HENRY SCHEIN®**

- Stay sharp and outlast regular wheels 10 to 1
- Long-lasting aggressive diamond
- Aggressive wheels outlast multiple ordinary wheels
- Far quieter trimming than carborundum wheels

10"

(365-1776-FPD) .....Ea. **679.99**

12"

(365-4147-FPD) .....Ea. **799.99**

**SAVE OVER 15%**



## Muslin Buffs

**HENRY SCHEIN®**

1" x 16 Ply

(900-4227-FPD) .....12 per Pkg. **19.99**

2" x 16 Ply

(900-4229-FPD) .....12 per Pkg. **26.99**

3" x 24 Ply

(900-4224-FPD) .....12 per Pkg. **29.49**

3" x 30 Ply

(900-4221-FPD) .....12 per Pkg. **35.49**


4" x 30 Ply

(900-4220-FPD) .....12 per Pkg. **47.99**

4" x 42 Ply

(900-4219-FPD) .....12 per Pkg. **62.79**

**SAVE OVER \$150**



## Blue Wonder Diamond Wheels

**Whip Mix**

- Made of plaster, stone, or epoxy
  - Trim quickly and efficiently and are self-cleaning
  - Outperform conventional trimmer wheels.
- Suited for high-volume labs

12" (557-1218-FPD).....Ea. **849.99**

**SAVE OVER 15%**



## MT3 Klettifix

**Renfert USA**

Klettifix 80-Grit Refill Kit

(228-6617-FPD).....5 per Pkg. **86.79**

#18031200, Refill Kit Grit 120

(228-5876-FPD).....5 per Pkg. **86.79**



Denture Teeth Promotion

## Another reason to smile



Portrait® IPN



Classic®

Buy Denture Teeth

**Buy 1 - Trade 1 - Get 1**

Valid September 28<sup>th</sup> - December 27<sup>th</sup>, 2025

Promo code USA: (Q4) LQ42025

Minimum \$5,000 purchase per account (based on MSRP value).  
Limit five redemptions per account per quarter.

Offer valid from September 28<sup>th</sup> to December 27<sup>th</sup>, 2025. Promo code must be provided with order. Offers available in the U.S. only. Only one code may be applied to a single order line. All discounts taken from retail prices (MSRP) only. Not valid on previous purchases, and cannot be combined with volume discounts, accounts with special pricing, or other product/order level promotions or offers. Not redeemable for cash. No substitutions. No rain check issued. Maximum five (5) redemptions per customer. Dentsply Sirona reserves the right to cancel any promotion at any time. To the extent required by law, customers must fully and accurately disclose the amount of these discounts in any cost report or claim for reimbursement. Submitted to any applicable government healthcare program and comply with any request to provide documentation of the discounts to representatives of the Secretary of Department of Health and Human Services and State agencies. Promotion must be Redeemed by December 27<sup>th</sup>, 2025.

USE  
COUPON  
CODE

FPD

# More For Your Money & A **GIFT** Too!



## Purchase

15 Discs

## Receive

2 Discs\*\*  
1 Pink Zirconia Bur  
1 TCS apron



## Purchase

10 Packs of Catridges

## Receive

1 Pack of Catridges\*  
1 Pink Zirconia Bur  
1 TCS apron



## Purchase

A Digital Furnace

## Receive

2 Pack of Catridges\*  
1 TCS apron

\*Medium size, st. pink shade, your choice of resin: Unbreakable™, iFlex™, Karadent™

\*\*Resin and shade of choice



\*To redeem promotion, email invoice(s)  
to [sales@tcsdentalinc.com](mailto:sales@tcsdentalinc.com) within **5 days** of purchase.  
\*Limit 2 per customer. Dropships only. Offer ends 12/27/2025.

**SAVE** OVER **\$100**



## Easy Cure 2® Light-Curing Tray

**HENRY SCHEIN®**

Cures virtually all tray materials. Compact design with a digital timing control. LED lights will not generate heat so there is no cooling fan. Last curing settings are stored when powered off and then displayed during next use. Large curing chamber with height adjustment insert can accommodate 2 trays. Safety switch deactivates light when door is open. Additional features include a digital touch-pad control, variable cycle times, and polished steel interior surfaces for more even curing. 1-year warranty.

- 110 volts, 50/60 Hz
- Lamp: 3W, 4UV LED, and 2 blue LED

(570-2721-FPD) ..... Ea. **485.99**

**SAVE** OVER **15%**



## Acraweld Repair Material

**HENRY SCHEIN®**

- Features porosity-free curing
- Easy mixing, and short dough time
- Provides a versatile cure and a high-polish finish

Powder Only..... 1 Lb. **45.99**

Specify:

Light Fibered (101-3959-FPD)

Light Characterized (101-4923-FPD)

Plain Pink (101-5054-FPD)

Dark Fibered (101-3910-FPD)

#188 Original Light Characterized (101-1493-FPD)

Clear (100-9772-FPD)

Powder Only..... 5 Lb. **179.99**

Specify:

Light Fibered (101-1105-FPD)

Light Characterized (101-3466-FPD)

Plain Pink (101-3633-FPD)

#188 Original (101-6706-FPD)

Clear (101-5011-FPD)

Liquid Only—Clear (101-8333-FPD) ..... 8 oz. **26.99**

(101-7550-FPD) ..... 32 oz. **49.99**

**SAVE** OVER **15%**



## Custom Tray Material

**HENRY SCHEIN®**

- Economical prefabricated sheets to eliminate mixing and mess
- Good working time and quick cure time
- Easy trimming characteristics reduce error and effort

Trays—2mm ..... 50 per Box **69.99**


Specify:

Pink (102-7977-FPD) Blue (102-4955-FPD)

Baseplate—1.4mm

(102-8184-FPD) ..... 50 per Box **69.99**

**SAVE** OVER **15%**



## Plastic Sheets

**HENRY SCHEIN®**

4" x 4"

(100-4342-FPD) ..... 1000 per Box **39.99**

5" x 5"

(101-3759-FPD) ..... 1000 per Box **16.99**

Contains: 1000 sheets (0.001mm thickness).

**SAVE** OVER **15%**



## Heat Seal Bags

**Buffalo Dental**

Heat-sealable bags are excellent for dentures. Treated surface permits writing.

5" x 6" Dispenser Roll

(365-4920-FPD) ..... 500 per Roll **41.99**



Buy 2 Polymax,  
Get 1 x D4304 Isolafilm 1 Liter Bottle  
or 1 x D504310 Silfix 4 x 50 ml Double  
Cartridges\*

**SAVE** **OVER \$900**



### Polymax 5 Dreve Dentamid

Solid pressure-polymerization unit for cold and hot polymerization of dental resins. Suitable for 1 large articulator, 3 brass flasks, or 4 Castdon flasks. Has a 90-minute timer with a temperature range of up to 95°C that can be adjusted continuously.

Complete Unit—Red  
(475-0095-FPD) ..... Ea. **3,259.99**  
Contains: 1 of each: pressure hose & water drain hose.

\*To redeem no charge items: Purchase product(s) between 9/28/2025 - 12/27/2025. E-mail invoice or proof of purchase to Julie Halvorson at Julie.Halvorson@dreve-america.com or sales@dreve-america.com. Invoice or proof of purchase must be sent no later than 30 days after purchase date.

**SAVE \$100**



### Denture Curing Unit Keystone Industries

Boilout/curing unit. Stainless steel construction promotes easy cleaning and prevents corrosion. Thermostatically controlled for variable holding temperatures. Cures up to 4 double flask compresses (8 denture arches). 1-year warranty.

Complete Unit  
(242-7090-FPD) ..... Ea. **469.99**

**SAVE** **UP TO \$250**



### Red Wing Lathe—¼ hp Handler Manufacturing

- Standard bench lathe ideal for office or laboratory ¼ HP, 115V/60Hz

Model #26A  
Low 4½" Base  
Polishing and grinding lathe. For use on top of the laboratory bench. Shipped with standard dental, tapered shafts, each lathe includes a right- and left-chuck remover, which provides ease of removal of standard dental chuck. 115 V/110 V and 4½" shaft height.

(612-0344-FPD) ..... Ea. **479.99**  
Model #16B

Precision Quick Chuck changer on #26A lathe. Combines a #26A motor with a chuck/changer.

(612-2539-FPD) ..... Ea. **989.99**  
Contains: 1 adapter, 1¼" arbor band, 1 tapered spindle, 2 stone chucks, 2 mandrels, & 1 collet wrench.

**SAVE** **OVER \$700**



### FlexPress™ Myerson

Automated Digital Injection System Injector is capable of injecting different appliances or appliance combinations. Easily process DuraFlex™, DuraCetal®, and VisiClear® cases using the preprogrammed settings. Can program additional settings. Integrated heating unit and precision digital software create consistent injection pressure and heating profiles, which result in the uniform flow of thermoplastic materials. 3-year warranty.

Advanced Kit  
(544-0228-FPD) ..... Ea. **4,869.99**

Contains: 1 injector, 3–6 packs of each medium-size cartridge (Medium Pink DuraFlex™, DuraCetal® A2 & VisiClear®), 1 injection flask, 1–1 oz. thermoplastic model separator (TMS), 10 TMS swabs, 10 TMS trays & 1–1 oz. release compound.

**SAVE** **OVER \$550**



### Neycraft™ Duplicating Unit Ransom & Randolph, LLC

This innovative and user-friendly duplicating unit offers precise and reliable operation as well as a power-saving on/off feature for an extended motor life. Unit will run a complete cycle and return to holding temperature within 3 hours. Includes a programmable timer for delayed starts, and the ability to adjust features easily. Comes preprogrammed and ready to use.

3-Gallon Unit  
(898-0741-FPD) ..... Ea. **4,899.99**

Use Promo Code FPD. Effective through 12/27/2025.

**zahndental.com • 1-800-496-9500**

California customers please use website when ordering for Prop 65 information.



If you think this  
has a lot of  
**dental lab solutions,**  
you should see  
our website:



**zahndental.com**




Access to over 100,000  
dental lab products



Wirelessly scan products  
to be ordered



Online payments and account  
management tools



Order confirmation, shipment  
tracking, and reporting