

EQUIPMENT

PRODUCT GUIDE 2025



Your **One Source** for a **Thriving Practice**

Practices are unique and constantly evolving. **Rely on Us** for customized solutions that empower optimal patient & practice outcomes.

www.henryschein.com/medical | 1.800.772.4346 8am-8pm (et)

California customers please use website when ordering for Prop 65 information.



WELCOME TO THE HENRY SCHEIN MEDICAL EQUIPMENT GUIDE

The right equipment is the backbone of a well-functioning practice. At **Henry Schein Medical**, we’re committed to providing advanced solutions that support efficiency, comfort, and quality care.

Inside, you’ll find a comprehensive range of equipment to meet the needs of modern health care settings. From exam tables, procedure chairs, and scales to diagnostic devices, monitoring systems, and imaging equipment, we have the trusted solutions you depend on. Our portfolio also includes sterilization systems, lighting, storage solutions, and laboratory essentials—all designed to help your practice run smoothly and effectively.

Every product is backed by Henry Schein’s dedication to quality, service, and support—so you can trust that the equipment you choose will deliver lasting performance and value.

Explore the Equipment Guide and discover solutions that streamline workflows, enhance the patient experience, and keep your practice operating at its best.



■ INDEX:

AEDs	26	Electrosurgery	32–34	Patient Monitoring	23–25
Capital Laboratory	38–46	Exam Room Equipment	3–9	Practice Green	47
Carts	11	Imaging & Ultrasound	27–30	Scales	14–15
Cold Storage	36–37	Instrument Reprocessing		Scales & Body Composition	13
Diagnostic Equipment	16–22	& Sterilization	34–35	Surgical Equipment	31
DME, Mobility Aids & Transfer	12	Medical Lighting	10		



Midmark® 626 Barrier-Free® Examination Chair
MIDMARK

Many health care providers are not aware of current standards for proper blood-pressure measurement. If this isn't addressed, it leaves health systems at serious risk. 9.8% of the population is affected by overestimation of high blood-pressure errors. Multiply the cost of treatment for hypertension by the typical patient panel size and the cost of overtreatment could be as much as \$135,000 per each practicing physician. The impact of overmedication can also result in complications such as dizziness, fainting, exhaustion, or even patient falls. The Midmark 626 is the first exam chair of its kind designed to promote proper patient positioning for a more accurate blood-pressure measurement.

- Easily position patients for a better BP measurement
- Powered movement of the back section ensures the patient's back is supported
- A Barrier-Free® low chair-height allows the majority of patients, even female patients with heights in the 3rd percentile, to place their feet flat on the floor

#626-001, Base Model (133-0301)	ea
#626-002, With Heated Upholstery (133-0312)	ea
#626-003, With IQ Scale (133-0323)	ea
#626-004, With Heated Upholstery and Scale (133-0334)	ea



Midmark 647/646 Barrier-Free® Power Podiatry Treatment Chair
MIDMARK

Boasts a wider foot section for ease of use and patient comfort for all patients, with several standard features to improve efficiency and functionality in the exam room.

- Power positioning
- Power back, tilt and height
- Adjustable, removable armrest
- Foot control and touch control
- Covered glides

647 Chair (111-7263)	ea
646 Chair (111-4145)	ea
647 Chair with Programmer (111-4146)	ea
Footrest Cover (111-8514)	ea



Ritter® 224 Barrier-Free® Examination Chair
MIDMARK

The number of elderly and obese in the U.S. is steadily growing, and over one fifth of adults in America are disabled. Assisting patients with mobility limitations onto exam chairs has become a frequently performed task. Although caregivers working in health care facilities are in one of the most at-risk professions for musculoskeletal disorders, little has been done to evaluate and reduce risk. The Ritter 224 is designed to offer low and high positioning heights so you can provide the best care to all patients, whether treating the elderly, expectant mothers, or patients with disabilities.

#224-002, Without Clean Assist Roller System (133-0159)	ea
#224-003, With Clean Assist Roller System (133-0268)	ea



Ritter® 225 Barrier-Free® Examination Chair
MIDMARK

The number of elderly and obese in the U.S. is steadily growing, and over one fifth of adults in America are disabled. Assisting patients with mobility limitations onto exam chairs has become a frequently performed task. Although caregivers working in health-care facilities are in one of the most at-risk professions for musculoskeletal disorders, little has been done to evaluate and reduce risk. The Ritter 225 is designed to offer low and high positioning heights so you can provide the best care to all patients, whether treating the elderly, expectant mothers, or patients with disabilities.

#225-002, Without Clean Assist Roller System (133-0279)	ea
#225-003, With Clean Assist Roller System (133-0290)	ea



Midmark® 272 Exam Stool
MIDMARK

Provides optimal comfort, support, and stability. Star-shaped base design allows for closer access to the patient, and stable, easy maneuverability.

- Black, plastic composite, 5-leg caster base and hood
- Pneumatic seat-height adjustment

Exam Stool.....	ea
Specify:	
#272-001-845, Lunar Gray	(132-2295)
#272-001-814, Stone	(132-2291)
#272-001-855, Blue	(132-2306)
#272-001-857, Obsidian	(132-2310)
#272-001-856, Linen	(132-2307)
#272-001-851, Ore	(132-2302)
#272-001-813, Mist	(132-2289)
#272-001-844, Mineral	(132-2294)

We're reaching new heights going lower than ever before.

The ADA set the legal requirement. The US Access Board defined the standard.

We're bringing it to you. **Midmark is the first and only manufacturer in the market to have both a procedure chair and examination chair that comply with the US Access Board Standard.**

For more information, visit: henryschein.com/medical



Scan to view
patient transfer.

1 <https://www.federalregister.gov/documents/2024/07/25/2024-16266/standards-for-accessible-medical-diagnostic-equipment>

© 2025 Midmark Corporation, Versailles, Ohio USA

Getting to Know the Standard

Accessible equipment can make a crucial difference in patient outcomes and caregiver safety. Are your exam and procedure chairs compliant?

- 01 Seat Height:** 17" or lower and 25" or higher
- 02 Transfer Surface:** 21" wide and 17" deep
- 03 Base Clearance:** 26" wide or less
- 04 Transfer Supports:** <1 ½" from surface edge, capable of supporting transfer on each side, >15 ½" long, 6"- 18" high, 1 ¼" - 2" grip diameter, continuous gripping surface
- 05 Leg Supports:** Supports the patient's thigh, knee and calf





Choose different, Choose better.

From exclusive safety features to whisper-quiet mobility and intuitive ergonomics, every Brewer PLUS exam table is designed to support your team and protect your patients. Backed by the industry's only 5-year warranty, we help you deliver better care with equipment you can trust.

**BREWER**



FLEX High-Low PLUS Exam Table
BREWER

Designed to enhance patient accessibility, safety, and provider efficiency. With its advanced high-low functionality, ergonomic design, and integrated safety features, this table supports a wide range of medical examinations while ensuring patient comfort.

- Adjustable height from 18" to 37"
- Fully compliant
- Integrated ComfortGrip Patient Assist Handles offer 700-lb weight capacity to ensure stability and durability for various patient needs
- Powered backrest and footrest enhance patient comfort and optimize provider's workflow efficiency
- Integrated safety grab bars offer additional patient support when accessing or exiting the table
- Optional EasyGlide Mobility System makes room reconfigurations and cleaning easier and more effective

Table..... ea

Specify:

- | | |
|--|------------|
| #5801PLUS-PR46, Blue Fog | (147-2583) |
| #5801PLUS-MV821, Cool Gray | (147-2593) |
| #5801PLUS-AF906, Chateau Grey | (147-2592) |
| #5801PLUS-PR92, Gunmetal | (147-2595) |
| #5801PLUS-AF945, Flannel | (147-2594) |
| #3500687PLUS, EasyGlide PLUS Mobility Option | |
| Accessory | |
| (147-2249) |ea |



Access High-Low PLUS Power Exam Table
BREWER

Designed to deliver unmatched comfort, convenience, and durability, this table is the ideal choice for modern clinical environments. Engineered to meet the diverse needs of both patients and healthcare professionals with advanced features like ComfortGrip® patient assist handles and ample storage options. 3-year warranty.

Key features:

- Adjustable height from 18" to 37", fully ADA Access Board compliant
- Extendable table length from 60" to 78" with 18" SafeGlide® leg extension length
- Mayo Pass-Through Work Surface for multifunctional use
- 6 cubic feet of storage divided between a front drawer (19¾" x 19¼" x 4") and side drawers (22½" x 20¼" x 8")
- 450-lb weight capacity

- | | |
|---------------------------|----------|
| #6000PLUS-PR46, Blue Fog | |
| (147-3124) | ea |
| #6000PLUS-PR57, Clamshell | |
| (147-3138) | ea |
| #6000PLUS-PR92, Gunmetal | |
| (147-3136) | ea |

Upgrade today

10 chair models that meet Enforceable ADA Standards

Our ADA-compliant 830 Series Chairs are our most popular and versatile model, improving patient care for a variety of specialties.

Contact your local Henry Schein Medical Representative to understand how MTI helps you get compliant.



Specialized Equipment
for Specialty Healthcare
Medical Technology Industries



830 Series Procedure Chair MTI

The highly versatile 830 Series chairs are ADA compliant in all base versions (standard, swivel, or mobile) using the transfer support, providing easy patient ingress/egress, and equal access for patients with disabilities and caregivers. All 830 chairs come equipped with SmartTech™ technology to enhance safety, reduce energy consumption, and provide smoother movement, additional memory options, and more user preferences.

- 725-lb patient load capacity
- 4 programmable positions
- Programmable home and ADA-compliant home positions
- Power lift, back, tilt, and foot extension
- Hand and foot controls
- Oval articulating headrest
- 22" or 25" seat width
- 18¾"-39" seat height
- Leg rest features knee-break with 7" power foot extension
- Optional patient arms: removable, floating (compensating) arms, or seat rail or back rail mounted arms
- Removable upholstery: seamless or designer
- Tapered backrest for easier caregiver access
- Sturdy all-steel frame
- Multiple G2 accessory options

#829-001-80, Standard Base (22" seat width) (135-8121)	ea
#829-002-80, Swivel Base (22" seat width) (139-4043)	ea
#829-003-22, Mobile Base (22" seat width) (140-6171)	ea

#830-001-80, Standard Base (25" seat width) (133-8866)	ea
#830-002-80, Swivel Base (25" seat width) (136-9044)	ea
#830-003-22, Mobile Base (25" seat width) (133-9127)	ea



463 and 464 Exam/Procedure Chair MTI

The most feature-packed, ADA compliant* chair on the market. Standard features include ComfortSync™ coordinated footrest, EasyView™ chair-back mounted switches, removable slide-back floating armrests, push-button 333° swivel base and oval articulating headrest.

- 725-lb patient load capacity
- 4 programmable positions
- Programmable home and ADA-compliant home positions
- Power lift, recline, and tilt (no tilt on 463)
- 21" seat width
- 18¾"-37" seat height (2" lower than any other ENT chair on the market)
- Removable upholstery: seamless or designer
- Tapered backrest for easier caregiver access
- Sturdy, all-steel frame

#463-002-80, Swivel Base, No Tilt (136-8741)	ea
#464-002-80, Swivel Base, With Tilt (136-8744)	ea

*ADA compliant with transfer support accessory for easy patient ingress/egress



#7820 Pediatric Treatment Table
CLINTON INDUSTRIES

Treatment table with built-in digital scale features 2-door storage with 2 adjustable shelves.

- Removable center pad for easy cleaning
- Top rimmed on all 4 sides
- Tanita digital scale, with measurements in 1-oz increments
- Built-in folding infantometer with adjustable folding ends
- Built-in 14" paper dispenser and cutter
- Convenient 36" work height

43"L x 36"H x 24"W
(122-0918) ea

Quick-ship program



Medic Portable Treatment Table
IMPACT ATHLETIC

These portable treatment tables feature heavy-duty construction, onboard storage and all-terrain wheels. The Medic is perfect for a budget with a 500-lb weight capacity, while the Medic XL offers almost double the surface area and a 750-lb weight capacity.

Dimensions:
Medic: 43"L x 24"W x 33"H
Medic XL: 73"L x 24"W x 33"H

#411700, Medic
(450-0031) ea
#411701, Medic XL
(450-0032) ea



Straight Line Treatment Table*
CLINTON INDUSTRIES

- Features an all-steel frame, a laminate shelf, and an adjustable backrest
- Paper dispenser included
- #3020-27: 31"H x 27"W x 72"L
- #3020-30: 31"H x 30"W x 72"L
- Available for FedEx® Ground and UPS® shipping to save you money

#3020-27
(122-0003) ea

#3020-30
(122-0005) ea

*Available with color-coordinated laminate and vinyl. Gray shelving will be shipped unless otherwise specified.



BOSS™ Table
OAKWORKS

Each table is complete with 4 field feet and a professional carry case.

#PKG4458RB, 29", Royal Blue
(541-0044) ea
#PKG4458BK, 29", Black
(541-0045) ea
#PKG4457RB, 31", Royal Blue
(241-0113) ea
#PKG4457BK, 31", Black
(241-0114) ea

Field Feet
(625-0032) 4/set

Carry Case for 29" Table
(241-2949) ea

EXAM ROOM EQUIPMENT



Tech-Med® Upholstered Blood Drawing Chair
DUKAL

- High-quality chair with tool-free assembly and padded seat for ultimate patient comfort
- One control knob allows the arms to be raised and lowered effortlessly
- Optional flip-up arm can be adjusted to varying angles with a single knob
- Arms can be raised and lowered in infinite positions and rotate left or right
- Adjustable levelers on all 4 legs allow the chair to sit solidly on an uneven floor
- Vinyl is mildew-resistant, sulfide-stain resistant and static dissipative
- Weight capacity: 500 lb
- Overall dimensions: 38"H x 34"W x 28"D
- Chair dimensions: Seat height to floor: 19", inside arm width: 21", seat depth: 17", height of adjustable arm: 28"-35"

#4382-F
(116-8624) ea



Mayo Stand with Mobile Base
TECH-MED

- Chrome-plated frame with secure height adjustments from 28" to 48"
- Mobile base includes locking casters
- Removable stainless steel tray: 12½" x 18"
- Weight capacity: 11 lb
- Base width: 20"

#4365
(100-9541) ea



Tape Measure
TECH-MED

Linen-like. White plastic; push-button retractor.

72"L x ¼"W
(100-2267) ea



Snellon Eye Test Wall Charts
TECH-MED

Nonreflecting, matte-finish eye charts that feature brass-eyelet reinforced holes. Designed for 20-ft test distance (unless otherwise specified). Individually polybagged.

Dimensions: 22" x 11", except where noted
Snellen
(100-4836) ea
Illiterate
(591-6940) ea
Kindergarten
(591-7664) ea
Snellen, Plastic
(104-3686) ea
Pocket Size, 6¾" x 3½"
(103-0532) ea

EXPLORE OUR ROBUST
COLLECTION OF
**PRODUCT
GUIDES**
ONLINE NOW!



www.henryschein.com/productguides



Welch Allyn Green Series™ 300 General Exam Lights BAXTER

Utilizing a single white LED, the GS 300 provides a whiter, brighter light than halogen and consumes less energy. It's the perfect combination of intense light output and cool operation in a compact design. Features touchless on/off, multiple mounting options for floor space maximization, low cost of ownership due to elimination of bulb replacement and lower energy consumption, and 5-year limited warranty to include the LED light source and power supply.

Light with Mobile Stand (114-2692)	ea
Light with Table/Wall Mount (114-5632)	ea



Welch Allyn Green Series™ 600 Minor-Procedure Lights BAXTER

Featuring 3 LEDs and a wider head, the GS 600 is designed specifically for demanding viewing conditions. Whiter, brighter, and greener: 3 white LEDs offer greater intensity. Other features include touchless power controls and an adjustable light output. It also has a disposable sheath available to help reduce the risk of infection due to cross contamination. Rugged, ergonomic design allows for improved performance and long-lasting durability. Low cost of ownership due to elimination of bulb replacement and lower energy consumption. 5-year limited warranty to include the LED light source and power supply only.

#44600, With Mobile Stand (566-0088)	ea
#44610, With Table/Wall Mount (566-0089)	ea




Clarus 7™ LED Exam Light MTI

MTI's brightest exam light uses 7 LEDs and has a slightly larger 3½" light head, providing maximum illumination and operational ease. Can be used with either a caster base or floor stand. Physical dimensions: head diameter: 3½"; caster base height: 73½"; floor stand height: 71½"; caster base diameter: 23¾", floor stand diameter: 15". 2-year limited warranty.

- LED life: 40,000 hours
- Color temperature (CCT): 4,300° K
- Light intensity @ 16": 33,400 lux
- Light intensity @ 12": 48,500 lux
- Working distance: 16"
- Spot size at 10": 5"
- Spot size at 12": 6"
- Spot size at 16": 8"
- Spot size at 24": 11"
- CRI: 75
- Control box is dimmable on light head
- Electrical input: 100–240 V, 50–60 Hz, 20W, 12V @ 1.67A
- 8-ft power cord: 3-wire grounded cord set
- Listings: UL 153, UL 8750, CAN/CSA: C22.2, No. 12

Light with Caster Base (121-8211)	ea
Light with Floor Stand (121-8212)	ea



Welch Allyn Green Series™ Exam Light IV BAXTER

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas, the GS Exam Light IV provides you with 3 times the typical light output and features a highly flexible light pipe that is twice as long as conventional pipes. And with a compact, unobtrusive head, it ensures pinpoint accuracy and provides true tissue color rendition during an exam. Utilizing white LEDs, the GS Exam Light IV provides a whiter, brighter light than halogen and consumes less energy. 5-year limited warranty to include the LED light source and power supply only.

- One LED
- Flexible pipe
- Touchless on/off
- Easy to assemble

Light with Table/Wall Mount (114-2547)	ea
Mobile Stand For great office mobility and space usage. Attaches easily to unit. (114-4177)	ea
Disposable Light Sheaths (566-0989)	25/box



Ritter 250-003 LED Exam Light MIDMARK

Features LED technology for longer life, brighter intensity and reduced energy consumption. Provides up to 25,000 hr of lighting—more than 12 times longer than a traditional halogen exam-room light. Light head size: 4½" diam. x 4" long. Height: 36" long, flex length: 22".

- Light source: 12V, 3 LEDs, MR16 type (with glass cover). Illumination @ 18": 1,800 fc. Color temperature: 4,000°K. Life: 25,000 hours
- Electrical: 115V AC, 50–60 Hz, 0.12A. Power cord: 8 ft hospital-grade plug
- Mechanical: Light head size: 4"L x 4½" diameter
- Flex arm: 22"L x 0.585" diameter. Neck material: high carbon steel covered in vinyl
- Base tube: 36"L. Outer tube diameter: 1". Tube material: steel
- Standard base housing: 2¼"H x 13" diameter
- Weight: 12 lb
- Color: Pebble Gray

#250-003, LED Exam Light (132-3839)	ea
--	----



Emergency Rescue Cart DETECTO

Detecto's Rescue series medical carts are manufactured with a quality steel cart body and feature a powder-coat painted finish for enhanced durability and ABS countertops that are easy to clean. The carts are equipped with either a keyed or EMG breakaway central lock that secures all drawers simultaneously for added convenience. Available with either 5 or 6 drawers.

Dimensions:

With guide handle: 43.2"H x 29.9"W x 21.3"D

Without guide handle: 41.8"H x 27.8"W x 21.3"D

Key features:

- 2 parking locks and 1 steering lock on wheels
- 360° rotation for maneuvering tight corners
- ABS bumpers on all 4 corners
- ABS countertop surface
- Auto-return rail on drawers
- Available in loaded models stocked with accessories

#RC33669RED, 5 Drawer
(129-3899) ea

#RC333369RED, 6 Drawer
(129-1293) ea

#RC33669RED-L, Loaded, 5 Drawer
(131-1496) ea

#RC333369RED-L, Loaded, 6 Drawer
(146-5166) ea



Four-Drawer Crash Cart with Specialty Accessory Package

HARLOFF

- Breakaway lock bar
- 100 breakaway seals included
- 24" of vertical drawer space
- Stabilizing bumper frame with extended front casters and 4 corner bumpers
- 5" ball-bearing casters, all full swivel, 2 with brake
- Predrilled threaded holes for V-Series side-mounted accessories
- Ball-bearing, full-extension drawer slides
- Ships with wheels installed and ready to use
- Minimal assembly required at delivery
- Not made with natural rubber latex
- Durable, hard-baked, red powder-coat finish
- Steel construction
- Emergency package accessories include: Cardiac board and brackets, side-mounted drop shelf, IV pole, oxygen tank holder (D/E size cylinders), and powder-coated steel top rail/push handle
- 12-year warranty

#V24-4B-HSACC
(738-0008) ea



Standard-Duty Utility Carts— Stainless Steel

LAKESIDE MANUFACTURING

Made from stainless steel, they feature 3 stain-resistant shelves, 3½" casters, and a load capacity of 300 lb.

- Durable, unitized electronically welded stainless steel angle leg design
- Reinforced, hemmed front on 22-ga stainless steel shelves
- Legs are 20-ga stainless steel
- Sound-deadening panels under shelves restrict vibrations and noise
- Bumpers on legs and handles
- Cleans and sanitizes easily
- 300-lb capacity

33"H x 18¾"W x 30¾"L
(522-7623) ea

32⅞"H x 16¼"W x 27½"L
(522-0083) ea

DETECTO

DRIVING BETTER CARE, ONE CART AT A TIME.

Engineered for efficiency and safety, our medical carts deliver secure storage, customizable configurations, and smooth mobility — helping clinicians save time and focus on patient care.



RESCUE
medical carts



MOBILECARE
medical carts



PEDIATRIC
scale cart

Rescue Carts

- Rapid Delivery** – Offers quick-ship availability
- Vibrant Selection** – Multiple color options (red, blue, yellow, white)
- Fully Loaded** – Accessories pre-stocked in "loaded" models
- Compliance Certified** – TAA-approved for trusted procurement

Mobile Care Carts

- Smart Control** – Software that simplifies cart management
- Certified Trust** – TAA-approved for compliance confidence
- Secure Tracking** – Audit trails deter narcotic theft
- Total Protection** – Advanced locks and instant alerts

Pediatric Scale Carts

- Perfect Height** – Ergonomic design eases infant weigh-ins
- Built-In Storage** – Keep diapers and supplies within reach
- Custom Fit** – Designed for DETECTO 6745 pediatric scale
- Safe & Mobile** – Protective lip and smooth-rolling casters



BridgeAir™ Single-Patient-Use Lateral Transfer Air Mattresses

BRIDGE HEALTHCARE BY PELSTAR

Air-assisted breathable lateral transfer and repositioning mattress that significantly reduces the physical effort required to move patients while minimizing patient discomfort.

- Minimize staff and patient injuries related to lateral transfer and positioning
- Durable, lightweight, and safe to leave under patient when not in use
- Easy-to-clean surfaces to support hospital hygiene protocols and infection prevention
- Also available in half sizes, reusable sizes, and with optional bed connecting straps

Single Patient Use, Full

#PS34SPUL, 34" x 78"	10/case
(146-8028)	
#PS39SPUL, 39" x 78"	10/case
(146-8027)	
#PS50SPUL, 50" x 78"	5/case
(146-8032)	



SlideMate™ Disposable, Single-Patient Use and Reusable SlideTubes

BRIDGE HEALTHCARE BY PELSTAR

A patient transfer and repositioning solution designed to reduce caregiver strain and enhance patient safety.

- Low-friction slide sheet reduces patient transfer force by up to 50%
- Easy to use and cost effective to prevent staff injuries
- Durable single-patient-use quality to last the patient's entire length of stay

Single Patient Use

#SMP-33SPU, 40" x 33" SlideMate™ Pocket	100/case
(145-4030)	
#SMP-29SPU, Original, 58" x 29"	50/case
(143-9003)	
#SMP-48SPU, Bari, 58" x 48"	50/case
(145-4028)	
Reusable	
#SMBH-5829R, Original, 58" x 29"	ea
(145-4031)	
#SMBH-5848R, Bari, 58" x 48"	ea
(145-4032)	
#SM-DISP, Wall-Mounted Dispenser	10/case
(145-4033)	



Sport Wheelchair, 18" with Fixed Arms

DRIVE MEDICAL

- Powder-coated silver-vein steel frame
- Embossed vinyl upholstery is durable, lightweight, attractive, and easy to clean
- Urethane tires mounted on composite wheels provide durability and low maintenance
- 8" front casters, adjustable in 3 positions
- Padded armrests provide additional comfort
- Carry pocket on backrest provides additional convenience
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable-length riggings
- Plastic foot plates with heel loops on footrests
- Comes with push-to-lock wheel locks

#SSP218FA-SF	
(116-8651)	ea



Bariatric Sentra EC Heavy-Duty Wheelchair

DRIVE MEDICAL

Crafted with a robust, carbon-steel frame with powder-coated, silver-vein finish for enduring appeal and slip-resistant protection. Its double-embossed vinyl upholstery offers a blend of durability, lightweight comfort, and easy maintenance. Transition effortlessly with the standard dual axle, while the urethane tires on composite wheels ensure long-lasting performance with minimal upkeep. Equipped with 8" front casters adjustable in 3 positions, and a choice of swing-away footrests or elevating leg rests with tool-free adjustable length riggings, the Sentra EC offers versatility and convenience tailored to your needs.

- Weight capacity: 450 lb
- Maintenance free, sealed, precision bearings
- Padded armrests provide additional comfort
- Aluminum foot plates
- Comes with push-to-lock wheel locks

#STD24ECDDA-ELR	
(116-8252)	ea



Steerable Knee Walker with Basket

DRIVE MEDICAL

Model #796

- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation, and ulcers of the foot
- New, 2-piece padded leg pad provides comfort
- Deluxe braking system
- Removable front basket
- Durable, steel frame with silver vein finish
- Tool-free height adjustment
- 8" casters are ideal for indoor/outdoor use

(124-2066)	ea
------------	----



Wheelchair, 22" With Removable Arms

DRIVE MEDICAL

- Weight capacity: 450 lb
- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Double-embossed vinyl upholstery is durable, lightweight, attractive, and easy to clean
- Standard dual axle provides easy transition of seat to hemi-level
- Urethane tires mounted on composite wheels provide durability and low maintenance
- Chrome hand rims are standard
- 8" front casters, adjustable in 3 positions
- No-maintenance, sealed, precision bearings
- Padded armrests provide additional comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable-length riggings
- Aluminum foot plates
- Comes with push-to-lock wheel locks


#STD22ECDDA-SF	
(115-5651)	ea

HELPING YOU MAKE BETTER HEALTH DECISIONS

—With Insight, Efficiency, and Compliance

Meet the Alpha

NEW!
mBCA Alpha
Body Composition
Analyzer
HS Item #148-9987



NEW!
XLine Hold
ADA Compliant Bariatric
Wheelchair Scale
HS Item # 149-2054



Scale-up Line
ADA Compliant Mobile
Handrail Scale
HS Item # 145-9845

seca 777 Physician Scale

SECA

With Eye-Level Display

Combines precision, durability, and safety for high-traffic clinical environments. Featuring a built-in measuring rod with a 4"–90" range, it ensures accurate height and weight readings for a wide range of patients. Designed with an oversized, low-profile, full-metal platform and a tilt-proof safety design, the seca 777 provides stability and reliability without plastic components. Its USB interface enables seamless PC connectivity for efficient data integration. Backed by an 8-year warranty—4 times the industry standard—and compatible with the Panda Bear Bao Bao accessory, this scale is built for long-lasting performance in busy healthcare settings.

- Built-in measuring rod with 4"–90" measuring range
- Oversized, low-profile, full-metal platform (no plastic) with tilt-proof safety design
- 8-year warranty (4 times the industry standard)
- USB interface for PC connectivity
- Compatible with Panda Bear Bao Bao accessory
- Capacity: 550 lb/250 kg
- Graduation: 0.2 lb/100 gm

#7771821004, Scale
(139-2675) ea

#4770000009, Optional Integrated Handrail
Accessory for seca 777, 787 and 797
Physician Scales
(138-4356) ea

#4590000009, Panda Bear Accessory for seca
Physician Scales, Scale-up Line, and Body
Composition Analyzers
(138-5120) ea



seca 554 Medical Body Composition Analyzer SECA

An effective training program needs validated data to determine a baseline, set attainable goals, and motivate athletes with valuable insights. seca BIA solutions provide medical and sports professionals precise data into body fat percentage, skeletal muscle mass, body water and many more parameters by using bioimpedance analysis (BIA). Their BIA assessment is a noninvasive, 24-second analysis that provides insights from the inside out. Provide athletes with real data and track progress to support their training and fitness goals.

- seca's enhanced BIA technology is the result of a long-term collaboration with medical experts and researchers in the field of body composition
- Validated as the only BIA solution for all components with the most accurate reference method (gold standard). This includes the complex comparison with a whole-body MRI for muscle mass and the 4C model for fat mass or fat-free mass
- Rapid 24-second assessment of vital body composition parameters, e.g., percentage of body fat, visceral adipose tissue, skeletal muscle mass, and body water
- Educate athletes with custom visuals that show trend analysis and progress over time
- Integrated handrail ensures reliable, repeatable results
- Capacity: 800 lb/300 kg

#ONIHMINNTN, Analyzer
(138-8511) ea

#ONIHMIUTNN, Analyzer with Ultrasonic Height
Measurement
(139-5321) ea



Note: Monthly cloud-based software fee applies.

Meet Accessibility Standards with *Certified* ADA Compliant Scales

The **only** ADA compliant scales in the market that are independently certified by a third party, supporting healthcare facilities in providing accessible, compliant, and inclusive care.

2101KL Digital Platform Scale

HS Item # 120-6338



1100KL Digital Platform Scale with Extra-wide Handrails

HS Item # 111-5112



2620KL Digital Wheelchair Scale

HS Item # 148-4830



Health o meter®
Professional

For more information, contact your local Henry Schein Representative, or visit henryschein.com/medical


Eye-Level Digital Physician Scale with Digital Height Rod **HEALTH O METER**

This high-capacity 660 lb/300 kg scale features a market-leading built-in digital height rod and the most connectivity options in its class.



- Functions: lb/kg conversion, lb/kg lockout, BMI, zero, tare, pre-tare, hold/release, auto zero, auto off, reweigh, recall, and Everlock®
- Low-profile nonskid platform with reinforced base for added stability
- State-of-the-art digital height rod
- Optional 500HB handlebars (137-2640) to ensure patient safety
- Wireless connectivity options available

#502KL
(145-6111)ea



Digital Bariatric Platform Scales Pro Plus™ **HEALTH O METER®**

- Capacity 700 lb/310 kg
 - Platform: 15¾" x 22"
 - ¾" LCD with 180° swivel movement
 - Functions: BMI, EMR capable, Hold/Release button, adjustable backlight, Auto-Off, Auto-Zero, LB KG conversion, Zero Out/Tare, date/time display, updatable software, modular design, smart digital load cells
 - 120 V adapter included; connectivity through USB
- #1100 KL
(111-5112)ea



Eye-Level Digital Physician Scale **HEALTH-O-METER**

High-capacity 660 lb/300 kg primary-care scale provides the optimum in quality, reliability, and ease of use.

- Functions: lb/kg conversion, lb/kg lockout, BMI, zero, tare, pre-tare, hold/release, auto zero, auto off, reweigh, recall, and Everlock®
- Low profile, nonskid platform with reinforced base for added stability
- Integrated height rod with 30"-84¼" measuring range
- Wireless connectivity options available

#500KL, Scale
(106-6626)ea
#500HB, Handlebar Accessory
(137-2640)ea



Digital Pediatric Scale with Extra-Wide Tray **HEALTH-O-METER**

High-performance digital pediatric scale with motion-sensing technology to capture the weights of even the most active babies.

- Capacity 44 lb/20 kg
- EMR connectivity
- Functions: Auto-Off, Auto-Zero, LB/KG conversion, and Zero-Out/Tare
- Additional features for #553KL: measuring tape, 1" LCD readout, EMR capable, and Hold/Release button
- Includes 6 AA batteries (#553KL); one 9V (#549KL)

#549KL, Scale
(110-0989)ea
#553KL, Scale with EMR
(104-7675)ea



Digital Electronic Eye-Height Scale

HENRY SCHEIN®

Compact, eye-height digital scale that features an integrated mechanical height rod, large backlit LCD display, and BMI calculator. Weight capacity up to 500 lb (225 kg). Low profile for easy access—platform is only 2½" high. Sta-Bilt stabilizing software ensures an accurate measure, even if the patient is moving. Autohold feature allows the weight to continue to be displayed after the patient has stepped off the scale. Integrated USB interface allows data to be transferred to a computer. Includes 2-year warranty.

22" x 15" x 53"

(112-7154) ea




Mechanical Beam Scale with Height Rod

HENRY SCHEIN®

Standing mechanical beam scale with built-in height rod. Weight capacity of 450 lb (200 kg). 1-year warranty.

22¼" x 18¼" x 58¼"

(112-7157) ea



Digital Wheelchair Scale

HENRY SCHEIN®

- Single, stable platform allows scale access with wheelchair
- Easily folded to minimize footprint and moved with casters
- Individual wheelchair weight can be programmed for easy recall
- 32" x 33" x 2½", 1000-lb maximum capacity
- 2-year limited warranty

(570-0597) ea

Digital Waist-High Scale

HENRY SCHEIN®

Compact, waist-high digital scale that features a large, backlit LCD display and shows lb or kg. Capacity up to 500 lb (225 kg). Low profile for easy access—platform is only 2½" high. Wheels for easy movement. Contains BMI calculator and clock. Includes 2-year warranty.

22" x 15" x 36"

(112-7155) ea



Digital Pediatric Table Scale

HENRY SCHEIN®

Small countertop pediatric scale with a curved design to support and protect the infant. Sta-Bilt stabilizing software ensures an accurate measure, even if the infant is moving. Large backlit LCD display measures in lb and oz. Weight capacity of 45 lb. Integrated tape measurer to determine the baby's size. Autohold feature allows the weight to be displayed, even after the baby has been taken off the scale. Integrated USB interface allows data to be transferred to a computer. Includes a 2-year warranty.

16¾" x 27¾" x 4"

(112-7156) ea



Digital Hand Rail Scale Height Rod


HENRY SCHEIN®

- Handrails allow patients to hold steady during weighing
- Low-profile design for easy access, platform is only 2¾" high
- Contains BMI calculator and clock
- 28" x 23½" x 2¾", 1000-lb maximum capacity
- 2-year limited warranty

(570-0598) ea



3M™ Red Dot™ Monitoring Electrodes
3M MEDICAL PRODUCTS
These multipurpose monitoring electrodes with foam backing and sticky gel are suited for all applications.
#2560, Foam, Diaphoretic, 1.6" x 1.38"
(4 cm x 3.5 cm)
(777-9475) 50/bag



Diagnostic Tab Electrodes
HENRY SCHEIN®
• Conductive adhesive hydrogel for excellent adhesion and clear, reliable tracings
• Silver/silver chloride (Ag/AgCl) sensing system is defibrillation-recoverable
• Not made with natural rubber latex
• DEHP-free, mercury-free
Electrode Resting Tabs
(570-1173) 100/pkg



Adult ECG Electrodes
CONMED
Clear Neoderm® tape allows visual monitoring of the patient's skin condition at all times. Perforated tape permits skin ventilation and moisture evaporation for increased patient comfort. Conforms to patient's body without pulling on the skin. Will not absorb fluids. Foam electrodes are ideal for general-purpose monitoring. Closed-cell foam keeps fluids from penetrating. Pull tab makes it easy to apply and remove the electrode.
Cleartrace™—Tape
(147-5056) 3 Pouches per Pkg.
Suretrace®—Foam
3 per Pouch
(147-5554) 3 Pouches per Pkg.
#1800-030
(520-6007) 30 Pouches per Pkg.
#1870-030
(115-2647) 30 Pouches per Pkg.



Defibrillator Pads
HENRY SCHEIN®
• Adult/pediatric multifunctional electrodes
• Peel-pouch
• 60" lead wire
• 3-year shelf life
• Not made with natural rubber latex
Defibrillator Pads pair
Specify:
Adult, Zoll (112-7167)
Adult, Physio (112-7163)
Pediatric, Zoll (112-7164)
Pediatric, Physio (112-7165)
Pediatric, Philips (112-7168)



Cleartrace™ and Cleartrace™
RTL Electrodes
CONMED
The combination of the low-profile perforated tape, GEL 2000X, and the round shape and size make the Cleartrace™ and Cleartrace™ RTL electrodes an ideal solution for general-purpose monitoring.
Features:
• Backing is made of a clear NeoDerm® tape, which allows visual monitoring of the patient's skin condition at all times
• Perforated tape permits skin ventilation and moisture evaporation for increased patient comfort
• Conforms to patient's body without pulling on their skin and will not absorb fluid
Cleartrace Adult Tape ECG
#1700-030
(681-3881) 30/pkg
#1700-050
(121-7460) 1000/case
#1700C-003, Cleartrace RTL Adult Tape ECG,
3/pouch
(115-4728) 600/case



Welch Allyn RetinaVue 700 Imager BAXTER

As the latest addition to the RetinaVue care delivery model family, the RetinaVue 700 Imager makes retinal exams simple and affordable for primary-care applications. This easy-to-use smart camera offers an automated retina imaging experience, featuring auto-alignment, auto-focus and auto-capture. It provides the user with flexible options, including EMR connectivity and WIFI for efficient Information transfer and seamless integration with existing clinical workflows.

- 2.5-mm small-pupil capability requires little-to-no chemical dilation
- 60° field of view
- HIPAA/HITECH-compliant RetinaVue Network software has end-to-end data encryption
- RetinaVue Network software is hosted on the Microsoft Azure cloud
- Small, handheld design is perfect for performing retinal exams in busy clinics and mobile health-care applications
- Can help improve HEDIS, Star, and HCC risk factors by closing gaps for patients living with diabetes

#RV700-B
(136-5631) ea



Welch Allyn® Universal Desk Charger BAXTER

- Integrated light indicators highlight that power is flowing to the unit and that handles are charging
- Designed for Welch Allyn Lithium Ion, Lithium Ion Plus battery handles, as well as 3.5 V or 2.5 V NiCad, Lithium Ion, and PocketScope® handles
- Optional wall mount kit to affix on wall or cabinet
- Now featuring iExaminer® Pro SmartBracket storage to help the user quickly apply to the MacroView Plus otoscope to capture and share images of the ear

#719-DSK, Universal Desk Charger
(140-1794) ea

#71-XXXLDE, Universal Desk Charger 3.5 V
(140-1796) ea

#71-PM3LDE-US, Desk Charger Set, PanOptic/Mac
(140-2385) ea



Connex® Spot Monitor BAXTER

Features an easy-to-use, vivid touch-screen display and provides accurate vital signs measurement, including blood-pressure averaging, spot checking, interval monitoring and custom scoring across patient populations. Device connects wirelessly to an EMR and accurately sends vitals to the patient chart from the point of care to help improve efficiency.

- Designed for the physician-office environment with an easy-to-use, vivid touch-screen display
- More accurate hypertension diagnosis without slowing down your workflow thanks to automated blood-pressure averaging
- Improved capabilities for assessing respiratory conditions such as COPD, asthma, and community-acquired pneumonia (CAP), thanks to integrated professional-grade pulse oximetry
- Increase staff efficiency and eliminate data-entry errors with wireless transfer of vital signs directly into an EMR

7100 Base Models with Wired Connectivity
#71XX-B, SureBP Non-invasive Blood Pressure
(566-0464) ea

#71XT-B, SureBP, SureTemp
(672-0176) ea

#71XE-B, SureBP, Braun Thermoscan PRO 6000 Ear Thermometer
(124-6103) ea

#71WX-B, SureBP, Nonin SpO₂
(672-0177) ea

#71WT-B, SureBP, Nonin SpO₂, SureTemp
(566-0466) ea

#71WE-B, SureBP, Nonin SpO₂, Braun Pro6000
(124-5328) ea

7300 Bluetooth Connectivity Physician Office Models
#73XT-B, SureBP, SureTemp
(566-0465) ea

#73XE-B, SureBP, Braun PRO 6000
(124-8200) ea

#73WX-B, SureBP, Nonin SpO₂
(672-0178) ea

#73WT-B, SurePB, Nonin SpO₂, SureTemp
(672-0179) ea

#73WE-B, SureBP, Nonin SpO₂, Braun PRO 6000
(124-4801) ea

#73CT-B, SureBP, Covidien SpO₂, SureTemp
(672-0175) ea

#73CE-B, SureBP, Covidien SpO₂, Braun PRO 6000
(125-8900) ea

#73MT-B, SureBP, Masimo SpO₂, SureTemp
(566-0467) ea

#73ME-B SureBP, Masimo SpO₂, Braun PRO 6000
(124-7063) ea



Welch Allyn Green Series™ 777 Integrated Wall System BAXTER

#777-PM2WAS-US, With PanOptic Basic LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-1989) ea

#777-PM3WAS-US, With PanOptic Plus LED Ophthalmoscope, MacroView Plus LED Otoscope for iExaminer, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-2495) ea

#777-SM2WAS, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-1801) ea

#777-SM2WAX, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, and Ear Specula Dispenser
(140-1799) ea



Welch Allyn Spot Vital Signs® 4400 Device BAXTER

With SureBP® Noninvasive Blood Pressure and SureTemp Plus Thermometer

- SureBP® technology can record blood pressure in just 15 seconds
- Blood-pressure averaging helps ensure accurate hypertension detection
- SureTemp® Plus oral/axillary thermometry
- Option for Nonin® SpO₂ technology
- 7" color touchscreen for a simplified workflow
- Options for documenting patient observations (like respiration rate and body mass index)
- EMR connectivity helps the care team avoid documentation errors
- Ships with one adult blood pressure cuff

#44XT-B
(137-3839) ea



YOU'RE ON THE FRONT LINE OF POPULATION HEALTH.

In a complex healthcare system, community health centers represent the first line of defense against the health issues patients face today. Helping to manage chronic conditions and identify signs of disease early are important.

We provide smart tools to help you identify, diagnose and manage the prevalent health issues you see every day.



Diabetes



Hypertension



Heart
Disease



Vision Loss

Scan here to learn more about the front line care portfolio.
Please contact your Henry Schein representative to learn more.

Rx only. For safe and proper use of the products mentioned herein,
please refer to the appropriate Operator's Manual or Instructions for Use.
Spot Vision Screener is not a prescriptive device.

Baxter, Connex, CP 150, ELI, Diagnostic Cardiology Suite, RetinaVue, Spot, Spot Vital Signs,
and Welch Allyn are trademarks of Baxter International Inc or its subsidiaries.

US-FLC45-250008 [v1.0] 01/2025





ELI 280 RESTING ECG

BUR280-C1X,
BUR280-81X,
BUR280-W1X,
ELI280-AAA-AAFBD,
ELI280-CAA-AAFBD,
ELI280-CAA-AAFBT,
ELI280-AAA-AAFBT



CP 150 RESTING ECG

CP150A-1ENB,
CP150AS-1ENB



DIAGNOSTIC CARDIOLOGY SUITE

CC-RXX-WAXX, CC-RXX-WADX,
CC-RXX-AAXX, CC-RXX-AADX,
DCSS-CXX, DCSS-NXX



WELCH ALLYN SPOT VITAL SIGNS 4400 DEVICE

44WT-B, 44XT-B



WELCH ALLYN CONNEX SPOT MONITOR

ALL MODELS



WELCH ALLYN SPOT VISION SCREENER

VS100S-B



WELCH ALLYN RETINAVUE 700 IMAGER

RV700-B

EasyOne® Air Spirometer
HENRY SCHEIN®

The all-in-one revolution in lung-function testing. Built for health-care providers large and small, EasyOne Air puts the power of the most advanced spirometry testing tools in the palm of your hand. Revolutionary TrueFlow ultrasound technology provides unprecedented accuracy and reliability. Equipped with ready-to-use EMR connectivity and a user-friendly interface, EasyOne Air makes the best in lung care available anywhere.

- Portable and PC-based solution
- Large, color touch screen
- Wireless connectivity to the PC via Bluetooth
- Instant interpretation based on ATS/ERS guidelines
- Rechargeable battery
- Pediatric incentive screens available with EasyOne Connect
- Calibration-free
- Compatible with EasyOne Connect

#2500-2HS, Spirometer (570-0683)	ea
#2050-1HS, Spirette Mouthpieces (570-0383)	50/case
#2050-5HS, Spirette Mouthpieces (570-0384)	200/case



Digital Ear Thermometer
HENRY SCHEIN®

- For professional or personal use
- Instant reading results in 1 second
- Memory feature recalls last reading
- Fever alarm
- Low battery indicator
- Replaceable lithium battery
- °F or °C measuring mode
- Thermometer includes 15 disposable probe covers

Thermometer (900-6438)	ea
Probe Covers (900-6439)	150/box

Pro Pocket Aneroid Sphygmomanometer
HENRY SCHEIN®

A value-priced aneroid that meets the needs of the most demanding professionals. Precision-crafted movement is rated to 30,000 cycles and backed by a 20-year calibration warranty. Extended graduation marks eliminate parallax error. White numbers on contrasting black scale with red indicator needle for easier viewing. Proprietary Size Guide-Marking System to enhance measurement ease and prevent miscuffing. Chrome-plated brass deflation valve. Leatherette zippered carrying case. Inspected, assembled, and packaged in the USA. Not made with natural rubber latex.



Infant (112-6073)	ea
Child (112-6071)	ea
Adult (112-6072)	ea
Large Adult (112-6075)	ea
Thigh (112-6074)	ea



Fingertip Pulse Oximeter
HENRY SCHEIN®

- Dual-color LED makes monitoring easy
- Display modes include SpO₂ (range 70%–100%) and pulse rate (range 30–235 bpm) along with pulse bar and plethysmogram
- 2 AAA batteries provide 30 hours of continuous use
- Visual and audible low-battery alarm
- Suitable for both adult and pediatric patients
- Electromagnetic compatibility

(570-3222)	ea
------------	---------

Green Series™ 777 Integrated Wall System
HENRY SCHEIN®

All models include a BP aneroid, an ear specula dispenser, and a SureTemp Plus thermometer.

- #777-PM2WAS-HS, With PanOptic Basic LED Ophthalmoscope and MacroView Basic LED Otoscope
(570-2666).....ea
- #777-PM3WAS-HS, With PanOptic Plus LED Ophthalmoscope and MacroView Plus LED Otoscope for iExaminer
(570-2665).....ea
- #777-SM2WAS-HS, With Coaxial LED Ophthalmoscope and MacroView Basic LED Otoscope
(570-2667).....ea
- #777-SM2WAX-HS, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, and Ear Specula Dispenser (Does not include SureTemp Plus Thermometer.)
(570-2668).....ea



MEDSOURCELABS

End-Tidal CO₂ Capnography

Breath-by-breath sampling of ETCO₂ while simultaneously delivering oxygen



Non-Contact Infrared Thermometer #IR 300 MEDSOURCE INTERNATIONAL

Provides instant temperature readings without worry of physical contact or cross contamination. Measuring distance of 0.4"–2". Temperature displays in <2 seconds. Color digital screen changes color to indicate normal, medium, or high temperature. Automatic power-off in <18 seconds. Memory function records temperature settings for up to 50 groups on 2 AAA batteries. Batteries included. 1-year device warranty.

Model #IR 300
(137-8300) ea



Sprague Stethoscope MEDSOURCE INTERNATIONAL

Features a double-tube configuration that provides improved acoustics. This set includes an adult and pediatric diaphragm and 3 bells (1", 1¼", and 1½") for accurate detection of high, middle, and low frequency sounds in adults and children. This stethoscope is 22" long and has an adjustable, chrome-plated brass binaural.

Stethoscope.....ea

Specify:
#MS-STBK, Black (499-2487)
#MS-STBLU, Blue (499-5138)
#MS-STG, Green (499-5139)
#MS-STR, Red (499-5140)



Fingertip Pulse Oximeter MEDSOURCE INTERNATIONAL

Offers 6 display modes, including SpO₂, Pulse Rate, and Pleth Waveform. This easy-to-read, 2-color OLED display makes monitoring clear and simple. The automatic switch-off feature preserves battery life. Suitable for adults and pediatric patients.

#MS-74002
(499-4744)ea



Blood Pressure Unit MEDSOURCE INTERNATIONAL

This high-quality BP unit features an accurate, no-needle stop gauge, a durable, easy-to-adjust nylon cuff with hook-and-loop fasteners and a standard pressure-release valve. It was created for the cost-conscious care provider who is not willing to compromise on accuracy and quality. Each unit is sold complete with matching carrying case.

- Range: 0–300 mmHg
- Durable nylon cuff
- Carrying case included
- Accurate, no-needle stop gauge

Blood Pressure Unit.....ea

Specify:
Infant (499-2556) Large Adult (499-2489)
Pediatric (499-2555) Thigh (499-5520)
Adult (499-2488)

The *e-sphyg*³ Digital Blood Pressure Monitor

Designed for busy professionals who value both **accuracy and affordability**, the e-sphyg 3 is the feature-packed, clinical-grade NIBP solution you've been waiting for.

- Systolic blood pressure
- Mean arterial pressure
- Diastolic blood pressure
- Irregular heartbeat
- Pulse Rate

1x **Single measurement mode**
with pulse rate and mean arterial pressure (MAP) displayed in as little as 30 seconds.

3x **Customizable, 3-reading average mode**
with individual measurement display

Manual mode
for auscultation with a stethoscope

High-capacity Rechargeable battery
STANDARD



USB RECHARGEABLE 3.5V Portable EENT Instruments

- C sized handle with bayonet mount for easy head changes
- USB rechargeable 18650 Li-ion battery, 2600 mAh, 9.62Wh
- CNC machined from solid brass and chrome plated in black
- Adtronic electronic rheostat for infinite control
- 18-inch universal USB type C charging cord included



10 hours lighting time when paired with LED lamps • 3.5 hours lighting time when paired with Xenon lamps



ADView 2 Diagnostic Station AMERICAN DIAGNOSTIC CORPORATION

The ADView 2 is the only truly modular diagnostic station that grows with your needs. Start with blood pressure and choose temperature or pulse oximetry options at the time of purchase; we'll ship your ADView 2 fully configured and ready to go. Or add new modules and accessories in the field, anytime, as you need them. Our user-expandable, connect-and-go construction makes upgrading effortless.

- Clinical-grade Advantage™ BP technology from SunTech Medical
- Contemporary, compact design with an integrated handle and easy-to-read display
- Manual blood pressure mode for real-time verification
- Customizable average mode options
- Next-generation EMR-ready with HL7 messaging for electronic data transfer
- Network-enabled with built-in Ethernet or optional Wi-Fi
- Lithium-ion power for longer-lasting battery life and a lighter weight
- Built-in memory with up to 99 charting sessions
- New Pediatric Mode to ensure accurate measurements for younger patients
- Firmware upgradable

#9005BP, ADView 2 BP (126-9887)	ea
#9005BPTO, ADView 2 BP/Temp (126-8542)	ea
#9005BPS, ADView 2 BP/SpO ₂ (127-1394)	ea
#9005BPSTO, ADView 2 BP/SpO ₂ /Temp (126-7166)	ea
#9005TO, Temperature Module, Oral/Axillary (127-1403)	ea
#9005TR, Temperature Module, Rectal (127-1404)	ea
#9005SM, SpO ₂ Module, Masimo (127-1401)	ea
#9005W, Wall Mount with Basket (127-1396)	ea
#9005D, Desktop Caddy (127-1400)	ea
#9005CKA, Adult BP Cuff Kit (Adult, Small Adult, Large Adult) (127-1472)	ea
#9005CKP, Pediatric BP Cuff Kit (Child, Small Adult, Adult) (127-1471)	ea
#9005CK, Full BP Cuff Kit (Child, Small Adult, Adult, Large Adult, Thigh) (127-1473)	ea
#900513TBR1MB, Thigh BP Cuff (127-1475)	ea
#9005SP, SpO ₂ ChipOx Finger Sensor, Adult (127-1470)	ea
#9000TP, Thermometer Probe Covers (499-6964)	500/pkg

Advantage™ 6015N Digital Wrist Blood Pressure Monitor AMERICAN DIAGNOSTIC CORPORATION

Featuring Smart Logic™ technology and advanced memory capabilities, this monitor measures systolic and diastolic pressure as well as pulse. The Smart Logic™ technology determines the ideal inflation level.

Key features:

- Pressure accuracy: +/- 3 mm Hg
- Pulse range: 40–200 BPM
- Stores the last 99 readings for 2 users
- Wrist cuff fits sizes from 5.3" to 7.7"
- Automatic shut-off feature
- Made without natural rubber latex
- Includes 2 AAA batteries and a compact storage case
- 5-year warranty

(216-1353)



Advantage™ 2200 Fingertip Pulse Oximeter AMERICAN DIAGNOSTIC CORPORATION

This portable device is perfect for use at home, during sports activities, or in aviation. It offers precise readings with an SpO₂ accuracy of +/-2% and a heart rate accuracy of +/- 5 BPM.

Key features include:

- 2 adjustable display modes
- Bright red LED display
- Compatible with finger thickness of 0.31" to 0.79"
- Includes a carry lanyard and instruction booklet
- Powered by 2 AAA batteries (included)
- 1-year limited warranty

(124-6142)



Contact your Henry Schein
Rep for more information!

Above the Rest



X8

X10

X12

(138-9631)	X8 Patient Monitor w/ G2
(138-9695)	X8 Patient Monitor w/ Printer
(138-9709)	X8 Patient Monitor w/ Roll Stand
(138-9710)	X8 Patient Monitor w/ Roll Stand and Printer
(138-9826)	X10 Patient Monitor w/ G2
(138-9825)	X10 Patient Monitor w/ Printer
(138-9827)	X10 Patient Monitor w/ Roll Stand
(138-9828)	X10 Patient Monitor w/ Roll Stand and Printer
(138-8976)	X12 Patient Monitor w/ G2
(138-8900)	X12 Patient Monitor w/ Printer
(138-9992)	X12 Patient Monitor w/ Roll Stand
(138-9993)	X12 Patient Monitor w/ Roll Stand and Printer

Edan helps maximize patient care through reliable purchasing and EMR solutions. Edan G2 technology is compatible with generic consumables in order to maximize efficiency across your supply chain. All monitors are out of box ready with accessories.



M3 Vital Signs Monitor EDAN

The M3 Vital Signs Monitor is a light, compact, and portable device, ideal for spot-checks, triage, or the out-patient department. Its cost-effectiveness makes it available to all. Parameters include NIBP (manual or automatic), SpO₂, and quick oral temp.
#23.02.302050, Patient Monitor with NIBP and SpO₂
(974-0026) ea

#03.02.30212, Patient Monitor with NIBP, SpO₂, and Oral Temp
(974-0041) ea
#M3C.P.MT, Patient Monitor with NIBP, SpO₂, and Oral Temp, Printer, and Rolling Stand
(974-0187) ea
#VS.RT, Rolling Stand
(127-1081) ea

iM70 Patient Monitor EDAN

The iM70 Patient Monitor is embedded with the latest advancements in design to provide accurate monitoring across different applications that range from general care to anesthesia monitoring. Its high-resolution, 12" touch screen clearly displays all patient data (NIBP, ECG, SpO₂, respiration, temperature, EtCO₂ optional, anesthetic agents optional).

#iM70.P, Patient Monitor with Printer
(974-0178) ea
#iM70.CS.P, Patient Monitor with CO₂ and Printer
(974-0181) ea
#iM70_OR, Patient Monitor with Multi-Gas Analyzer and Printer
(974-0131) ea
#PM.RT, Rolling Stand
(127-1083) ea



iM60 Patient Monitor EDAN

The iM60 monitor spans accurate monitoring across different care environments, acuity levels, and clinical requirements. All data is clearly displayed on a 10½" color touch screen while providing real-time readings on NIBP, ECG, SpO₂, respiration, temperature, and optional EtCO₂.

#iM60.P, Patient Monitor with Printer
(974-0177) ea
#iM60.CS.P, Patient Monitor with CO₂ and Printer
(974-0180) ea
#iM60.P.MT, Patient Monitor with Printer and Rolling Stand
(974-0185) ea
#iM60.CS.P.MT, Patient Monitor with CO₂, Printer, and Rolling Stand
(974-0186) ea



iM80 Patient Monitor EDAN

The iM80 patient monitor is equipped with a massive, 15" color touch screen to clearly display patient readings and waveforms. Applications range from general cases to anesthesia monitoring. Parameters include: NIBP, ECG, SpO₂, respiration, temperature, EtCO₂ optional, anesthetic agents optional.

#iM80.P, Patient Monitor with Printer
(974-0179) ea
#iM80.CS.P, Patient Monitor with CO₂ and Printer
(974-0182) ea
#iM80_OR, Patient Monitor with Multi-Gas Analyzer and Printer
(974-0135) ea
#IM80.RT, Rolling Stand
(127-1139) ea



Mindray Ambulatory Care Solutions

www.mindray.com

Ambulatory Surgery Center and Clinic



The **A4 Advantage** complements clinical workflow with customizable ventilation capabilities and an intuitive touchscreen user-interface. It is well-suited for ambulatory surgery centers and offers outstanding functionality in a comprehensive, cost-effective anesthesia delivery system.



ePM 10M/12M/15M continuous, multi-parameter patient monitors are ideally suited for the ambulatory surgery center where flexibility, efficiency, and value are essential.

VS 8 is a reliable, easy-to-use vital signs monitor ideal for **spot check** vitals in an office setting. Offering **fast** new TrueTemp and True BP configurations.

All products offer an industry-leading 3-Year Warranty

To learn more about our entire portfolio of innovative solutions, contact your Henry Schein Representative.

mindray
healthcare within reach

ePM 10A Configured Patient Monitor **MINDRAY**

Exceptional user interface and enhanced parameter functionality, delivering a monitoring solution to fit multiple care environments. Immediate and measurable benefits of the smartphone-like user interface that delivers 98% of all common monitoring functions with just 2 gestures, reducing time clinicians spend navigating monitor functions, and more time on what matters most: patient care. Advanced monitoring options include 3- and 5-lead ECG, SpO₂, NIBP, respiration, and temperature supporting multiple clinical workflows. ePM 10A has a compact, ergonomic design that is perfectly suited for healthcare facilities managing a steady flow of patients daily. Add optional CO₂ capability to aid clinicians to identify signs of patient deterioration due to respiratory distress. Manufactured with industry-leading plastics that stand up to the most stringent medical-grade cleaners. The concealed grab-and-go handle eliminates hard-to-clean crevices and eases clinical workflow while on the move. Flexible mounting options for bedside and transport.

- Dimensions: 10.6" x 9.9" x 6.3"
- Weight: 8.8 lb
- Warranty: 3 years

Integrated
Sidestream CO₂, Recorder.....Ea

Specify:
Masimo SpO₂ (139-2003) Nellcor SpO₂ (139-2005)

Recorder.....Ea

Specify:
Masimo SpO₂ (139-1920) Nellcor SpO₂ (139-1927)



ePM 12MA Patient Monitor **MINDRAY**

Exceptional user interface and enhanced parameter functionality, delivering a monitoring solution to fit multiple care environments. Immediate and measurable benefits of the smartphone-like user interface that delivers 98% of all common monitoring functions with just 2 gestures, reducing time clinicians spend navigating monitor functions, and more time on what matters most: patient care. Advanced monitoring options are integrated: 3- and 5-lead ECG, SpO₂, NIBP, respiration, and temperature-supporting multiple clinical workflows. The ePM 12MA features a semimodular design with a plug-and-play module structure that adds monitoring functions such as CO₂ and multigas, where and when needed. The 12MA eGM multigas module automatically identifies and measures the inspired and expired values of O₂, CO₂, N₂O and 5 anesthetic agents, providing a compact, yet comprehensive OR monitoring solution. Manufactured with industry-leading plastics that stand up to the most stringent medical-grade cleaners. Concealed grab-and-go handle eliminates hard-to-clean crevices and eases clinical workflow while on the move. Flexible mounting options for bedside and transport.

- Dimensions: 12.2" x 11.4" x 6.7"
- Weight: 10.6 Lb.
- Warranty: 3 years

Semi-Configured/Modular.....Ea.

Specify:
Masimo SpO₂ (139-2126) Nellcor SpO₂ (139-2128)

ePM AG Module
(139-3276).....Ea.

CO₂ Module Adult/Pediatric
(137-8484).....Ea.





SIX HEALTH INDICATORS. ONE HANDHELD SOLUTION.

Rad-G® with Temperature

- > Measures multiple key health indicators in a single, rugged, portable device
- > Leverages the same technology trusted by hospitals (Masimo SET®) for accurate monitoring on the go, with an integrated non-contact temperature sensor



Oxygen Saturation



Pulse Rate




Perfusion Index



Pleth Variability Index



Respiration Rate from the Pleth



Temperature

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

PLCO-006965/PLMM-12507A-0823 PLLT-12038A
© 2024 Masimo Corporation. All rights reserved.



Masimo Set® LNCS® Sensors and Cables **MASIMO**

LNCS, or low-noise cabled sensors, are compatible with Masimo SET® pulse oximeters and Nellcor sensor-compatible instruments without OxiMax®.

#1863, LNCS DC-I Adult Reusable Sensor (114-1951)	ea
#2653, LNCS DBI, Adult, Soft, 3 ft (499-9823)	ea
#1864, LNCS DCIP Pediatric Reusable Sensor (116-0986)	ea
#1859, LNCS Amtx Adult Sensor (114-1911)	20/box
#2317, LNCS Amtx3, Adult, 3 ft (700-1279)	20/box
#1860, LNCS Pmtx Pediatric Sensor (499-6226)	20/box
#2318, LNCS Pmtx3, Pediatric, 3 ft (700-1280)	20/box
#2918, LNCS E1®, Adult, Ear, 3 ft (700-1785)	10/box

LNCS® Disposable, Adhesive SpO₂ Sensors with Replacement Tape

#2328, LNCS, Infant, 18" (116-0769)	20/box
#2319, LNCS, Infant3, 3 ft (119-3255)	20/box
#2329, LNCS Neo-L Neonatal Sensor (116-0770)	20/box
#2320, LNCS, Neonatal/Adult, 3 ft (119-1696)	20/box

LNCS® RED 20-Pin Patient Cables for Use with Masimo rainbow® SET® Devices

#2055, Red LNC-04 LNCS Patient Cable, 4 ft (114-5834)	ea
---	----

LNC Patient Cables for Use with Masimo SET® Devices

#2364, LNC1, 1 ft (125-7805)	ea
#2017, LNC-4 LNCS Patient Cable (700-1281)	ea



Rad-57® Handheld Pulse CO-Oximeter® **MASIMO**

#9216-U, Rad-57®, SpO ₂ Only (700-0668)	ea
#3734, Rad-57® Kit with SpCO® Option, Adult Sensor (700-1444)	ea
#3735, Rad-57® Kit with SpCO® Option (700-1303)	ea
#3736, Rad-57® with SpCO® and SpMet® Options (700-1308)	ea
#2406, rainbow RC4, 4 ft (118-6479)	ea
rainbow® Reusable Sensors for rainbow® Devices (SpO ₂ , SpCO®/SpMet®)	
#2696, rainbow DCI, Adult, 3 ft (700-1273)	ea
#2697, rainbow DCIP, Pediatric, 3 ft (499-9933)	ea

ZOLL Medical's interactive AEDs guide you through the process of performing high-quality CPR with voice, text, and image prompts. Rigorous, automatic self-testing ensures they are always ready. Each AED has unique features to provide the best solution for your needs.

ZOLL®



AED Plus®

- * Real CPR Help® technology
- * Industry leading five-year battery and pad life
- * One-piece pad design

(866-7559) Semi-Automatic Unit

Package includes: AED, electrodes, batteries, carry case.



Powerheart® G5 AED

- * Intellisense® CPR Feedback (ICPR)
- * Non-polarized pads (can be placed in either position)
- * Dual language functionality

(139-5101) Semi-Automatic Unit

Package includes: AED, electrodes, battery.



ZOLL R Series®

- * Professional defibrillator/monitor
- * OneStep Complete™ electrodes with integrated three-lead and Real CPR Help® technology

(127-9979) Professional Defibrillator

Package includes: Monitor/defibrillator, battery, 2 sets of OneStep electrodes.

ZOLL AED Plus®

ZOLL

Bringing Lifesaving Defibrillation to Everyone

The only AED with CPR feedback and an ECG waveform. By supporting CPR as well as defibrillation, this device can help every unresponsive victim. Features Real CPR Help®, a tool that is able to actually see what user is doing and provides audio and visual prompts to help them do it with confidence and clarity.

- Provides a one-piece electrode for accurate placement
- Uses consumer lithium batteries that are easy to replace, long-lasting, and inexpensive

Semi-Automatic Unit

(866-7559)ea

Contains: AED, electrodes, batteries, carry case, ECG, wall-mount bracket & operator's guide.

Fully Automatic Unit

(119-7910)ea

Contains: AED, electrodes, batteries, carry case, ECG, wall-mount bracket & operator's guide.

CPR-D padz®—Adult Defibrillator Pad

(253-2611)ea

Pedi padz® II—Pediatric Electrode Pad

(203-7054)ea

Type 123 Lithium Batteries

(610-4769) 10/pkg

Soft Carrying Case

(106-1937)ea

3D Wall Sign

(110-5384)ea

AED Plus Simulator

(499-7385)ea



Powerheart® G5 AED

ZOLL

Features Intellisense™ CPR, which provides corrective feedback to ensure delivery of high-quality compressions. One-button, dual-language functionality and comprehensive self-tests give rescuers the confidence to act.

- Intellisense™ CPR feedback provides real-time voice and text corrective prompts for CPR quality
- Rescue Ready® conducts regular, automatic self-tests on pads, battery, and internal circuitry
- RescueCoach™ guides users through each critical step of the rescue with user-paced voice/text prompts
- Dual language change from primary to secondary language in the push of a button

#G5S80C-S, Semi-Automatic

(139-5101)ea

#G5A80C-S, Fully Automatic

(139-5102)ea



AED Pro® Automated External Defibrillator

ZOLL

An external defibrillator that provides advanced capabilities for BLS and ALS users with innovative Real CPR Help® technology for rate and depth of chest compressions. See-Thru CPR® capability filters CPR chest compressions, displaying the underlying cardiac rhythm. Unit is unmatched in its ruggedness and durability and is ideal for any environment.

- Available as semiautomatic only or semiautomatic with manual override for ALS professionals
- Built-in internal memory and external USB memory technology allows expanded data storage capacity and maximum flexibility for data management and transfer
- Electrodes and battery compatibility with ZOLL M Series®, E Series®, and R Series® professional defibrillators
- ECG monitoring with 3-lead cable and high-quality display
- Unrivaled battery options (long-life disposable or rechargeable) that meet specific needs

(869-6270) ea



GE HealthCare

Patient Care Solutions Monitoring

Anesthesia

Carestation 620

Compact
Smart
Simple
Familiar



Carestation 650

Compact
Smart
Simple
Familiar



Kardia 12L



Portable
Efficient
Easy



Contact Us

For more product information
or to reach your local rep,
please click the QR
link to the right.



Portrait VSM

Care with
Confidence



B105 / B125 / B155

Simple, Durable, Versatile



Maternal Infant Care

Corometrics 170

NST Fetal Monitor



Bilisoft

Phototherapy
System



Cardiology

MAC 5 / MAC 5 Lite

Simple
Precise
Secure



MAC 7

Intuitive
Connected
Secure



MAC VU360

Smart
Integrated
Secure



CASE Stress

Clinical Excellence
Signal Quality
Seamless Workflow





Vscan Extend* R2 Dual Probe DICOM Handheld Ultrasound GE HEALTHCARE

A general-purpose diagnostic ultrasound imaging system that can enable qualified and trained health care professionals to visualize and measure anatomical structures and fluid. Its pocket-sized portability and simplified user interface enable integration into examination and training sessions. Information can be used for basic or focused assessments, procedural guidance, and in conjunction with other medical data for clinical diagnostic purposes during routine periodic monitoring and triage assessments for adult and pediatric patients. Unit supports applications for cardiac, abdominal, renal, obstetrics and gynecology, urology, fetal, adult cephalic, peripheral vascular imaging (arteries, veins, lower extremity, and carotid), small organs including thyroid and musculoskeletal (long bone, hip, shoulder, elbow, and knee joints), evaluation of presence of fluid, thoracic/lung (pleural motion/sliding and line artifacts), pediatrics, imaging guidance for needle/catheter placement (paracentesis, pericardiocentesis, thoracentesis, and amniocentesis), and procedure guidance for arterial or venous vessels (central lines and upper extremity). DICOM** configuration provides the ability to export images wirelessly to DICOM servers or to shared Windows*** network folders, along with wireless access to the GE Marketplace for download and installation of Vscan Extend apps. It also allows the transfer of images or videos via USB cable connection to a PC. The product is complemented with a Vscan* web portal, including online access to product and clinical information. In some countries, 1–60-gm bottle of gel may also be included in the standard configuration.

#H8038VP, DICOM Base Package (133-4462) ea

Contains: 1 control unit with standard components, 1 global AC adapter with interchangeable region-specific plugs, 1 USB cable, 2 rechargeable batteries, 1 external battery charger, 2 micro-SD memory cards (for potential service needs only), 1 soft case & 1-hour session of clinical applications training to be provided via telepresence.

#H8038VF, DICOM POCUS Focus Package (133-4461) ea

Contains: 1 control unit with standard components, 1 global AC adapter with interchangeable region-specific plugs, 1 USB cable, 2 rechargeable batteries, 1 external battery charger, 2 micro-SD memory cards (for potential service needs only), 1 soft case, 1-hour session of clinical applications training to be provided via telepresence & 1 subscription to the Point of Care Ultrasound FocusClass.

#H8038VS, DICOM with SonoSim 365 Package (133-4464) ea

Contains: 1 control unit with standard components, 1 global AC adapter with interchangeable region-specific plugs, 1 USB cable, 2 rechargeable batteries, 1 external battery charger, 2 micro-SD memory cards (for potential service needs only), 1 soft case, 1-hour session of remote clinical applications training & 1 SonoSim 365 for GE Healthcare training solution.

*Vscan Extend and Vscan are trademarks of General Electric Company.

**DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

***Windows is a registered trademark of Microsoft Corporation in the United States and/or other countries.



Vscan Air™ GE HEALTHCARE

This wireless, handheld ultrasound system enables whole-body scanning and delivers crystal-clear images.

- The power of a high-performance ultrasound machine in a small, lightweight wireless solution
- Optimized for the most common adult and pediatric physical exams to support a broad range of applications
- The dual probe's 2-sided design allows for both shallow and deep exams

#H8031VA, Vscan Air Base Package (140-0608) ea

#H45611AG, Vscan Air Protective Carrying Case (140-0615) ea



Accessories for Vscan™ Extend* R2 Dual Probe DICOM Handheld Ultrasound GE HEALTHCARE

#H41212RY, Global AC Adapter Secondary or replacement global AC adapter with 100–240 V, 50/60 Hz, and global plugs that can be connected with Vscan Extend, or with separate external battery charger. (135-2055) ea

#H41212SB, Robust Case, Complete Set Rugged design for transporting Vscan Extend with standard configuration and all accessories. (135-2058) ea Contains: 1 case, 1 global AC adapter with interchangeable region-specific plugs, 1 USB cable, 2 micro-SD memory cards, 1 soft case, 1 user manual CD & 1 hard-copy quick card.

#H41212SH, Robust Case, Scanner Set Rugged design for transporting Vscan Extend. (135-2060) ea Contains: 1 case, 1 optional extra battery, 1–60-g bottle of gel & 1 USB cable.

#H41212SA, External Battery Charger Adapter unit to charge Vscan Extend battery taken out from device. Global AC adapter not included. (135-2057) ea



Sheathes™ Vaginal/Rectal Ultrasound Probe Covers

MEDICAL RESOURCES

1¼"D x 8"L, Straight, Rolled, Latex	
#10001, Non-Sterile, Seamless	
(100-6867)	100/box
#10020, Non-Sterile, Seamless	
(103-1597)	600/bag
#10080, Non-Sterile, Individually Wrapped	
(103-9008)	100/box
#10040, Sterile, with 1 Elastic, Seamless	
(103-5591)	50/box
1"D x 9"L, Non-Sterile, Tapered, Not Made With Natural Rubber Latex	
#20001, Semi-Rolled on Paper Back	
(103-9004)	100/box
#20080, Semi-Rolled on Paper Back, Individually Wrapped	
(103-7983)	60/box
#25001, Rolled	
(421-4286)	100/box
#25080, Rolled, Individually Wrapped	
(358-2697)	100/box
1"D x 11.8"L, Tapered, Rolled, Extra-Long, Not Made With Natural Rubber Latex	
#25340, Sterile, With 1 Elastic	
(113-4410)	50/box
#27380, Pre-Gelled, Non-Sterile, Individually Wrapped	
(106-5607)	30/bag
¾"D Tapered to 2.9"W x 11.8"L	
#27690, Non-Sterile, Rolled with Pull Tab, Extra-Long, Individually Wrapped, Pre-Gelled, Not Made With Natural Rubber Latex	
(114-1035)	30/bag
1"D x 8"L, Non-Sterile, Tapered, Bidirectionally Rolled, Not Made With Natural Rubber Latex	
#42001	
(402-5928)	100/box
#40001, With 3" Debris Shield, Bulk	
(103-1311)	40/box



Blue Ultrasound Gel

Multipurpose, hypoallergenic, medium-viscosity ultrasound gel that is bacteriostatic for patient satisfaction. Excellent conductive agent with superior lubricating and clinical features for the most common professional diagnostic and imaging procedures. Simple operation and disposal. Long shelf life.

8½-oz Bottle, Blue	
(900-4094)	ea
5-liter Bottle, Blue	
(900-4095)	ea



Aquasonic® 100

PARKER LABORATORIES

60-g Tube, Blue	
(100-7634)	ea
0.25-liter Bottle, Blue	
(101-1042)	ea
1-liter Bottle with Dispenser, Blue	
(100-5242)	ea
5-liter SONICPAC® with Refillable Dispenser, Blue	
(100-1833)	ea



Clear Ultrasound Gel

Multipurpose, hypoallergenic, medium-viscosity gel is clear, dye-free, and bacteriostatic for patient satisfaction. It is also an excellent conductive agent with superior lubricating and clinical features, making it suitable for the most common professional diagnostic and imaging procedures. For use in therapeutic ultrasound, diagnostic imaging, and other external procedures that require a thick medium with excellent lubricating and glide characteristics. Will spread evenly without running or drying under standard clinical applications.

8½-oz Bottle, Clear	
(900-4352)	ea
5-liter Bottle, Clear	
(900-4353)	ea



Aquasonic Clear®

Acoustically correct for the broad range of frequencies used. Formula is hypoallergenic, bacteriostatic, and nonirritating.

- Water soluble
- Fragrance-free
- Colorless
- Formaldehyde-free
- Not a spermicide

8.5-oz Dispensers	
(727-9534)	12/box
1-liter Bottle with Dispenser Cap, Clear	
(101-2417)	ea
5-liter Bottle with Dispenser, Clear	
(101-9866)	ea
2-oz Tube, Clear	
(101-0693)	ea
5-liter Econopaks, Clear	
(426-2081)	4/case

NEW and Platinum Certified Preowned Systems


X-ray, CT, MRI, PET/CT, C-Arms, Mammography, and Ultrasound



A reliable resource for imaging professionals for nearly 70 years.



- Equipment
- Medical Supplies
- Mobile Imaging
- Service
- Parts



Contact Henry Schein to Learn More!

Acclimate® Gel Warmer

Designed to provide safe gel warming for EcoVue® and EcoVue® HV Ultrasound Gels to help enhance patient comfort. The warmer is available in two configurations — one for single-use packets and one for multi-use bottles.

- LCD display allows for temperature regulation (95°F-105°F)
- Two-stage thermal cut-off feature protects the warmer from overheating
- Suitable for wall mounting or countertop use

Item #	Description
(142-2675)	Acclimate® Gel Warmer, Bottle Warmer
(142-2681)	Acclimate® Gel Warmer, Packet Warmer
(136-5581)	EcoVue® Sterile Ultrasound Gel, 20 g Packet
(136-5586)	EcoVue® Non-Sterile Ultrasound Gel, 20 g Packet
(135-4814)	EcoVue® Non-Sterile Ultrasound Gel, 250 ml (8.5 oz.) Bottle
(139-8370)	EcoVue® HV Sterile Ultrasound Gel, 20 g Packet
(139-8366)	EcoVue® HV Non-Sterile Ultrasound Gel, 20 g Packet
(139-8375)	EcoVue® HV Non-Sterile Ultrasound Gel, 250 ml (8.5 oz.) Bottle

Bottle Warmer holds (2) 8.5 oz. bottles | Packet Warmer holds (14) 20 g single-use packets

For more information, contact your local Henry Schein Medical Representative or visit [henryschein.com/medical](https://www.henryschein.com/medical).

Acclimate, EcoVue, HR HealthCare, and the HR HealthCare logo are trademarks of HR HealthCare. © 2025 HR HealthCare. All rights reserved. LBL-0132 Rev.0 (09/2025)





BRASSELER
USA® SURGICAL

IMPROVE YOUR BOTTOM LINE

Save 40-60%*
over the OEM
when you switch
to Brasseler



The Right Product

We provide non-surgeon preference, high volume replacement products



The Right Price

Achieve significant cost reduction when you choose Brasseler orthopedic cutting accessories



High Quality

Made in the USA of high quality surgical grade materials, adhering to strict quality standards



OEM Compatible

All of our products are compatible with your current power systems



Easy Conversions

Dedicated Brasseler representative to support clinical cases and guide you through the conversion process



Follow Us!

Visit our website at BrasselerSurgical.com

To order call 800.535.6638

California customers please use website when ordering for Prop 65 information.

©2025 Brasseler USA. All rights reserved.

*Based on analysis of publicly available retail pricing. Data on file and available upon request.



50 Years of dependable high speed drills for Oral and Maxillofacial Surgeons



Surgairtome Two (148-3277)



Oral Max™ Elite High Speed Drill (140-5068)



Hall MicroFree™ High Speed Drill (131-8761)

HALL®
POWERED INSTRUMENTS

CONMED

To experience an easier day in the OR, contact your
Henry Schein Representative today for more information.



Hyfrecator® Pencil Sheaths CONMED

For use with all models of Hyfrecator
electrosurgical units.

#7-796-19, Sterile, 3'	25/box
(520-9814).....	
#7-796-18, Non-sterile, 8"	100/box
(686-4301).....	



Hyfrecator® 2000 CONMED

Electrosurgical Unit

Dual-microprocessor technology makes it one of the safest units for in-office procedures. Controlled Tissue Response™ minimizes tissue damage, resulting in excellent clinical outcomes. Over \$100 in accessories included at no extra charge.

- Procedural versatility: With a multitude of applications in such diverse specialties as dermatology, plastic surgery, gynecology, dental and oral surgery, ophthalmology, otolaryngology, proctology, and urology, this unit lets you perform a wide variety of desiccation, fulguration, and coagulation office procedures. Includes 50 blunt and 50 sharp, non-sterile disposable hyfrecator tips
- Precise power settings: Digital 2/10 incremental readouts in the critical 0–5-watt low-power range enhance your skill, accuracy, and control when performing delicate procedures
- Convenient fingertip control: ergonomic, handheld pencil enables you to remotely activate the unit and make rapid or slow power adjustments

Complete Unit

(520-9134).....ea
Contains: 1 remote control, 1–#7-800-5 reusable switching pencil, 2 Electrolase® disposable electrodes (1 of ea: #7-100-12 sharp & #7-101-12 blunt), 2–#7-796-19 sterile sheaths, 2–#7-796-18 non-sterile sheaths, 1–#7-796-20 wall-mounting kit & 1 operator's manual.

Aspen Surgical® Electrosurgery Solutions

Electrosurgery



SURGI-CENTER PRO™



OR PRO 300™



SMOKE SHARK® II
WITH BOVIE® ORCA®

Contact your Henry Schein Sales Representative at
www.henryschein.com/medical to purchase yours today!

122-9913 | OR PRO A3350 ELECTROSURGICAL GENERATOR | 1/EA

124-9073 | SURGI-CENTER PRO A2350 ELECTROSURGICAL GENERATOR | 1/EA

119-3868 | SMOKE SHARK® II SMOKE EVACUATOR | 1/EA

139-6511 | BOVIE® ORCA® INTEGRATED SMOKE EVACUATION PENCIL | 10/BOX

aspensurgical®

©2024 Aspen Surgical Products, Inc., ALL RIGHTS RESERVED M-SM-ADV-1167129, Rev I



Aaron® Surgi-Center Pro Electrocautery Generator

ASPEN SURGICAL

The Aaron® Surgi-Center Pro is a 200-watt, high-power operating-room generator. The unit incorporates the power and safety of a premium OR generator at a significant savings to the medical facility. The system delivers 12 monopolar and bipolar RF energy modalities to provide a high level of surgical precision.

- Offers pure cut, 4 modes of blend, 3 modes of coagulation, bipolar coagulation, and Bovie FCFS® technology
- Designed with state-of-the-art safety features
- Backed by a 4-year manufacturer's limited warranty

#A2350
(124-9073)..... ea



Aaron® Smoke Shark™ Smoke Evacuator

ASPEN SURGICAL

Designed to be easy to use and is extremely effective at smoke and particle removal at the operative site. While designed for use with all electrocautery generators, it is also excellent for use in controlling laser smoke.

- Easy to operate
- Quiet and effective smoke evacuation
- Long-life, 18-hour filter
- Low cost per procedure
- Excellent particle capture

#SE02
(119-3868)..... ea



High-Frequency Desiccator

HENRY SCHEIN®

Designed and engineered to be the most reliable and durable desiccator available today, this unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat. Post-procedure and long-term cleanup are simplified by its nonabsorbing, nonstaining, high-impact plastic case.

- 40-watt high-frequency desiccator
- Adjustable in 1/10-watt increments in the first 10 watts, in monopolar and bipolar
- More accurate power with Universal Power Control 100–240-volt AC
- Less to do and less to buy makes this the most cost-effective and user-friendly unit on the market
- This desiccator is a ground reference unit, which means no grounding pad is required
- 4-year warranty

#A942
(570-0583) ea



Electrosurgical Generator and High-Frequency Desiccator

HENRY SCHEIN®

This electrosurgical generator is a high-frequency desiccator with cutting ability. One half of the unit is the Henry Schein High Frequency Desiccator.

The other half is an electrosurgical generator that features 50 watts of cutting power and coag capabilities—plenty of power for most cutting procedures performed in the doctor's office.

- 40-watt high-frequency desiccator
- 2 machines in 1 convenient package
- Adjustable in 1/10-watt increments in the first 10 watts, in monopolar and bipolar
- No grounding pad is required during desiccation mode
- 4-year warranty

#A952
(570-0582) ea

Desiccator Tips

HENRY SCHEIN®

Utilize the standard 3/32" (2.4 mm) stainless steel shaft. All molded electrodes feature safety-grip insulators combining patient and user safety with easy insertion into and removal from the surgical pencil.



Sharp Dermal Tips, Non-Sterile (Bulk Packaged) (112-7171).....	100/box
Sharp Dermal Tips, Sterile (Packaged Individually) (112-7172).....	50/box
Blunt Dermal Tips, Non-Sterile (Bulk Packaged) (112-7173).....	100/box
Blunt Dermal Tips, Sterile (Packaged Individually) (112-7175).....	50/box
Disposable Handpiece Sheaths, Sterile (112-7176).....	25/box
Sharp H-Type Tip, Sterile (Packaged Individually) (570-0888).....	50/box
Sharp H-Type Tip, Non-Sterile (Bulk Packaged) (570-0890).....	100/box
Blunt H-Type Tip, Sterile (Packaged Individually) (570-0889).....	50/box
Sharp H-Type Tip, Non-Sterile (Bulk Packaged) (570-0891).....	100/box

Midmark M11® Steam Sterilizer

MIDMARK

The M11 is engineered with user-friendly features to simplify compliance, streamline cycle operations, and enhance device care tracking. These improvements help minimize transcription and documentation errors, which can account for up to 56% of sterilization issues. Automated mechanical monitoring and chemical indicator reporting not only save time but also facilitate easier staff training, instilling confidence in the sterilization process.

Key features:

- 5" fingerprint-resistant touchscreen allows for easy navigation of cycle setup and other operations, even while wearing gloves. Users can choose between light and dark modes to suit their preferences.
- Stay prepared for audits with built-in reminders, notifications, and user authentication. The sterilizer offers unlimited storage for routine care events throughout its lifespan, along with automated cycle recordkeeping.
- Midmark sterilizers are built to last, with an upgraded chamber design that extends the device's life to 25,000 cycles. This can reduce the need for maintenance by both users and certified technicians by more than half.
- Height: 18.6" (47.2 cm), Width: 17.8" (45.2 cm)
- Chamber dimensions: 11" x 18" (27.9 cm x 45.7 cm)
- Water reservoir capacity: Full: 1.4 gal (5.3 liter); Usable: 0.8 gal (3 liter)

#M11-050, 115 Volt
(148-3221) ea

#M11-051, 230 Volt
(148-8941) ea



Vistacool™ Direct to-Drain Cooling System for Autoclave Wastewater

MIDMARK

Efficiently cools sterilizer wastewater before safely directing it to the drain. This non-electric, self-monitoring, and self-regulating system allows users to focus more on patient care and less on equipment maintenance.

Key benefits:

- Eliminates the need to handle hot wastewater tanks, preventing potential injuries and avoiding damage to cabinetry from spills or overflows
- Prevents scalding hot wastewater from damaging plumbing with built-in backflow prevention that meets or exceeds local plumbing codes
- Automates a time-consuming task with a simple, self-monitored system, reducing the need for manual intervention
- Dimensions (including fittings): 4"W x 13"L x 17"H (10.2 cm x 33.1 cm x 43.2 cm)
- Tank volume: 3 gal (11.4 liter)

#9A586002
(142-4954) ea



Cool Hand Tool

MIDMARK

Attaches to trays and cassettes, helping to minimize the risk of injury for the user when loading or unloading the sterilizer.

#9A307001
(934-6594) ea



Automated High-Level Disinfection (HLD) of ultrasound probes just got better.

The trophon®3 device is the ultimate solution for healthcare environments where patient safety and efficiency are top priorities and compliance is non-negotiable.

- ✓ Enhanced experience
- ✓ Streamlined workflow
- ✓ Built-in traceability

With the trophon ecosystem, every component works seamlessly to deliver an end-to-end reprocessing solution.



trophon³
**40%
FASTER**
than previous
trophon technology

Contact your Henry Schein representative for more information.



The trophon® family includes trophon® EPR, trophon®2, and trophon®3 devices which share the same core technology of 'sonically activated' hydrogen peroxide. © Nanosonics 2025. All rights reserved. 2226517-US-AD V01. September 2025.



trophon®2 Ultrasound Probe Disinfection Device NANOSONICS

This fully enclosed, automated high-level disinfection (HLD) system streamlines workflow with a simple, fast, and integrated solution. Designed specifically for patient examination rooms, the trophon2 enables efficient point-of-care reprocessing.

Key features include:

- Enhanced user compliance and audit readiness through automation, with comprehensive data capture throughout the reprocessing workflow
- AcuTrace® RFID technology that records operator, probe, and cycle data to ensure and demonstrate compliance
- A printer that generates labels with AcuTrace inputs and disinfection results
- The ability to download disinfection records to a USB
- An enclosed system with an easy-close door that seals automatically
- A safe-to-handle disinfection cartridge operable with gloves
- No special ventilation requirements
- Environmentally friendly design for safe and easy disposal
- Compatible with more than 1,300 ultrasound probes; approved and endorsed by 26 probe manufacturers

#N05000-US
(142-4635) ea



trophon Sonex-HL® Disinfectant NANOSONICS

This proprietary disinfectant liquid contains 35% hydrogen peroxide, formulated for effective high-level disinfection (HLD) of the entire ultrasound probe, including the shaft and handle. The disinfectant packaging features AcuTrace® technology, allowing batch and expiration details to be stored directly in the trophon2 device.

#N05002
(142-4628) ea



trophon Chemical Indicator NANOSONICS

This indicator provides an independent confirmation of the success of each trophon disinfection cycle. The chemical indicator packaging includes AcuTrace® technology, enabling the recording and storage of batch and expiration details directly in the trophon2 device.

#N05003
(142-4646) ea



trophon Printer NANOSONICS

The trophon printer is a convenient traceability solution that prints high-quality disinfection labels for each cycle. The labels are suitable for medical records and have a minimum lifespan of 7 years.

#N00048-US, Printer
(142-4643) ea
#N00049-ROW, Printer Label Roll
(142-4648) ea

COLD STORAGE



Pharma-Vac Performance Series Vaccine Refrigerators **ACCUCOLD**

This line of refrigerators supports meeting regulations and guidelines from the CDC, State Agencies, World Health Organization, and other health ministries on vaccine storage and handling. Designed with a hydrocarbon refrigerant, they provide an energy-efficient solution to vaccine storage. Units are Pyxis, Omnicell, and AcuDose RX compatible.

- Certified to the NSF/ANSI 456 Vaccine Storage Standard
- Adjustable-cycle defrost operation ensures the unit avoids temperature spikes and ice buildup
- Battery-powered alarm sounds if the power connection is interrupted
- Digital controls with a large LED readout in Celsius or Fahrenheit for intuitive temperature management
- Multiple shelf positions spaced at ½" to accommodate virtually any size
- Ecofriendly hydrocarbon-based refrigerant for a greener footprint that meets current CARB regulations
- Environmentally friendly design without ozone-damaging chemicals

#ARG1PV456, 1-cu-ft Compact	
(143-5669)	ea
#ARG15PV456, 15-cu-ft Upright	
(143-5665)	ea
#ARG18PV456, 18-cu-ft Upright	
(143-5753)	ea



General-Purpose Freestanding Undercounter Refrigerator **AMERICAN BIOTECH SUPPLY**

Intended for general-purpose applications only where precise temperature control and maintenance are not essential and there are no requirements for temperature alarms.

- Swing door with keyed door lock
- 3 shelves (1 fixed/2 adjustable)
- Interior light
- HFC-free refrigerant (R600A), foam insulation, and packaging
- UL/C-UL listed
- Temperature range: 0°C to 10°C
- Capacity: 5 cu ft
- Exterior dimensions: 20¼"W x 22" D x 31½"H
- 1-year parts and labor warranty

#ABT-HC-UCFS-0504W	
(125-4697)	ea

This unit must have 4" clearance on sides and back for adequate ventilation.



Vaccine Temperature Data Logger



- High accuracy
- 24 million readings
- Maximum of 170 files (.csv)
- Internal (ambient) and external (fridge/freezer) sensors reading display
- Minimum and maximum reading memories
- Low and high alarm limits setting
- Flashing light alarm indicator
- °C or °F scale
- Real-time display
- Waterproof external sensor
- Low battery indicator
- Removable SD memory card or USB interface for data retrieval
- Optional external power

(570-2722)	ea
------------	----

COLD STORAGE



Premier Pharmacy Dual-Temp Refrigerator & Auto-Defrost Freezer Combination Unit

AMERICAN BIOTECH SUPPLY

- Separate microprocessor temperature controller allows for precise temperature control, verification, and recovery
- Interior fan-driven, forced-draft air-circulation system provides excellent temperature uniformity and recovery after door openings
- Equipped with digital temperature displays for convenient monitoring
- Critical samples and supplies are security backed with an extensive package of audible and visual alarms
- Freezer unit features an automatic defrost operating system
- 3 adjustable shelves for refrigerator, 2 adjustable shelves for freezer
- Shelves have back rails to prevent product from being pushed against back wall
- Capacity: 7 cu ft total (5.2 cu ft for refrigerator, 1.3 cu ft for freezer)
- 53 $\frac{1}{8}$ "H x 23 $\frac{3}{4}$ "W x 24 $\frac{1}{4}$ "D
- 1– $\frac{3}{8}$ " probe access port for refrigerator, 1– $\frac{3}{8}$ " probe access port for freezer
- Keyed door locks

#PH-ABT-HC-RFC7A

(125-4637)ea



Pharmacy Undercounter Built-In Refrigerator

AMERICAN BIOTECH SUPPLY

- 4.6-cu-ft capacity
- Features a microprocessor temperature controller and display that allows for precise temperature control, verification, and recovery (vs. mechanical-style thermostats)
- Acceptable temperature ranges for vaccine storage are 36°F to 46°F (2°C to 8°C) in the refrigerator and -58°F to 5°F (-50°C to -15°C) in the freezer
- Keyed door lock
- 3 shelves with guardrail on back to prevent product from being pushed against the back wall
- 115V, 60 Hz, $\frac{1}{8}$ hp, 1.74 amps
- 23 $\frac{3}{4}$ "W x 24 $\frac{1}{2}$ "D x 33 $\frac{3}{8}$ "H

#PH-ABT-HC-UCBI-0404

(125-4626)ea



Compact Vaccine Refrigerator

ACCUCOLD

- Designed and purpose-built for pharmacy, medication, and vaccine applications to support meeting CDC/VFC/VFA guidelines
- Meets UL-471 commercial standards
- Antimicrobial, silver ion-coated handle is powder-coated with naturally occurring silver ions to help reduce germ spread
- Electronic controller is equipped with a volt-free contact for alarm forwarding to an external remote system (contacts located in the external rear of the unit)
- Control panel requires password entry to adjust the parameters beyond the set point
- Digital controls with a large LED readout in Fahrenheit or Celsius for intuitive temperature management: from 35.5°F to 46.4°F (2°C to 8°C)
- Capacity: 1 cu ft (28 L)
- 15 $\frac{3}{4}$ "H x 14"W x 8"D

#ARG1PV

(136-3278)ea



Premier Pharmacy/ Vaccine Undercounter Refrigerator

AMERICAN BIOTECH SUPPLY

- Features a microprocessor temperature controller and display that allows for precise temperature control, verification, and recovery (vs. mechanical-style thermostats)
- Freestanding model ("UCFS") requires 4" of clearance around back and sides of the unit for proper ventilation. Improper installation will void the warranty and will lead to temperature maintenance issues, overheating, and possible failure
- Acceptable temperature ranges for vaccine storage are 36°F to 46°F (or 2°C to 8°C) in the refrigerator and -58°F to 5°F (or -50°C to -15°C) in the freezer
- 3 adjustable shelves with guardrail
- Bottom-footed shelf allows for air circulation
- Keyed door lock
- 115V, 60 Hz, 1/6 hp, 1.3 amps
- 23 $\frac{3}{4}$ "W x 24"D x 32 $\frac{3}{8}$ "H

#PH-ABT-HC-UCFS-0504

(138-5734)ea

Standardized, Easy Loading, Liquid, Ready-To- Use Reagents


Beckman Coulter Chemistry
and Immunoassay Systems



(#104-7912)

CAPITAL LABORATORY

HSI#	MFR#	Description	Size
(117-4246)	628017	Diluent f/DXH	Ea
(117-8884)	A16793	UniCel DXI Wash Buffer II	1/Bx
(133-9880)	73627	Aution 9EB Strips 100/Bottle	100/Btl
(126-7011)	B63284	Access TSH Reagent 3rd IS	200/Bx
(331-0052)	33000	Access B12	100/Kt
(331-0116)	33880	Access Free T4	100/Kt
(331-0134)	37200	Access Hybritech PSA	100/Kt
(108-5050)	A16792	Wash Buffer II	4Bt/Bx
(135-6108)	B36868	FP DXH 500 Series Cleaner	500mL/Bt
(861-0698)	A98856	Access 25(OH) Vit D Total	2x50/Bx
(133-8778)	B1017-435	Microscan Negative/Urine Combo	Ea
(331-0159)	81906	Substrate	4Bt/Bx
(117-3886)	A98032	Access Folate2x50 Determi	2x50/Ea
(331-0059)	33020	Access Ferritin	100/Kt
(133-9911)	660-3021	Urinalysis Test Tube 16x100mm	250/Bx
(122-9827)	B1026-10D	Microscan System Inoculation	60/Bx
(331-0089)	33560	Access Testosterone	100/Kt
(135-6101)	B36845	DxH 500 Diluent	Ea
(135-6106)	B36846	FP DXH 500 Series LYSE	Ea
(331-0111)	A13422	Access Free T3	100/Kt
(331-0109)	33830	Access Total T3	100/Kt
(135-6111)	B36872	FP DXH 500 Series Ctrl Kit Assembly	1/Bx
(142-8192)	OSR6178	Olympus Reagent Creatinine OSR	4x990/Bx
(747-2179)	OSR6295	Direct HDL Olympus	4/St
(331-0156)	81901	Reaction Vessels	16X98/Bx
(131-1913)	B84493	Access Sensitive Estradiol	1/Bx
(133-9472)	B1017-219	Panel Postv MicroScan Combo 45	20/Pk
(114-2975)	A48617	Access SHBG Reagent	2x50/Bx
(142-8186)	OSR6137	Bicarbonate Co2 Osr Reage	4X580/Bx
(331-0103)	33800	Access Total T4	100/Kt
(129-2467)	800-3236	Lamina iQ Series 200	2Bt/Ca
(861-0679)	B24838	Access 25(OH) Vit D Total	2x50/Bx
(336-9242)	OSR6189	Magnesium AU400 Test Kit	4x250/Bx
(331-0086)	33550	Access Progesterone	100/Kt
(117-4248)	628020	Diff Pack f/DXH	Ea
(331-0079)	33520	Access hFSH	100/Kt
(109-2666)	A16972	Access 2 Intact PTH 2x50Test	1/Bx
(109-7644)	A10826	Access DHEA-S	2x50/Bx
(331-0077)	33510	Access hLH	100/Kt
(122-1937)	A85264	Access Total BhCG Reagent	Ea
(331-0073)	33410	Access Ultrasensitive Insulin	100/Kt
(141-1364)	C48639	Panel Negative Urine Combo 89	20/Bx
(331-0095)	33600	Access Cortisol	100/Kt
(124-0099)	OSR6004	Alp Osr Reagent Olympus	4X310/Bx
(117-4247)	628019	Cell Lyse f/DXH	Ea
(129-8778)	A16793	UniCel DXI Wash Buffer Solutn	Ea
(331-0105)	33810	Access T-Uptake	100/Kt
(126-2413)	OSR6107	Alt Osr Reagent Olympus	4X980/Bx
(666-5676)	OSR6195	Direct HDL Reagent for Olympus	4/St
(116-3597)	AUH1012	Mid Standard Ise Reagent	4/Bx
(108-6459)	OSR60118	Olympus Triglyceride	4/Kt
(126-5601)	OSR6109	Ast Osr Reagent Olympus	4X980/Bx
(127-1976)	OSR6112	Bilirubin Total Osr Reage	4X590/Bx
(331-0038)	386167	Unicel Dxl 800 Reaction Ves	10bg/Bx
(109-2668)	A12985	Access TPO AB	2x50/Bx
(142-8179)	OSR6134	Blood Urea Nitrogen Osr Reagnt	4X620/Bx
(331-0081)	33530	Access Prolactin	100/Kt
(142-4989)	OSR6132	Total Protein Osr Reagent	4X750/Bx
(331-0138)	37210	Access Hybritech Free PSA	100/Kt
(109-3392)	A16364	Access EPO	100/Kt



VITROS® Syphilis Assay

**Accurate and early
infection detection**

The American College of Obstetricians and Gynecologists (ACOG) highlights that syphilis screening is crucial for pregnant women to prevent the transmission of congenital syphilis, which can pose significant prenatal risks. The VITROS Syphilis Assay helps clinicians accurately detect infections early to treat and prevent the spread of the disease.

Now available on the VITROS® 3600 Immunodiagnostic System and VITROS® 5600 and VITROS® XT 7600 Integrated Systems in the U.S.

**For more information,
visit henryschein.com/medical
or contact your local
Henry Schein Representative.**

New QuidelOrtho branding may not be available in all markets, subject to country-specific regulatory approval. Please confirm with your local commercial team.

© 2025 QuidelOrtho Corporation. All rights reserved. | PR-107569-NA-EN-US-v1

CAPITAL LABORATORY

HSI#	MFR#	Description	Size
(600-6259)	8379034	Vitros 250 Sodium Slides	5x50/Pk
(600-6266)	8157596	Vitros Potassium 250 Slides	5x50/Pk
(127-5604)	6844471	Vitros Chloride w/ Urine	5x50/Pk
(360-5433)	8392292	Vitros 250 Total Protein	250/Bx
(334-4037)	8196057	Vitros 250 Albumin	250/Bx
(340-7792)	8262396	Vitros 250 Slide Eco2	300/Bx
(101-0675)	1707801	Vitros Glucose 300 Slides	5x60/Pk
(370-9860)	1450261	Vitros Slides Calcium-60	300/Pk
(611-4331)	8102204	Vitros 250 BUN Slides-60	5x60/Bx
(365-6437)	8433815	Vitros 250 Slide AST	300/BX
(129-8122)	6844288	Vitros ALTV Reagent Slides	300/Bx
(858-1524)	1053180	Vitros Slides Alkaline Phos-60	300/Bx
(333-1286)	8159931	Vitros 250 Total Bilirubin	300/Bx
(858-4795)	6801895	Vitros dHDL Slides	5x60/Bx
(136-4879)	6844296	Vitros XT TBIL-ALKP 600 Slides	300/Bx
(982-4750)	1336544	Vitros Triglyceride 300 Slides	300/Bx
(101-0617)	1669829	Vitros Total Cholest Slides	300/Bx
(600-9890)	1142561	Vitros B12/Folate Reagent Pk3	3/Pk
(136-4633)	6844294	Vitros XT UREA-CREA 600 Slides	5/Bx
(136-4878)	6844295	Vitros XT GLU-CA 600 Slides	300/Bx
(370-9133)	1921204	Vitros Magnesium 90 Slides	90/Pk
(370-7492)	8383051	Vitros Slides BuBc-60	5x60/Bx
(330-2575)	8097990	Vitros 250 CRP Slides	90/Bx
(127-8167)	6844463	Vitros Reference Fluid 800	Ea
(370-0823)	8297749	Vitros Lipase 90 Slides	90/Bx
(131-0687)	6842781	Vitros HIV Combo Test	100/Bx
(779-7315)	8112724	Vitros Amylase 90 Slides	90/Bx
(602-1722)	1924547	Vitros Iron 90 Slides	90/Bx
(816-3901)	1612365	Vitros Bunc 90 Slides	90/Bx
(817-9840)	8478034	Vitros Ck 90 Slides	Bx
(136-4403)	6844292	Vitros XT TRIG-CHOL 600 Slides	5/Bx
(886-8412)	6802001	Vitros TIBC T.Iron Binding	1/Bx
(123-9761)	6842905	Vitros HbA1c Reagent	300/Bx
(370-0053)	1943927	Vitros Uric Acid 300 Slides	300/Pk
(370-4916)	8255093	Vitros 250 Magnesium Slide	5x60/Bx
(136-4211)	6844508	Vitros Cal Kit 32	Ea
(370-2858)	1513209	Vitros Phosphorus 300 Slides	300/Pk
(370-3612)	8329930	Vitros Triglycerides 90 Slides	90/Pk
(101-0879)	1882208	Vitros Cal Kit 1	4/Bx
(816-0714)	1515808	Vitros Slides Iron/TIBC-60	300/Bx
(264-9089)	6801740	Vitros mALB Micro Albumin	1/Bx
(318-5604)	6844471	Vitros 250 Chloride w/Urine SI	5x50/Pk
(370-3599)	1290709	Vitros Cal Kit 3	4/Bx
(370-6273)	1204668	Vitros Cal Kit 4	4/Bx
(318-6266)	8157596	Vitros 250 Potassium Slides	5x50/Bx
(739-2730)	6802584	Vitros 250 Creatinine IDMS	300/Bx
(739-6259)	8379034	Vitros 250 Sodium Slides	5x50/Bx
(318-5433)	8392292	Vitros 250 Total Protein Slide	250/Bx
(739-6266)	8157596	Vitros 250 Potassium Slides	5x60/Bx
(129-3222)	6844425	Vitros Reagent Pack	Ea
(739-5604)	6844471	Vitros 250 Chloride w/Urine SI	5x50/Pk
(318-4037)	8196057	Vitros 250 Albumin Slides	250/Bx
(318-2730)	6802584	Vitros 250 Creatinine IDMS	300/Bx
(146-1152)	6900052	Vitros Intact PTH II Reagent	Ea
(739-5433)	8392292	Vitros 250 Total Protein Slide	250/Bx
(600-9665)	6801896	Vitros Cal Kit 25	2/Bx
(318-7792)	8262396	Vitros 250 Eco2 Slides	300/Bx
(739-4037)	8196057	Vitros 250 Albumin Slides	250/Bx
(739-9860)	1450261	Vitros 250 Calcium-60 Slides	300/Bx
(318-4331)	8102204	Vitros 250 BUN-60 Slides	5x60/Bx



HARMONIZED SYSTEMS

CLINICAL CHEMISTRY, IMMUNOASSAY AND INTEGRATED SYSTEMS TO TRANSFORM YOUR LABORATORY



Clinical Chemistry | Immunoassay | Hematology | Transfusion | Molecular | Point of Care | Professional Services

CAPITAL LABORATORY

HSI#	MFR#	REAGENTS	
		Description	Size
(140-6514)	07P4820	Alinity I TSH	Ea
(143-4867)	04Z41-20	Alinity 6-Acetylmorphine	Ea
(141-8811)	04Z4520	Alinity C Oxycodone	Ea
(143-4726)	04Z42-20	Alinity Buprenorphine	Ea
(140-6441)	07P7020	Alinity I Free T4	200/Bx
(140-6511)	07P6920	Alinity i Free T3	Ea
(140-6553)	07P5020	Alinity I ESTRADIOL (200T)	Ea
(140-6454)	07P6821	Alinity I Testosterone II	Ea
(140-6989)	09P3720	Alinity Dehydroepiandrosterone Sulfate Reagent	200/Bx
(140-3411)	09P3820	Alinity I SHBG	Ea
(140-6501)	09P3521	Alinity I Anti-TPO	Ea
(140-6517)	08P4522	Alinity I Vitamin D	Ea
(140-6522)	08P0521	Alinity I ANTI-HCV	Ea
(141-8856)	09P5120	Alinity C Barbituates	Ea
(141-8859)	09P5420	Alinity C Cannabinoids	Ea
(141-8966)	09P6620	Alinity C Phencyclidine	Ea
(143-4722)	09P45-20	Alinity C Amphetamine/Methamphetamine	Ea
(143-4723)	09P58-20	Alinity C Ecstasy	Ea
(140-6478)	07P4920	Alinity Follicle Stimulating Hormone	200/Bx
(141-8880)	09P6520	Alinity C Opiates	Ea
(140-6494)	08P1021	Alinity I HBSAG 2	Ea
(140-6549)	07P9221	Alinity I TOTAL PSA US 200	Ea
(141-8882)	09P5520	Alinity C Cocaine	Ea
(140-6452)	08P3620	Alinity I Progesterone Reagent	200/Bx
(141-8858)	09P5220	Alinity C Benzodiazepine	Ea
(140-6497)	07P6520	Alinity I Ferritin Reagent	Ea
(140-6502)	07P6721	Alinity I B12	Ea
(137-8874)	04S1740	Alinity C Cartridge	Ea
(140-6583)	08P1421	Alinity I FOLATE (200T) US	Ea
(141-0278)	04T7520	Alinity I Insulin	Ea
(137-8598)	07P7220	Alinity C CO2	3000/Bx
(137-8855)	07P7520	Alinity C Ultra HDL	1400/Bx
(140-6533)	08P0721	Alinity I HIV COMBO	Ea
(140-6548)	07P9120	Alinity I LH	Ea

HSI#	MFR#	SOLUTIONS	
		Description	Size
(138-1489)	06P1368	Alinity I Wash Buffer Concentrated	Ea
(138-1485)	06P1401	Alinity I Vessel Reaction	Ea
(137-8745)	08P7840	Alinity C Alkaline Bath	Ea
(137-8744)	08P7740	Alinity C Acid Wash Solution	Ea
(137-8836)	08P9670	Alinity C Detergent A	Ea
(147-3667)	06P1170	Alinity Trigger Solution	4/Bx
(137-8863)	08P7640	Alinity C ICT Reference Solution	Ea
(137-8770)	08P4420	Alinity C Liquid UIBC	300/Bx
(147-3670)	06P1270	Alinity Pre-Trigger Solution	4/Bx
(137-8596)	07P5520	Alinity C Glucose	4000/Bx
(137-8609)	08P6001	Alinity C MMC Calibrator	Ea
(137-8811)	08P1925	Alinity C Magnesium	720/Bx
(138-1481)	01R3801	Alinity I Cups Sample	1000/Bx
(137-8601)	08P2020	Alinity C ALK Phos	4000/Bx
(137-8793)	01R6070	Alinity C Probe Wash	Ea
(137-8617)	07P5720	Alinity C Calcium	4000/Bx
(137-8616)	08P1620	Alinity C Urea Nitrogen	4000/Bx

Diagnostics that deliver, wherever care happens



CA-660
Coagulation



Atellica CI Analyzer
Integrated Chemistry and
Immunoassay



Immulite 2000 XPI
Immunoassay and Allergy



From large hospitals to non-acute care settings, Siemens Healthineers delivers a comprehensive portfolio of laboratory diagnostics solutions — from cutting-edge analyzers to a broad menu of high-quality assays. Wherever care happens, our diagnostics deliver trusted results that empower clinical decisions and improve patient outcomes.

SIEMENS
Healthineers

HSI#	MFR#	Description	Size
(147-7503)	11208706	Atellica IM TSH3ULII Reagent 130 Test	Ea
(133-0403)	10995569	Atellica IM Ferritin Reagent 90 Test	Ea
(133-0407)	10995588	Atellica IM fT4: Free Thyroxine Reagent 250 Tests	Ea
(110-8565)	10385206	IMMULITE 2000 Reaction Tube	1000/Bx
(133-0406)	10995585	Atellica IM fT3: Free Triiodothyronine Reagent 60 Test	Ea
(133-0455)	10995704	Atellica IM TSH3-UL Reagent 550 Test	Ea
(136-2543)	10995572	Atellica IM Folate Reagent	140/Bx
(133-0386)	10995538	Atellica IM Cor Reagent 50 Tests	Ea
(111-6646)	10385232	IMMULITE 2000/2500 Chemiluminescent Substrate 2x205mL For Bottle Liquid Volume	2Bt/Bx
(136-2560)	10995720	Atellica IM VitD(US) Reagent	500/Bx
(133-7638)	10995459	Atellica HIV Ag/Ab Reagent 100 Tests	1/Bx
(142-5404)	11537213	Atellica CH HDLC Assay	4/Bx
(142-5438)	11537222	Atellica CH Trig_2 Assay	4/Bx
(136-0902)	10995679	Atellica T3: Triiodothyronine Reagent 120 Test	Ea
(141-9407)	11417929	Atellica IM Solution Acid 2x1.5L	Ea
(141-9410)	11417930	Atellica IM Solution Base 2x1.5L	Ea
(132-9938)	10995576	Atellica IM Folate DTT/REL Kit	600/Bx
(133-7291)	10995663	Atellica IM PSA: Prostate Specific Antigen Reagent 500 Tests	1/Bx
(113-8515)	10381455	IMMULITE 2000 Insulin Reagent 200 Count Contains	Ea
(136-0907)	10995708	Atellica Testosterone II Reagent 100 Tests	Ea
(133-0459)	10995715	Atellica IM Vitamin B12 Reagent 500 Tests	Ea
(136-2550)	10995604	Atellica IM HBsII Reagent	200/Bx
(133-0398)	10995561	Atellica IM eE2: Enhanced Estradiol Reagent 100 Tests	Ea
(136-2591)	11203473	Atellica IM aHCV: Anti-Hepatitis C Antibody Reagent	200/Bx
(133-7032)	10995662	Atellica IM PSA Reagent	100/Bx
(361-4974)	10385393	IMMULITE 2000 Probe Wash	2/Bx
(110-7433)	10901844	IMMULITE 2000 Cortisol Reagent	Ea
(117-6599)	11129118	IMMULITE 2000 TSH: Thyroid Stimulating Hormone Reagent	Ea
(133-0435)	10995660	Atellica IM PRGE Reagent 90 Test	Ea

HSI#	MFR#	Description	Size
(133-0353)	10995471	Atellica IM BNP: Brain Natriuretic Peptide Reagent	Ea
(133-0519)	11097639	Atellica CH ASTPLc Reagent	Ea
(133-0396)	10995555	Atellica IM DHEA-S: Dehydroepiandrosterone Sulfate Reagent	Ea
(136-2555)	10995690	Atellica IM thCG Reagent	90/Bx
(133-0402)	10995568	Atellica IM Ferritin Reagent	Ea
(833-6025)	10385584	IMMULITE Canine - Control 2x2mL	EA
(133-0518)	11097638	Atellica CH ALTPLc Reagent	Ea
(136-2559)	10995719	Atellica IM VitD(US) Reagent	100/Bx
(113-9755)	10385229	IMMULITE 2000/2500 Probe Cleaning Kit	Ea
(133-0417)	10995615	Atellica IM tIgE: Total Immunoglobulin E Reagent	Ea
(133-0505)	11097618	Atellica CH RF: Rheumatoid Factor Reagent	Ea
(110-7451)	10380914	IMMULITE 2000 Vitamin B12 Reagent	Ea
(133-0405)	10995584	Atellica IM fT3: Free Triiodothyronine Reagent	Ea
(136-2538)	10995453	Atellica IM aHBs2 Reagent	200/Bx
(136-3967)	11097536	Atellica CH Hemoglobin A1c Reagent	Ea
(146-9481)	11537245	Atellica CH Calibrator Consumable	Bx
(146-9237)	11537243	Atellica CH Bilirubin Calibrator	6/Bx
(133-0388)	10995541	Atellica IM Reagent 100 Tests	Ea
(133-0420)	10995628	Atellica IM IRI Reagent 100 Tests	Ea
(112-4548)	10380911	IMMULITE 2000 Folic Acid Reagent	Ea
(113-3741)	10381622	IMMULITE fT4: Free Thyroxine Reagent	1/Bx
(133-0441)	10995672	Atellica IM SHBG: Sex Hormone Binding Globulin Reagent	Ea
(113-8547)	10381647	IMMULITE 2000 Thyroglobulin Reagent	Ea
(140-2819)	11202257	Atellica IM Alpha-Fetoprotein Reagent	Ea
(122-3526)	10901873	IMMULITE 2000 PROG: Progesterone Reagent	Ea
(146-9482)	11537223	Atellica CH Revised C-Reactive Protein	Bx
(136-2536)	10995449	Atellica IM aHBcM Reagent	100/Bx
(133-0443)	10995675	Atellica IM Syphilis Reagent	Ea
(133-0358)	10995478	Atellica IM BR Reagent	Ea
(133-0439)	10995670	Atellica IM Rub G Reagent	Ea
(133-0493)	11097530	Atellica CH AlbP Reagent	Ea

HORIBA

Make HORIBA Your **Trusted Partner** for Chemistry & Hematology in POL

Why choose HORIBA?

Compact, reliable, and efficient — HORIBA streamlines your lab with hematology and chemistry in one solution.



Lab Solutions Combo

Including:

- Pentra C400 Chemistry Analyzer
- Micros 60 Hematology Analyzer
- Lite^{DM} Data Management



Contact your Henry Schein Representative to learn more about the lab solutions combo.

CAPITAL LABORATORY

HSI#	MFR#	Description	Size
(106-5492)	1210906014	Alphalyse 360 Reagent 360mL For ABX/Micros 60/Micros CRP 200	1/Bt
(110-4172)	1220001932	Pentra C400/C200 TP: Total Protein Reagent	1/Bx
(338-9153)	1210903010	Pentraclean Cleaning Solution For Pentra	Ea
(338-5479)	1210901020	Pentra 60/80 Sample Diluent 20L	Ea
(108-9175)	1210906013	Lysebio Lytic Reagent 400mL For ABX Pentra 60	Ea
(117-9574)	1220001891	Segment Cuvette For Pentra 400 Clinical Chemistry Analyzer	450/Ca
(113-3181)	1220001901	Reference Reagent 280mL RTU For ABX Pentra	Ea
(110-4171)	1220001871	Pentra C400/C200 LDH: Lactated Dehydrogenase Reagent	1/Bx
(110-4170)	5301101890	Carbon Dioxide Level 1-2 Control 8x3mL RTU For Pentra C400/C200 Liquid Volume	8/Bx
(139-7173)	1300054414	ABX Pentra N MultiControl Multi Multi Control Kit 10x5mL Vial	10/Bx
(140-3380)	1300054415	ABX Pentra P MultiControl Multi Control Kit 10x5mL Vial	10/Bx
(690-8400)	5300000276	Minocal Multi-Analyte Hematology Calibrator 2x2.5mL	Kt
(131-1931)	1300010862	ABX Micros ESV 60 Minipack LMG Reagent	1/Bx
(128-5209)	1300017438	Precitest Quality Control pH 7.4 Solution 15mL For ABX/C200/C400 Chemistry Analyzer	Ea
(128-4930)	1203601326	Lamp For Pentra C400	Ea
(117-4989)	1229001072	Liquid Complete Filter For Pentra 400	Ea
(129-6879)	1300032565	mALB: Microalbumin Calibrator For ABX Pentra 3x2mL Vial	Ea
(133-4779)	5300000944	Transferrin Reagent Assay	Ea
(108-2474)	5300100135	CL020 Hematology Linearity For Micros 45/ Micros	14VI/Bx
(133-4776)	5300000945	Multi Analyte Calibrator For Transferrin/RF/C3/C4	Ea
(128-4959)	1209119847	Probe Assembly Needle Sampler For Pentra C400 Sample	Ea
(117-9276)	5300001371	Micros 60 Lite Hematology Analyzer Complete	Ea
(129-7842)	5360600241	ABX Micros Hematology Analyzer ES60	Ea
(146-5703)	5300001371	Micros 60 Lite Hematology Analyzer Complete	Ea
(690-0015)	1220001626	Pentra C400/C200 ALP: Alkaline Phosphatase Reagent	Ea
(690-0085)	1220001755	Pentra C400/C200 Cleaning Solution 30mL	Ea
(118-8518)	1220001933	Pentra C400/C200 Creatinine CP Reagent	Ea
(435-1352)	1210403010	Miniclean Mini Clean Reagent 1L For ABX/Micros	Ea
(806-8290)	1210802010	Minidil LMG Diluent 10L Cube For Micros	1Ea/Ca
(338-8145)	1210906003	Basolyse II Reagent 1L For ABX/Pentra	Ea
(690-0074)	1220001717	Pentra C400/C200 Solution I Standard 280mL	Ea
(690-0029)	1220001639	Pentra C400/C200 Total Bilirubin Reagent	Ea
(110-4928)	1220001645	Pentra C400/C200 Carbon Dioxide Reagent 2x20mL	Ea
(338-4034)	1210206010	Eosinofix Plus Hematology Reagent 1L For ABX Pentra 60/Pentra XL 80 Plasma (Sodium Citrate)	Ea
(147-6447)	5300100312	Minotrol Low/ Normal/ High Control 12x2.5mL	Ea
(690-0044)	1220001667	Pentra C400/C200 Glucose Hexokinase Reagent	Ea
(690-0016)	1220001627	Pentra C400/C200 ALT: Alanine Aminotransferase Reagent	Ea
(690-0018)	1220001629	Pentra C400/C200 AST: Aspartate Aminotransferase Reagent	Ea
(690-0031)	1220001641	Pentra C400/C200 Urea Glutamate Dehydrogenase Reagent	Ea
(114-5825)	1210401005	Minoclair Bleach Reagent 500mL For ABX/ Pentra/ Micros Series	Ea
(133-7605)	5300000206	Zebra Printer Thermal Label	1/Pk
(690-0084)	1220001754	Pentra C400/C200 Deproteinizer Reagent	Ea
(690-0267)	1220001954	Calcium Reagent 90mL	90mL
(109-2216)	1220001664	Pentra C400/C200 ALB: Albumin Reagent	Ea
(117-6041)	5300100143	Reagent Grade Water For Pentra C400/C200 20liter	20liter
(690-0091)	1220001769	Pentra C400/C200 Etching CP Solution 25mL	Ea
(690-0026)	1220001636	Pentra C400/C200 HDL Direct CP Reagent	Ea
(690-0034)	1220001646	Pentra C400/C200 Magnesium Reagent 2x25mL	2X25mL
(690-0075)	1220001718	Pentra C400/C200 Solution II Standard 100mL Bottle Liquid	Ea
(690-0030)	1220001640	Pentra C400/C200 Triglyceride CP Reagent	Ea
(147-6440)	5300100300	Minocal Calibrator For Whole Blood Hematology 1x2.5mL	Ea
(690-0025)	1220001635	Pentra C400/C200 Direct Bilirubin Reagent	Ea
(690-0024)	1220001634	Pentra C400/C200 Cholesterol CP Reagent	Ea
(690-0046)	1220001670	Pentra C400/C200 Uric Acid CP Reagent	Ea
(110-4929)	1220001648	Pentra C400/C200 Carbon Dioxide Calibrator 3x3mL RTU Liquid Volume	Ea
(690-0022)	1220001632	Pentra C400/C200 CK NAC CP Reagent	Ea
(690-0000)	1220001611	Pentra C400/C200 CRP Reagent	Ea
(110-1239)	5300000502	Pentra Diffrol Multi-Analyte Tri-Level Control 12x3mL Liquid Volume	12/Bx
(690-0028)	1220001638	Pentra C400/C200 dLDL: Low Density Lipoprotein Reagent 100 Count Bottle	Ea
(690-0042)	1220001665	Pentra C400/C200 Phosphorus CP Reagent 100 Count	Ea
(138-5135)	1220001971	Cleaner Cassette For Pentra 400 90mL ISE	Ea
(690-0019)	1220001630	Pentra C400/C200 GGT: Gamma Glutamyl Transferase Reagent 250 Count Bottle CP	Ea
(690-0017)	1220001628	Pentra C400/C200 AMYL: Amylase Reagent 125 Count Bottle	Ea
(348-4269)	5300000490	Diffrol Multi-Analyte Level 1-3 Control 3x3mL	3/Kt
(690-0021)	1220001631	Pentra C400/C200 Lipase CP Reagent 100 Count Bottle	Ea
(690-0036)	1220001652	Multi-Cal Multi-Analyte Calibrator For Pentra C400/C200	10X3mL
(690-0115)	1221034626	Reagent Cup For Pentra 400 Clinical Chemistry Analyzer	100/Bag

Care for Your Patients and the Planet



For a Healthier **TOMORROW**

By integrating environmental sustainability into health care practices worldwide, we can support practitioners to enhance both patient care and the environment, fostering a health care system that positively contributes to a healthy planet.

www.henryschein.com/practicegreen





To Order: **1-800-772-4346** 8am–8pm (et)
henryschein.com/medical

*California customers please use website
 when ordering for Prop 65 information.*

Due to manufacturer price increases, prices are subject to change without notification. Offer valid on our full array of products. Upon any request by Medicare, Medicaid, Tricare or other payer, you must report the value of any benefit received from a discount program (e.g., points redemptions, gifts or other special awards). ©2025 Henry Schein, Inc. No copying without permission.

Not responsible for typographical errors. No adjustments from prior sales. Offer not to be combined with other promotions and/or special contract pricing. Void where prohibited. Local regulatory requirements may apply to use or installation of certain products. It is your responsibility to understand and comply with any such requirements prior to purchase, use or installation of products.

[†]For informational purposes only.

Source: ReimbursementCodes.com and/or www.cms.gov. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.

PRSRT STD
 U.S. POSTAGE
PAID
 HENRY SCHEIN, INC.



A FRIENDLY REMINDER TO RECYCLE—
 Thank you for helping us care for the environment.



Equipment Finance Made Easy

**CUSTOMIZED FINANCING SOLUTIONS TO EXPAND,
 UPGRADE, AND ENHANCE YOUR PRACTICE**

OUR FINANCING BENEFITS INCLUDE:

- ✓ 100% Financing*
- ✓ Up to 6-month Payment Deferrals
- ✓ Customized Repayment Programs
- ✓ Same-Day Approvals and Documentation
- ✓ Simple, Fast Application Process

PLUS

Take advantage of Section 179 and get a possible tax deduction for up to 100% of your purchase on your 2025 taxes for new equipment, technology, or other qualifying purchases.**



Contact us today to schedule a consultation **to learn more
 about how we can support your growth and success**

Scan the QR code to learn more

Disclaimer:

* Rates and programs are subject to change without notice. All transactions are subject to the satisfaction of underwriting guidelines, credit approval by third-party lenders, and documentation requirements, and not all applicants will qualify. Certain other restrictions and additional terms and conditions may apply. ** Neither Henry Schein, Inc. nor Henry Schein Financial Services provides tax advice. Please consult with a qualified professional tax advisor to determine your eligibility for a Section 179 tax deduction. Must purchase and place into service by December 31, 2025. ©2025 Henry Schein, Inc. Neither Henry Schein, Inc. nor Henry Schein Financial Services provides financial advice. Please consult your financial advisor. Neither Henry Schein, Inc. nor Henry Schein Financial Services is a bank and neither represents itself as such, nor conducts banking activities. Henry Schein Financial Services may receive a marketing fee from the vendor for products/services purchased. "World's Most Ethical Companies" and "Ethisphere" names and marks are registered trademarks of Ethisphere LLC.