

FROM HERE

TO HERE

WEAR IT WITH
CONFIDENCE.

by **HENRY SCHEIN**[®]

ORDER NOW AT
HENRYSCHEIN.COM/CRITERION

Gloves Comparison Chart

GLOVE NAME	TYPE	COLOR	STERILE STATUS	NOT MADE WITH NATURAL RUBBER LATEX	SURFACE TEXTURE	THICKNESS AT PALM (MM/MIL)	ADDITIONAL FEATURES	ITEM CODES
NITRILE								
CRITERION™ Nitrile Exam N100 Gloves	Exam	Blue	Non-Sterile	Yes	Finger Textured	0.08/3.15	Chemotherapy tested to ASTM D6978 standard.	100/Box: XS (111-8533), S (111-8535), M (111-8536), L (111-8537) 90/Box: XL (111-8539)
CRITERION™ Nitrile Exam N200 Gloves	Exam	Periwinkle Blue	Non-Sterile	Yes	Finger Textured	0.06/2.36	Chemotherapy tested to ASTM D6978 standard.	200/Box: XS (900-7437), S (900-7438), M (900-7439), L (900-7440) 180/Box: XL (900-7441)
	Exam	Sky Blue	Non-Sterile	Yes	Finger Textured	0.06/2.36	Chemotherapy tested to ASTM D6978 standard.	200/Box: XS (570-2097), S (570-2098), M (570-2099), L (570-2100) 180/Box: XL (570-2101)
	Exam	Black	Non-Sterile	Yes	Finger Textured	0.06/2.36	Fentanyl tested. Chemotherapy tested to ASTM D6978 standard.	200/Box: XS (570-3487), S (570-3488), M (570-3489), L (570-3490) 180/Box: XL (570-3491)
CRITERION™ Nitrile Exam N300 Gloves	Exam	Ice Blue	Non-Sterile	Yes	Finger Textured	0.05/1.97	Chemotherapy tested to ASTM D6978 standard. Eco-friendly Packaging.	300/Box: XS (112-6829), S (112-6830), M (112-6831), L (112-6832) 250/Box: XL (112-6833)
	Exam	Electric Blue	Non-Sterile	Yes	Finger Textured	0.05/1.97	Fentanyl tested. Chemotherapy tested to ASTM D6978 standard. Eco-friendly Packaging.	300/Box: XS (570-3540), S (570-3541), M (570-3542), L (570-3543) 250/Box: XL (570-3544)
	Exam	Black	Non-Sterile	Yes	Finger Textured	0.05/1.97	Fentanyl tested. Chemotherapy tested to ASTM D6978 standard. Eco-friendly Packaging.	300/Box: XS (570-3482), S (570-3483), M (570-3484), L (570-3485) 250/Box: XL (570-3486)

GLOVE NAME	TYPE	COLOR	STERILE STATUS	NOT MADE WITH NATURAL RUBBER LATEX	SURFACE TEXTURE	THICKNESS AT PALM (MM/MIL)	ADDITIONAL FEATURES	ITEM CODES
NITRILE (CONT.)								
CRITERION™ Nitrile Exam Colloidal Oatmeal Gloves	Exam	White	Non-Sterile	Yes	Finger Textured	0.05/1.97	Colloidal Oatmeal Coated.	<u>200/Box:</u> XS (112-6220), S (112-6221), M (112-6222), L (112-6223), XL (112-6224)
CRITERION™ Nitrile Exam Low Derma Gloves	Exam	Slate Blue	Non-Sterile	Yes	Finger Textured	0.05/1.97	Fentanyl tested. Chemotherapy tested to ASTM D6978 standard. Low Dermatitis Potential.	<u>200/Box:</u> XS (570-0630), S (570-0631), M (570-0632), L (570-0633) <u>180/Box:</u> XL (570-0634)
CRITERION™ Nitrile Exam Aloe + Vitamin E Gloves	Exam	Magenta	Non-Sterile	Yes	Finger Textured	0.05/1.97	Aloe Vera + Vitamin E Coated. Polymer Coated.	<u>200/Box:</u> XS (112-7143), S (112-7144), M (112-7146), L (112-7147)
	Exam	Lavender	Non-Sterile	Yes	Finger Textured	0.05/1.97	Aloe Vera + Vitamin E Coated. Polymer Coated.	<u>200/Box:</u> XS (570-3545), S (570-3546), M (570-3547), L (570-3548), XL (570-3549)
CRITERION™ Nitrile Mint Scented	Exam	Blue	Non-Sterile	Yes	Finger Textured	0.08/3.15	Mint Scented.	<u>100/Box:</u> XS (570-3535), S (570-3536), M (570-3537), L (570-3538), XL (570-3539)
POLYCHLOROPRENE								
CRITERION™ Chloroprene Exam Ultrathin Gloves	Exam	Lime	Non-Sterile	Yes	Finger Textured	0.05/1.97	Polymer Coated.	<u>100/Box:</u> XS (112-7306), S (112-7307), M (112-7308), L (112-7309), XL (112-7310)
CRITERION™ Chloroprene Exam Gloves	Exam	Green	Non-Sterile	Yes	Finger Textured	0.08/3.15	Polymer Coated.	<u>100/Box:</u> XS (112-6857), S (112-6858), M (112-6859), L (112-6860) <u>90/Box:</u> XL (112-6861)
	Exam	Pink	Non-Sterile	Yes	Finger Textured	0.08/3.15	Polymer Coated.	<u>100/Box:</u> XS (570-3532), S (570-3533), M (570-3534), L (570-3572) <u>90/Box:</u> XL (570-3573)
CRITERION™ Polychloroprene Surgical Gloves	Surgical	White	Sterile	Yes	Palm to Fingertip Textured	0.15/5.91	Polymer Coated. Extended Cuff.	<u>50 Pair/Box:</u> 5.5 (112-7079), 6 (112-7080), 6.5 (112-7081), 7 (112-7082), 7.5 (112-7083), 8 (112-7085), 8.5 (112-7086), 9 (112-7088)
CRITERION™ Polyisoprene Surgical Gloves	Surgical	White	Sterile	Yes	Palm to Fingertip Textured	0.21/8.27	Polymer and Polyurethane Coated. Extended Cuff.	<u>50 Pair/Box:</u> 5.5 (570-1119), 6 (570-1120), 6.5 (570-1121), 7 (570-1122), 7.5 (570-1123), 8 (570-1124), 8.5 (570-1125), 9 (570-1126)
LATEX								
CRITERION™ Latex Aloe & Vitamin E Gloves	Exam	White	Non-Sterile	No	Finger Textured	0.08/3.15	Aloe Vera + Vitamin E Coated. Polymer Coated. Chlorinated.	<u>100/Box:</u> XS (570-3448), S (570-3449), M (570-3450), L (570-3451) <u>90/Box:</u> XL (570-3452)
CRITERION™ Latex Exam Polymer-Coated Gloves	Exam	White	Non-Sterile	No	Finger Textured	0.08/3.15	Polymer Coated. Chlorinated.	<u>100/Box:</u> XS (570-3412), S (570-3413), M (570-3414), L (570-3415) <u>90/Box:</u> XL (570-3416)
CRITERION™ Latex Exam Chlorinated Gloves	Exam	White	Non-Sterile	No	Palm to Fingertip Textured	0.11/4.33	Chlorinated.	<u>100/Box:</u> XS (102-5354), S (102-8296), M (102-5355), L (102-7763), XL (102-5357)
CRITERION™ Latex Extra Grip Gloves	Exam	White	Non-Sterile	No	Palm to Fingertip Textured	0.10/4.33	Polymer Coated. Non Chlorinated.	<u>100/Box:</u> XS (570-3527), S (570-3528), M (570-3529), 7 (570-4498), 7.5 (570-4499) <u>90/Box:</u> XL (570-3531)
CRITERION™ Latex Surgical Gloves	Surgical	White	Sterile	No	Smooth	0.14/5.51	Polymer Coated. Extended Cuff. Non Chlorinated.	<u>50 Pair/Box:</u> 5.5 (570-4495), 6 (570-4496), 6.5 (570-4497), 7 (570-4498), 7.5 (570-4499), 8 (570-4500), 8.5 (570-4501), 9 (570-4502)
VINYL								
CRITERION™ Vinyl Exam Gloves	Exam	Clear	Non-Sterile	Yes	Smooth	0.08/3.15	Smooth. Accelerator Free.	<u>100/Box:</u> XS (570-0087), S (570-0088), M (570-0089), L (570-0090) <u>90/Box:</u> XL (570-0091)

CRITERION™ NITRILE EXAM N100 GLOVES

**BUILT FOR
DURABILITY,
ENGINEERED
FOR COMFORT.**



SAME GREAT PRODUCTS, NEW PACKAGING



For more information about the full Criterion line of personal protection equipment scan here or visit HENRYSCHEIN.COM/CRITERION

CRITERION™ NITRILE EXAM N100 GLOVES

- Durable nitrile construction
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Beaded cuff. 3.15-mil/0.08mm thickness at palm



BLUE COLOR

X-Small	100/box	111-8533
Small	100/box	111-8535
Medium	100/box	111-8536
Large	100/box	111-8537
X-Large	90/box	111-8539

CHEMOTHERAPY TESTED TO ASTM D6978 STANDARD

Caution: gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of chemicals used. This glove has been tested for resistance to permeation of various chemotherapy drugs per ASTM D6978, "standard practice for assessment of resistance of medical gloves to permeation by chemotherapy drugs."

CHEMOTHERAPY DRUGS & CONCENTRATION*	MINIMUM BREAKTHROUGH DETECTION TIME
Cisplatin (1.0mg/mL)	>240min
Cyclophosphamide (Cytoxan) (20mg/mL)	>240min
Dacarbazine (DTIC) (10.0mg/mL)	>240min
Doxorubicin Hydrochloride (2.0mg/mL)	>240min
Etoposide (Toposar) (1.0mg/mL)	>240min
5-Fluorouracil (50.0mg/mL)	>240min
Glutaraldehyde (80.0mg/mL)	>240min
Methotrexate (100.0mg/mL)	>240min
Paclitaxel (Taxol) (6.0mg/mL)	>240min
Vincristine Sulfate (1.0mg/mL)	>240min

*When chemotherapy drugs are present, gloves selection should be based on the specific types of drugs used. Users should review drug labeling or material safety data sheets for the drugs being used to determine an adequate level of protection.

SAME GREAT PRODUCTS, NEW PACKAGING



1. UPDATED PACKAGING

Makes it easy to know what's inside the box with attractive imagery and bold design elements.



2. INFO ZONE

Appears in the bottom right of all packaging. The icons indicate key product information: PPE standards, size, color and quantity.



3. REORDER PANEL

Helps to simplify the reordering process by clearly stating the full product name and reorder item number.

CRITERION™ NITRILE EXAM N200 GLOVES

**THE PERFECT PAIR:
DURABILITY AND
PERFORMANCE**



SAME GREAT PRODUCTS, NEW PACKAGING



For more information about the full Criterion line of personal protection equipment
scan here or visit [HENRYSCHEIN.COM/MEDICAL](https://www.henryschein.com/medical)

CRITERION™ NITRILE EXAM N200 GLOVES

- Durable nitrile provides reliable protection and tactility
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Beaded cuff. 2.36-mil / 0.06mm thickness at palm



PERIWINKLE BLUE COLOR

● X-Small	200/box	900-7437
● Small	200/box	900-7438
● Medium	200/box	900-7439
● Large	200/box	900-7440
● X-Large	180/box	900-7441

SKY BLUE COLOR

● X-Small	200/box	570-2097
● Small	200/box	570-2098
● Medium	200/box	570-2099
● Large	200/box	570-2100
● X-Large	180/box	570-2101

BLACK COLOR

● X-Small	200/box	570-3487
● Small	200/box	570-3488
● Medium	200/box	570-3489
● Large	200/box	570-3490
● X-Large	180/box	570-3491

CHEMOTHERAPY TESTED TO ASTM D6978 STANDARD

Caution: gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of chemicals used. This glove has been tested for resistance to permeation of various chemotherapy drugs per ASTM D6978, "standard practice for assessment of resistance of medical gloves to permeation by chemotherapy drugs."

CHEMOTHERAPY DRUGS & CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME
Carboplatin (10.0mg/ml)	>240min
Cisplatin (1.0mg/mL)	>240min
Cyclophosphamide (Cytosan) (20mg/mL)	>240min
Dacarbazine (DTIC) (10.0mg/mL)	>240min
Doxorubicin Hydrochloride (2.0mg/mL)	>240min
Etoposide (Toposar) (20.0mg/mL)	>240min
Fluorouracil (50.0mg/mL)	>240min
Ifosamide (50mg/ml)	>240min
Methotrexate (25mg/mL)	>240min
Mitomycin C (0.5mg/mL)	>240min
Mitroxantrone (2.0mg/mL)	>240min
Paclitaxel (Taxol) (6.0mg/mL)	>240min
Vincristine Sulfate (1.0mg/mL)	>240min

Warning: Do not use Carmustine (BCNU) (3.3mg/mL) and Thiotepa (10mg/mL)
When chemotherapy drugs are present, gloves selection should be based on the specific types of drugs used. Users should review drug labeling or material safety data sheets for the drugs being used to determine an adequate level of protection.

FENTANYL-TESTED

Tested using ASTM D6978 test method, which evaluates glove permeation resistance to fentanyl. The results show that gloves have resistance to fentanyl for at least 240 minutes.

Fentanyl breakthrough time (100 mcg/2ml): >240 minutes

CRITERION™ NITRILE EXAM N300 GLOVES

ULTRATHIN FORMULATION

DELIVERED IN
ECO-FRIENDLY
PACKAGING



SAME GREAT PRODUCTS, NEW PACKAGING



For more information about the full Criterion line of personal protective equipment,
scan here or visit HENRYSCHEIN.COM/CRITERION

CRITERION™ NITRILE EXAM N300 GLOVES

- Ultrathin formulation enhances tactile sensitivity and dexterity
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Beaded cuff. 2.0-mil/0.05mm thickness at palm
- Eco-friendly packaging, with 300 gloves in a box, reduces waste and saves space



ICE BLUE COLOR

X-Small	300/box	112-6829
Small	300/box	112-6830
Medium	300/box	112-6831
Large	300/box	112-6832
X-Large	250/box	112-6833

ELECTRIC BLUE COLOR

X-Small	300/box	570-3540
Small	300/box	570-3541
Medium	300/box	570-3542
Large	300/box	570-3543
X-Large	250/box	570-3544

BLACK COLOR

X-Small	300/box	570-3482
Small	300/box	570-3483
Medium	300/box	570-3484
Large	300/box	570-3485
X-Large	250/box	570-3486

CHEMOTHERAPY TESTED TO ASTM D6978 STANDARD

Caution: gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of chemicals used. This glove has been tested for resistance to permeation of various chemotherapy drugs per ASTM D6978, "standard practice for assessment of resistance of medical gloves to permeation by chemotherapy drugs."

CHEMOTHERAPY DRUGS & CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME
Cisplatin (1.0mg/mL)	>240min
Cyclophosphamide (Cytoxan) (20mg/mL)	>240min
Dacarbazine (DTIC) (10.0mg/mL)	>240min
Doxorubicin Hydrochloride (2.0mg/mL)	>240min
Etoposide (Toposar) (20.0mg/mL)	>240min
Fluorouracil (50.0mg/mL)	>240min
Methotrexate (25mg/mL)	>240min
Mitomycin C (0.5mg/mL)	>240min
Paclitaxel (Taxol) (6.0mg/mL)	>240min
Vincristine Sulfate (1.0mg/mL)	>240min

Warning: Do not use Carmustine (BCNU) (3.3mg/mL) and Thiotepa (10mg/mL)
When chemotherapy drugs are present, gloves selection should be based on the specific types of drugs used. Users should review drug labeling or material safety data sheets for the drugs being used to determine an adequate level of protection.

FENTANYL-TESTED (BLACK GLOVES ONLY)

Tested using ASTM D6978 test method, which evaluates glove permeation resistance to fentanyl. The results show that gloves have resistance to fentanyl for at least 240 minutes.

Fentanyl breakthrough time (100 mcg/2mL): >240 minutes