



# DASH<sup>®</sup> Rapid PCR System

Rapid PCR Made Simple<sup>™</sup>

## Lab-Quality Results in 15 Minutes at the Point of Care

- ✓ **SIMPLE**  
Intuitive, easy-to-operate user experience
- ✓ **RAPID**  
Delivers results in 15 minutes, with only one minute of hands-on time
- ✓ **AFFORDABLE**  
Low-cost solution enables PCR testing at the point of care
- ✓ **ACCURATE**  
Lab-quality results enhance confidence in treatment decisions

## Test and Treat in a Single Patient Visit



PATIENT ARRIVES AND  
SAMPLE IS COLLECTED



RUN DASH<sup>®</sup> SARS-COV-2  
& FLU A/B TEST



RESULTS IN JUST  
15 MINUTES



CLINICIAN PRESENTS RESULTS  
AND DISCUSSES TREATMENT PLAN

## DASH<sup>®</sup> Rapid PCR System Capabilities



Cartridges are stored at  
room temperature



Designed to support a wide  
variety of sample types



Intuitive operation by non-  
laboratory personnel



No calibration or preventive  
maintenance required



Expansive planned test menu, including  
both viral and bacterial pathogens




Built to be multiplex-capable  
in a single cartridge



## Why Choose the DASH® Rapid PCR System?

- ✓ Rapid PCR test results at the point of care
- ✓ Lab-quality sensitivity and specificity
- ✓ Anterior nasal swab specimens
- ✓ Qualitative, intuitive patient report
- ✓ Wi-Fi and cloud connectivity enabled



## DASH® SARS-CoV-2 & Flu A/B Test

FDA Cleared & CLIA Waived | 15-Minute Runtime | PCR at the POC

### Results From a Prospective Clinical Study Against an FDA-Cleared PCR Test

| POSITIVE AGREEMENT<br><i>(Sensitivity)</i> |       | NEGATIVE AGREEMENT<br><i>(Specificity)</i> |       |
|--------------------------------------------|-------|--------------------------------------------|-------|
| SARS-CoV-2                                 | 95.2% | SARS-CoV-2                                 | 99.5% |
| Flu A                                      | 94.3% | Flu A                                      | 98.1% |
| Flu B                                      | 97.3% | Flu B                                      | 99.2% |

Source: DASH® SARS-CoV-2 & Flu A/B Test package insert.

The DASH® SARS-CoV-2 & Flu A/B Test is affordably priced to support high-quality PCR testing at the point of care. The test utilizes the existing CPT code 87636-QW.